

Artwork below shown at approximately 145% of actual size



Department of Agriculture
STATE OF HAWAII

LICENSED

PERIOD 2014-2016 LIC. NO.

8013.62



Department of Agriculture
STATE OF HAWAII

LICENSED

PERIOD 2011-2013 LIC. NO.

Kills *C. difficile* Spores

DISPATCH[®] BY



Bleach Germicidal Cleaner

- Contains 6500 ppm Sodium Hypochlorite
- Ready-to-Use 1:10 Dilution of Bleach and Detergent
- Cleans, Disinfects & Deodorizes in One Step
- For Hospital, Veterinary & Food Contact Hard, Non-porous Surfaces
- For Instrument Pre-cleaning & Disinfection of Non-critical Instruments

Kills 99.99% of bacteria, viruses[†] & fungi in 1 MINUTE on hard, non-porous surfaces

Kills *Mycobacterium bovis* - BCG, (TB) in 30 SECONDS

Kills *Clostridium difficile* spores in 5 MINUTES*

ACTIVE INGREDIENT:	
Sodium Hypochlorite.....	0.65%
OTHER INGREDIENTS:.....	99.35%
TOTAL:.....	100.00%

KEEP OUT OF REACH OF CHILDREN

CAUTION: READ COMPLETE DIRECTIONS AND PRECAUTIONS ON PACKAGE REAR PANELS.

32 FL OZ (1 QT) 946 mL

NI-15592

R0884-1.1

Artwork below shown at approximately 145% of actual size

OPEN HERE
FOR DIRECTIONS

Multipurpose, broad spectrum cleaner disinfectant for use in healthcare and other areas where control of cross contamination is required, especially during outbreaks of *Clostridium difficile* and other Healthcare Associated Infections (HAIs) listed on this label. Meets OSHA Bloodborne Pathogen Standards. May be used for instrument precleaning and disinfection of noncritical instruments. **Surfaces:** May be used on hard, nonporous surfaces such as stainless steel, glazed ceramic tile, hard plastics, nonporous vinyl, painted surfaces, Plexiglas, laminated plastic countertops, fiberglass (sealed) surfaces, glazed porcelain and glass. **APPLICATIONS: Medical** – hard, nonporous surfaces such as autoclaves, blood glucose monitors, external surfaces of medical equipment, IV poles, exam tables, bed railings, external surfaces of dialysis machines, gurneys, mammography equipment, plastic mattress covers, phlebotomy trays, reception counters, external surfaces of respiratory equipment, stretchers, external surfaces of ultrasound transducers and probes, wash basins and wheelchairs. **Dental** – chairs, counters, operatory hard surfaces and light lens covers. **Veterinary** – animal equipment, cages, fountains, vet clinic surfaces. **Other** – changing tables, doorknobs, hand railings, highchairs, playpens. **Food service** – both food and non-food contact surfaces. **USE SITES: Medical** – ambulance, central service, critical care, dialysis clinics, emergency rooms, exam rooms, home health care, hospitals, intensive care, laboratories, neonatal nursery, nurses' stations, nursing homes, operating rooms, orthopedic clinics, outpatient clinics, patient rooms, physical therapy rooms, physician offices, radiology, recovery rooms, respiratory therapy and surgery. **Dental** – dental clinics, labs. **Veterinary** – animal housing, veterinary clinics/hospitals. **Other** – blood banks, correctional facilities, day care, funeral homes, health clubs, recreational centers, restrooms and schools. Efficacy tests have demonstrated that this product is an effective bactericide, virucide, tuberculocide, and fungicide in the presence of organic soil (5% blood serum).

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling. **FIRST AID: If in eyes:** Hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor or going for treatment. **PHYSICAL OR CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers or toilet bowl cleaners. This product is not intended for use on fabrics.

STORAGE AND DISPOSAL: Store this product in a cool, dry area away from direct sunlight and heat. Refillable container. Refill empty DISPATCH spray container only from DISPATCH refill bottle. Cleaning before refilling is the responsibility of the refiller. Take care to update expiration date on refillable bottle and replace bottle if it becomes worn or damaged. Do not reuse this container for any other purpose. Rinse container with water prior to disposal. Cleaning the container before final disposal is the responsibility of the person disposing of the container. When discarding container, place in trash collection or take to recycling center. Do not contaminate water, food or feed by storage and disposal.

EPA Reg. No. 56392-7
EPA Est. Nos. 48181-IL-1, 50757-WI-1, 62801-MI-1,
8251-OH-1, 5813-GA-2
Actual EPA Est. No. near batch code.

Manufactured for:
©2010 Clorox Professional Products Company
1221 Broadway, Oakland, CA 94612
Visit us at www.cloroxprofessional.com 800-234-7700
® Trademark of Clorox Professional Products Company



R0884-1.1
NI-15592

Artwork below shown at approximately 145% of actual size

HINGED AREA, NO PRINT ALLOWED

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

To Clean, Disinfect and Deodorize in One Step: Spray 6"–8" from surface until surface is completely wet or apply with a clean cloth, mop or sponge saturated with DISPATCH until surface is completely wet. (When applying with a saturated cloth or sponge gloves should be worn.) To disinfect, let stand for one minute. To kill TB, let stand for 30 seconds at room temperature (20°C–25°C). Wipe with a clean, damp cloth or paper towel or allow to air dry. For food contact surfaces, rinse with potable water. For all other surfaces no rinsing is required. Gross soil must be removed prior to disinfecting. **To Clean:** Follow directions above, but no standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel. ***Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores.** **Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks, or eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Spray 6" – 8" from surface until surface is completely wet. Let stand for 5 minutes. Wipe with a clean, damp cloth or paper towel or allow to air dry. **To Clean and Disinfect in a Veterinary Application:** Use to clean and disinfect hard, nonporous surfaces such as feeding and watering equipment, cages, utensils, instruments, kennels, stables, catteries, etc. Remove all animals and feed from premises, animal transportation vehicles, crates, etc. Remove all litter, droppings and manure from walls, floors and surfaces of facilities occupied or traversed by animals. Empty all feeding and watering equipment. Preclean all surfaces with soap or detergent and rinse with water. Saturate surfaces with DISPATCH and let stand for one minute. To kill TB, let stand for 30 seconds at room temperature. Rinse all surfaces with potable water. Ventilate buildings and other closed spaces. Do not house animals or employ equipment until treated surfaces have been thoroughly rinsed with water and allowed to dry. Thoroughly scrub all treating, feeding and watering appliances with soap or detergent, and rinse with potable water before reuse. **To Preclean Instruments Prior to Terminal Sterilization/High-Level Disinfection: As a precleaning spray** – Place instruments into a suitable container. Spray DISPATCH onto instruments to thoroughly wet all surfaces. Let stand for up to 10 minutes. Rinse instruments. **As a precleaning immersion solution** – Fill appropriate size container with a sufficient amount of DISPATCH to completely submerge instruments. Place instruments into the container of DISPATCH, cover and allow to soak for up to 10 minutes. Remove and rinse, and follow with an appropriate cleaning and disinfecting process. Change solution daily. **As a manual instrument cleaner** – Thoroughly prerinse dirty instruments under running water to remove visible, gross debris. Immerse prerinsed instruments into an appropriate size container filled with DISPATCH. Scrub instruments using a stiff bristle brush until visibly clean. Submerge instruments while scrubbing. Rinse instruments thoroughly. Change solution daily. Follow with an appropriate disinfection process. **Cleaning of critical and semicritical devices must be followed by an appropriate terminal sterilization/high-level disinfection process.** R0884-1.1

Artwork below shown at approximately 145% of actual size

To Disinfect Noncritical, Precleaned Instruments:

Instruments must be thoroughly precleaned to remove excess organic debris, rinsed and rough dried. (Clean and rinse lumens of hollow instruments before filling with DISPATCH or before immersion.) **Immersion method** – Using a soaking tray, immerse instruments into DISPATCH and let stand for one minute. To kill TB, immerse for 30 seconds at room temperature. Remove and wipe with a clean, damp cloth or paper towel and allow to air dry. Change solution for each use. **Spray method** – Spray all surfaces of instruments with DISPATCH until thoroughly wet. Let stand for one minute. To kill TB, let stand for 30 seconds at room temperature. To kill *Clostridium difficile* (spores), let stand for five minutes. Wipe with a clean, damp cloth or paper towel and allow to air dry.

†Kills HIV, HBV and HCV on precleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV). **Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV and HCV on Surfaces/Objects Soiled with Blood/Body Fluids.**

Personal Protection: When handling items soiled with blood/body fluids, personal protective equipment such as disposable gloves, gowns, masks and eye coverings should be used.

Cleaning Procedure: Blood/body fluids must be thoroughly cleaned from surfaces and objects before disinfecting. **Contact Time:** Let stand for one minute. **Disposal of Infectious Materials:** Blood, body fluids and cleaning materials must be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semicritical medical devices prior to sterilization or high-level disinfection.

Kills bacteria, viruses† & fungi in 1 MINUTE:

Bacteria:

- *Acinetobacter baumannii*
- *Enterobacter aerogenes*
- *Enterococcus faecium* (VRE)
- *Escherichia coli* ESBL
- *Klebsiella pneumoniae*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Staphylococcus aureus*
- *Staphylococcus aureus* (MRSA)
- *Streptococcus pyogenes*

Fungi:

- *Trichophyton mentagrophytes*
(Athlete's Foot Fungus)

Viruses:

- †Avian Influenza A Virus
- †Canine Parvovirus
- †Feline Panleukopenia Virus
- †Hepatitis A Virus (HAV)
- †Hepatitis B Virus (HBV)
- †Hepatitis C Virus (HCV)

- †Herpes Simplex Virus (HSV-2)
- †HIV-1 (AIDS Virus)
- †Influenza A H1N1
- †Influenza A Virus
- †Norovirus
- †Poliovirus Type 1 (Mahoney)
- †Rhinovirus
- †Rotavirus

Kills *Mycobacterium bovis* - BCG, (TB) in 30 SECONDS

Kills *Clostridium difficile* spores in 5 MINUTES*

HINGED AREA, NO PRINT ALLOWED

R0884-1.1

LICENSED

Artwork below shown at approximately 110% of actual size

PERIOD 2014-2016 LIC. NO.

8013.62

LICENSED

PERIOD 2011-2013 LIC. NO.

Kills *C. difficile* Spores

DISPATCH[®] BY



Bleach Germicidal Cleaner

- Contains 6500 ppm Sodium Hypochlorite
- Ready-to-Use 1:10 Dilution of Bleach and Detergent
- Cleans, Disinfects & Deodorizes in One Step
- For Hospital, Veterinary & Food Contact Hard, Non-porous Surfaces
- For Instrument Pre-cleaning & Disinfection of Non-critical Instruments

Kills 99.99% of bacteria, viruses[†] & fungi in 1 MINUTE on hard, non-porous surfaces

Kills *Mycobacterium bovis* - BCG, (TB) in 30 SECONDS

Kills *Clostridium difficile* spores in 5 MINUTES*

ACTIVE INGREDIENT:
Sodium Hypochlorite 0.65%
OTHER INGREDIENTS: 99.35%
TOTAL: 100.00%

KEEP OUT OF REACH OF CHILDREN

CAUTION: READ COMPLETE DIRECTIONS AND PRECAUTIONS ON PACKAGE REAR PANELS.

128 FL OZ (1 GAL) 3.78 L

R0885-1.1

NI-15594

Artwork below shown at approximately 110% of actual size

Multipurpose, broad spectrum cleaner/disinfectant for use in healthcare and other areas where control of cross contamination is required, especially during outbreaks of *Clostridium difficile* and other Healthcare Associated Infections (HAIs) listed on this label. Meets OSHA Bloodborne Pathogen Standards. May be used for instrument precleaning and disinfection of noncritical instruments. **Surfaces:** May be used on hard, nonporous surfaces such as stainless steel, glazed ceramic tile, hard plastics, nonporous vinyl, painted surfaces, Plexiglas, laminated plastic countertops, fiberglass (sealed) surfaces, glazed porcelain and glass. **APPLICATIONS:**
Medical – hard, nonporous surfaces such as autoclaves, blood glucose monitors, external surfaces of medical equipment, IV poles, exam tables, bed railings, external surfaces of dialysis machines, gurneys, mammography equipment, plastic mattress covers, phlebotomy trays, reception counters, external surfaces of respiratory equipment, stretchers, external surfaces of ultrasound transducers and probes, wash basins and wheelchairs.
Dental – chairs, counters, operatory hard surfaces and light lens covers. **Veterinary** – animal equipment, cages, fountains, vet clinic surfaces. **Other** – changing tables, doorknobs, hand railings, highchairs, playpens. **Food service** – both food and non-food contact surfaces. **USE SITES: Medical** – ambulance, central service, critical care, dialysis clinics, emergency rooms, exam rooms, home health care, hospitals, intensive care, laboratories, neonatal nursery, nurses' stations, nursing homes, operating rooms, orthopedic clinics, outpatient clinics, patient rooms, physical therapy rooms, physician offices, radiology, recovery rooms, respiratory therapy and surgery.
Dental – dental clinics, labs. **Veterinary** – animal housing, veterinary clinics/hospitals. **Other** – blood banks, correctional facilities, day care, funeral homes, health clubs, recreational centers, restrooms and schools. Efficacy tests have demonstrated that this product is an effective bactericide, virucide, tuberculocide, and fungicide in the presence of organic soil (5% blood serum).

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

To Clean, Disinfect and Deodorize in One Step: Spray 6"–8" from surface until surface is completely wet or apply with a clean cloth, mop or sponge saturated with DISPATCH until surface is completely wet. (When applying with a saturated cloth or sponge gloves should be worn.) To disinfect, let stand for one minute. To kill TB, let stand for 30 seconds at room temperature (20–25°C). Wipe with a clean, damp cloth or paper towel or allow to air dry. For food contact surfaces, rinse with potable water. For all other surfaces no rinsing is required. Gross soil must be removed prior to disinfecting. **To Clean:** Follow directions above, but no standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel. ***Special Instructions for Cleaning Prior to Disinfection against Clostridium difficile spores. Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks, or eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Spray 6"–8" from surface until surface is completely wet. Let stand for 5 minutes. Wipe with a clean, damp cloth or paper towel or allow to air dry. **To Clean and Disinfect in a Veterinary Application:** Use to clean and disinfect hard, nonporous surfaces such as feeding and watering equipment, cages, utensils, instruments, kennels, stables, catteries, etc. Remove all animals and feed from premises, animal transportation vehicles, crates, etc. Remove all litter, droppings and manure from walls, floors and surfaces of facilities occupied or traversed by animals. Empty all feeding and watering equipment. Preclean all surfaces with soap or detergent and rinse with water. Saturate surfaces with DISPATCH and let stand for one minute. To kill TB, let stand for 30 seconds at room temperature. Rinse all surfaces with potable water. Ventilate buildings and other closed spaces. Do not house animals or employ equipment until treated surfaces have been thoroughly rinsed with water and allowed to dry. Thoroughly scrub all treating, feeding and watering appliances with soap or detergent, and rinse with potable water before reuse. **To Preclean Instruments Prior to Terminal Sterilization/High-Level Disinfection: As a precleaning spray** – Place instruments into a suitable container. Spray DISPATCH onto instruments to thoroughly wet all surfaces. Let stand for up to 10 minutes. Rinse instruments. **As a precleaning immersion solution** – Fill appropriate size container with a sufficient amount of DISPATCH to completely submerge instruments. Place instruments into the container of DISPATCH, cover and allow to soak for up to 10 minutes. Remove and rinse, and follow with an appropriate cleaning and disinfecting process. Change solution daily. **As a manual instrument cleaner** – Thoroughly prerinse dirty instruments under running water to remove visible, gross debris. Immerse prerinsed instruments into an appropriate size container filled with DISPATCH. Scrub instruments using a stiff bristle brush until visibly clean. Submerge instruments while scrubbing. Rinse instruments thoroughly. Change solution daily. Follow with an appropriate disinfection process. **Cleaning of critical and semicritical devices must be followed by an appropriate terminal sterilization/high-level disinfection process. To Disinfect Noncritical, Precleaned Instruments:** Instruments must be thoroughly precleaned to remove excess organic debris,

rinsed and rough dried. (Clean and rinse lumens of hollow instruments before filling with DISPATCH or before immersion.) **Immersion method** – Using a soaking tray, immerse instruments into DISPATCH and let stand for one minute. To kill TB, immerse for 30 seconds at room temperature. Remove and wipe with a clean, damp cloth or paper towel and allow to air dry. Change solution for each use. **Spray method** – Spray all surfaces of instruments with DISPATCH until thoroughly wet. Let stand for one minute. To kill TB, let stand 30 seconds at room temperature. To kill *Clostridium difficile* (spores), let stand for five minutes. Wipe with a clean, damp cloth or paper towel and allow to air dry. **†Kills HIV, HBV and HCV on precleaned environmental surfaces/objects previously soiled with blood/body fluids** in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV). **Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV and HCV on Surfaces/Objects Soiled with Blood/Body Fluids. Personal Protection:** When handling items soiled with blood/body fluids, personal protective equipment such as disposable gloves, gowns, masks and eye coverings should be used. **Cleaning Procedure:** Blood/body fluids must be thoroughly cleaned from surfaces and objects before disinfecting. **Contact Time:** Let stand for one minute. **Disposal of Infectious Materials:** Blood, body fluids and cleaning materials must be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal. This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semicritical medical devices prior to sterilization or high-level disinfection.

Kills bacteria, viruses & fungi in 1 MINUTE:

Bacteria:

- *Acinetobacter baumannii*
- *Enterobacter aerogenes*
- *Enterococcus faecium* (VRE)
- *Escherichia coli* ESBL
- *Klebsiella pneumoniae*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Staphylococcus aureus*
- *Staphylococcus aureus* (MRSA)
- *Streptococcus pyogenes*

Fungi:

- *Trichophyton mentagrophytes* (Athlete's Foot Fungus)

Viruses:

- †Avian Influenza A Virus
- †Canine Parvovirus
- †Feline Panleukopenia Virus
- †Hepatitis A Virus (HAV)
- †Hepatitis B Virus (HBV)
- †Hepatitis C Virus (HCV)
- †Herpes Simplex Virus (HSV-2)
- †HIV-1 (AIDS Virus)
- †Influenza A H1N1
- †Influenza A Virus
- †Norovirus
- †Poliovirus Type 1 (Mahoney)
- †Rhinovirus
- †Rotavirus

Kills Mycobacterium bovis - BCG, (TB) in 30 SECONDS

Kills Clostridium difficile spores in 5 MINUTES*

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION:

Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling. **FIRST AID: If in eyes:** Hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor or going for treatment. **PHYSICAL OR CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers or toilet bowl cleaners. This product not intended for use on fabrics.

STORAGE AND DISPOSAL:

Store this product in a cool, dry area away from direct sunlight and heat. Nonrefillable container. Do not reuse or refill this container. When discarding container, place in trash collection or take to recycling center. Do not contaminate water, food or feed by storage and disposal.

EPA Reg. No. 56392-7
EPA Est. Nos. 48181-IL-1
50757-WI-1
62801-MI-1
8251-OH-1
5813-GA-2

(Actual EPA Est. No. near batch code.)

Manufactured for:

©2010 Clorox Professional Products Company

1221 Broadway, Oakland, CA 94612

Visit us at www.cloroxprofessional.com

800-234-7700

® Trademark of Clorox

Professional Products Company



NI-15594

R0885-1.1