



8013.68



Kills
C. difficile^s
 spores in
 3 minutes

HEALTHCARE®

Bleach Germicidal Cleaner

Limpiador Germicida con Blanqueador

- Bactericidal, Virucidal[§], Fungicidal, Sporicidal*
- Ready-to-Use 1:10 Dilution of Bleach and Detergent
- Meets OSHA Bloodborne Pathogens Standard
- Aids in Reducing Cross-Contamination Between Treated Surfaces

KEEP OUT OF REACH OF CHILDREN
CAUTION: READ COMPLETE DIRECTIONS
 AND PRECAUTIONS ON REAR PANELS.

MANTÉNGASE FUERA DEL ALCANCE DE LOS NIÑOS
PRECAUCIÓN: LEA LAS INSTRUCCIONES Y PRECAUCIONES
 COMPLETAS EN LOS PANELES DE ATRÁS.
 Las instrucciones y la información de precaución
 se pueden encontrar en español en el panel posterior.

ACTIVE INGREDIENT:
 Sodium Hypochlorite..... 0.65%
 OTHER INGREDIENTS..... 99.35%
 TOTAL..... 100.00%
 Contains no phosphorous.

1.38 PT (22 FL OZ) 650 mL

NI-335012

Peel for
 Additional Directions

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Always refer to manufacturer's care instructions before using on equipment or devices.

To Refill: Remove: Remove trigger sprayer from empty bottle. **Pour:** Unscrew cap on refill and pour contents directly into empty bottle. **Use:** Replace trigger sprayer and use as you normally would. **To Clean, Disinfect and Deodorize Hard, Nonporous Surfaces:** Spray 6"-8" from surface until surface is thoroughly wet. Allow this product to remain wet for 1 minute. To kill *Clostridium difficile*^s spores: let stand for 3 minutes. Then wipe with a clean, damp cloth. For heavily soiled areas, a precleaning is required. Rinse with potable water for food-contact surfaces. **To Clean:** Follow directions above or apply with a clean cloth, mop, or sponge saturated with Clorox Healthcare® Bleach Germicidal Cleaner until surface is completely wet. (When applying with a saturated cloth or sponge, gloves should be worn.) No standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel.

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS.

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. **FIRST AID: If in Eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. If present, remove contact lenses after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor or going for treatment. **PHYSICAL OR CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners. This product not intended for use on fabrics.

QUESTIONS OR COMMENTS?

Call (800) 234-7700. For SDS information, please go to www.cloroxhealthcare.com.
 A list of this product's ingredients is available at www.IngredientsInside.com.

Mfd. for Clorox Professional Products Company,
 1221 Broadway, Oakland, CA 94612.

Clorox Healthcare is a registered trademark of
 Clorox Professional Products Company.
 © 2016 Made in USA of global ingredients.

EPA Reg. No. 56392-7. EPA Est. No. 8251-OH-1 (Q6);
 5813-GA-1 (A4), GA-2 (VG). Beginning of batch code
 indicates Est. No.



NI-33501

HINGED AREA, NO PRINT ALLOWED

STORAGE AND DISPOSAL:

Store this product in a cool, dry area, away from direct sunlight and heat.

Refillable container. Refill empty Clorox Healthcare® Bleach Germicidal Cleaner spray container only from Clorox Healthcare® Bleach Germicidal Cleaner refill bottle. Cleaning before refilling is the responsibility of the refiller. Take care to update expiration date on refillable bottle and replace bottle if it becomes worn or damaged. Do not reuse this container for any other purpose. Rinse container with water prior to disposal. Cleaning the container before final disposal is the responsibility of the person disposing of the container. When discarding container, place in trash collection or take to recycling center. Do not contaminate water, food or feed by storage and disposal.

Para Limpiar, Desinfectar y Desodorizar, superficies duras y sin poros: Rocíe de 6 a 8 pulgadas de la superficie hasta que esté bien mojada. Deje que este producto permanezca mojado durante 1 minuto. Para eliminar las esporas de *Clostridium difficile*⁵, deje reposar 3 minutos. Luego limpie con un paño limpio y húmedo. Para áreas muy sucias, se requiere una limpieza previa. Para superficies que entren en contacto con los alimentos, enjuague con agua potable. **Para Limpiar:** Siga las instrucciones de arriba o aplique con un paño, esponja o trapeador limpio y saturado con el Limpiador Germicida con Blanqueador Clorox Healthcare® hasta que la superficie esté completamente mojada. (Se deben usar guantes cuando se aplique con un paño o esponja saturada.) No requiere tiempo de reposo. Si se observan rayas, limpie con un paño o con una toalla de papel limpia y húmeda.

DECLARACIONES PREVENTIVAS: RIESGOS PARA LOS SERES HUMANOS Y ANIMALES DOMÉSTICOS.

PRECAUCIÓN: Causa irritación moderada a los ojos. Evite el contacto con los ojos o la ropa. Lávese bien con agua y jabón después de usar y antes de comer, beber, mascar goma, usar tabaco o usar el inodoro. **PRIMEROS AUXILIOS: Contacto con los Ojos:** Mantenga los ojos abiertos y enjuague con agua lentamente, con cuidado, durante 15–20 minutos. Si hay lentes de contacto, retírelos después de los primeros 5 minutos, luego continúe enjuagando los ojos. Llame a un médico o a un centro de control de envenenamiento para obtener más asesoría sobre el tratamiento. Cuando llame a un centro de control de envenenamiento o a un médico, o intente obtener tratamiento, tenga a la mano el envase o la etiqueta del producto. **RIESGOS FÍSICOS O QUÍMICOS:** No use este producto con amoníaco ni ácidos como vinagre, removedores de óxido o productos para limpiar el inodoro. Este producto no está hecho para usarse en telas.

HINGED AREA, NO PRINT ALLOWED

DIRECTIONS FOR USE continued: Special

Instructions for Cleaning Prior to Disinfection Against *Clostridium difficile* spores. Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks, and eye covering. **Cleaning**

Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Spray 6"–8" from surface until surface is completely wet. Let stand for 3 minutes. Wipe with a clean, damp cloth or paper towel or allow to air dry. **KILLS HIV, HBV AND HCV ON PRE-CLEANED, NONPOROUS SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV). **Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV and HCV on Surfaces/Objects Soiled with Blood/Body Fluids. Personal Protection:** When handling items soiled with blood or body fluids, use disposable gloves, gowns, masks, and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying this product. **Contact Time:** Let stand for one minute. For all other organisms, see directions for contact times. **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

PEEL

Multipurpose, broad spectrum cleaner disinfectant for use in healthcare and other areas where control of cross-contamination on surfaces is required. Meets OSHA Bloodborne Pathogens Standard.

Efficacy tests have demonstrated that this product is an effective bactericide, virucide[†] and fungicide in the presence of organic soil (5% blood serum).

ORGANISMS

Bacteria: 1 minute contact time

- *Acinetobacter baumannii*
- *Bordetella pertussis*
- *Campylobacter jejuni*
- Carbapenem-resistant *Klebsiella pneumoniae* (CRKP)
- *Enterobacter aerogenes*
- *Enterobacter cloacae* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- *Enterococcus hirae*
- ESBL (Extended Spectrum Beta Lactamase) producing *Escherichia coli*
- ESBL (Extended Spectrum Beta Lactamase) producing *Klebsiella pneumoniae*
- *Escherichia coli* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- *Klebsiella pneumoniae*
- *Klebsiella pneumoniae* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Multi-drug-resistant *Enterococcus faecium*
- *Proteus mirabilis*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Serratia marcescens*
- *Staphylococcus aureus*
- *Staphylococcus epidermidis*
- *Streptococcus pyogenes*
- *Stenotrophomonas maltophilia*
- *Streptococcus pneumoniae*
- Vancomycin-intermediate-resistant *Staphylococcus aureus* (VISA)
- Vancomycin-resistant *Staphylococcus aureus* (VRSA)
- Vancomycin-resistant *Enterococcus faecalis* (VRE)
- Vancomycin-resistant *Enterococcus faecium* (VRE)

Bacteria: 3 minute contact time

- * *Clostridium difficile*[‡] spores

Viruses[†]: 1 minute contact time

- *Adenovirus Type 2
- *Avian Influenza A (H3N2) Virus
- *Canine Parvovirus
- *Enterovirus Type D68
- *Feline Panleukopenia Virus
- *Hepatitis A Virus (HAV)
- *Hepatitis B Virus (HBV) (Duck Hep B Virus [DHBV] as surrogate)[†]
- *Hepatitis C Virus (HCV) (Bovine Viral Diarrhea Virus [BVDV] as surrogate)[†]
- *Herpes Simplex Virus Type 2 (HSV-2)
- *Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus)[†]
- *Influenza A H1N1 Virus
- *Influenza A Virus
- *Influenza B Virus
- *Measles Virus
- *MERS – Coronavirus
- *Norovirus (Feline Calicivirus as surrogate)
- *Poliovirus
- *Rhinovirus Type 37
- *Rotavirus
- *SARS-Associated Coronavirus

Fungi: 1 minute contact time

- *Candida albicans*
- *Trichophyton mentagrophytes* (Athlete's Foot Fungus)

Kills *Clostridium difficile*[‡] spores in 3 MINUTES in the presence of 3-part soil load.

[‡]Follow special instructions for cleaning prior to disinfection.
[†]Bloodborne pathogens

Compatible with a wide variety of surfaces.
 Check cloroxhealthcare.com/compatible.

SURFACES

May be used on hard, nonporous surfaces such as:

- Chrome
- Enamel
- Fiberglass (sealed) surfaces
- Formica[®]
- Glass
- Glazed ceramic tile
- Laminated surfaces
- Plastic^{††}
- Stainless steel
- Vinyl

Avoid contact with:

- Aluminum
- Brass
- Copper
- Corian[®]
- Fabric
- Galvanized steel
- Rubber
- Sealed granite
- Sealed marble
- Silver
- Wood

Medical hard, nonporous surfaces:

- Autoclaves^{††}
- Bed pans
- Bed railings
- Blood glucose monitors^{††}
- Exam tables
- High touch surfaces
- Mattress covers
- Medical equipment surfaces^{††}

Dental hard, nonporous surfaces:

- Counters
- Operatory hard surfaces

General hard, nonporous surfaces:

- Bathroom fixtures
- Countertops
- Doorknobs
- Exterior toilet surfaces
- Garbage cans
- Showers
- Sinks
- Tables

USE SITES

Medical:

- Dialysis clinics
- Emergency rooms
- Home healthcare
- Hospitals
- Intensive care units
- Isolation areas
- Nursing homes
- Outpatient clinics
- Operating rooms
- Patient restrooms
- Patient rooms

Dental:

- Dental clinics
- Labs

^{††}For any use site that includes electronics, refer to manufacturers' instructions before use.

^{‡‡}For unfamiliar surfaces, test in a small, inconspicuous area first. Consult medical equipment warranty information and equipment care and cleaning guides prior to selecting and using any surface disinfectant on plastic surfaces.

HINGED AREA, NO PRINT ALLOWED

HINGED AREA, NO PRINT ALLOWED



DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Always refer to manufacturer's care instructions before using on equipment or devices.

To Refill: Remove cap from empty bottle. **Pour:** Unscrew cap on refill and pour contents directly into empty bottle. **Use:** Replace cap and use as you normally would.

To Clean, Disinfect and Deodorize Hard, Nonporous Surfaces: Spray 6"-8" from surface until surface is thoroughly wet. Allow this product to remain wet for 1 minute. To kill *Clostridium difficile*[§] spores: let stand for 3 minutes. Then wipe with a clean, damp cloth. For heavily soiled areas, a precleaning is required. Rinse with potable water for food-contact surfaces. **To Clean:** Follow directions above or apply with a clean cloth, mop, or sponge saturated with Clorox Healthcare® Bleach Germicidal Cleaner until surface is completely wet. (When applying with a saturated cloth or sponge, gloves should be worn.) No standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel. **Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores.**

Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks, and eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Spray 6"-8" from surface until surface is completely wet. Let stand for 3 minutes. Wipe with a clean, damp cloth or paper towel or allow to air dry. **KILLS HIV, HBV, AND HCV ON PRE-CLEANED, NONPOROUS SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV). **Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV and HCV on Surfaces/Objects Soiled with Blood/Body Fluids.**

Personal Protection: When handling items soiled with blood or body fluids, use disposable gloves, gowns, masks, and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying this product. **Contact Time:** Let stand for 1 minute. For all other organisms, see directions for contact times. **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Multipurpose, broad spectrum cleaner disinfectant for use in healthcare and other areas where control of cross-contamination on surfaces is required. Meets OSHA Bloodborne Pathogens Standard.

Compatible with a wide variety of surfaces. Check cloroxhealthcare.com/compatible.

SURFACES: May be used on hard, nonporous surfaces such as chrome, enamel, fiberglass (sealed) surfaces, Formica® glass, glazed ceramic tile, laminated surfaces, plastic[‡], stainless steel, vinyl. **Avoid contact with:** aluminum, brass, copper, Corian®, fabric, galvanized steel, rubber, sealed granite, sealed marble, silver, wood. **Medical hard, nonporous surfaces:** autoclaves^{††}, bed pans, bed railings, blood glucose monitors^{††}, exam tables, high touch surfaces, mattress covers, medical equipment surfaces^{††}. **Dental hard, nonporous surfaces:** counters, operator hard surfaces. **General hard, nonporous surfaces:** bathroom fixtures, countertops, doorknobs, exterior toilet surfaces, garbage cans, showers, sinks, tables. **USE SITES: Medical** – dialysis clinics, emergency rooms, home healthcare, hospitals, intensive care units, isolation areas, nursing homes, outpatient clinics, operating rooms, patient restrooms, patient rooms. **Dental** – dental clinics, labs.

†† For any use site that includes electronics, refer to manufacturers' instructions before use.

‡ For unfamiliar surfaces, test in a small, inconspicuous area first. Consult medical equipment warranty information and equipment care and cleaning guides prior to selecting and using any surface disinfectant on plastic surfaces.

QUESTIONS OR COMMENTS? Call (800) 234-7700.

For SDS information, please go to www.cloroxhealthcare.com. A list of this product's ingredients is available at www.Ingredientsinside.com.

Mfd. for Clorox Professional Products Company, 1221 Broadway, Oakland, CA 94612.

Clorox Healthcare is a registered trademark of Clorox Professional Products Company.

© 2016 Made in USA of global ingredients.

EPA Reg. No. 56392-7. EPA Est. No. 8251-OH-1 (06); 5813-GA-1 (A4), GA-2 (V6).

Beginning of batch code indicates Est. No..



HEALTHCARE®

Bleach Germicidal Cleaner

Limpiador Germicida con Blanqueador

Kills *C. difficile*[§] spores in 3 minutes

- Bactericidal, Virucidal[‡], Fungicidal, Sporicidal*
- Ready-to-Use 1:10 Dilution of Bleach and Detergent

- Meets OSHA Bloodborne Pathogens Standard
- Aids in Reducing Cross-Contamination Between Treated Surfaces

KEEP OUT OF REACH OF CHILDREN

CAUTION: READ COMPLETE DIRECTIONS AND PRECAUTIONS ON REAR PANELS.

MANTÉNGASE FUERA DEL ALCANCE DE LOS NIÑOS

PRECAUCIÓN: LEA LAS INSTRUCCIONES Y PRECAUCIONES COMPLETAS EN LOS PANELES DE ATRÁS.

Las instrucciones y la información de precaución se pueden encontrar en español en el panel posterior.

ACTIVE INGREDIENT:
 Sodium Hypochlorite.....0.65%
 OTHER INGREDIENTS.....99.35%
 TOTAL.....100.00%

Contains no phosphorus.

1 QT (32 FL OZ) 946 mL

Efficacy tests have demonstrated that this product is an effective bactericide, virucide[‡] and fungicide in the presence of organic soil (5% blood serum).

ORGANISMS:

Bacteria: 1 minute contact time

- *Acinetobacter baumannii*
- *Bordetella pertussis*
- *Campylobacter jejuni*
- Carbapenem-resistant *Klebsiella pneumoniae* (CRKP)
- *Enterobacter aerogenes*
- *Enterobacter cloacae* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- *Enterococcus hirae*
- ESBL (Extended Spectrum Beta Lactamase) producing *Escherichia coli*
- ESBL (Extended Spectrum Beta Lactamase) producing *Klebsiella pneumoniae*
- *Escherichia coli* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- *Klebsiella pneumoniae*
- *Klebsiella pneumoniae* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Multi-drug-resistant *Staphylococcus aureus* (VISA)
- *Proteus mirabilis*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Serratia marcescens*

- *Staphylococcus aureus*
- *Staphylococcus epidermidis*
- *Streptococcus pyogenes*
- *Streptotrophomonas maltophilia*
- *Streptococcus pneumoniae*
- Vancomycin-intermediate-resistant *Staphylococcus aureus* (VISA)
- Vancomycin-resistant *Staphylococcus aureus* (VISA)
- Vancomycin-resistant *Staphylococcus aureus* (VISA)
- Vancomycin-resistant *Enterococcus faecalis* (VRE)
- Vancomycin-resistant *Enterococcus faecium* (VRE)

Bacteria: 3 minute contact time

- *Clostridium difficile*[§] spores

Viruses[‡]: 1 minute contact time

- Adenovirus Type 2
- Avian Influenza A (H3N2) Virus
- Canine Parvovirus
- Enterovirus Type D68
- Feline Panleukopenia Virus
- Hepatitis A Virus (HAV)
- Hepatitis B Virus (HBV) (Duck Hep B Virus [DHBV] as surrogate)[†]
- Hepatitis C Virus (HCV) (Bovine Viral Diarrhea Virus [BVDV] as surrogate)[†]
- Herpes Simplex Virus Type 2 (HSV-2)
- Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus)[†]
- Influenza A H1N1 Virus
- Influenza A Virus
- Influenza B Virus
- Measles Virus
- MERS – Coronavirus
- Norovirus (Feline Calicivirus as surrogate)
- Poliovirus
- Rhinovirus Type 37
- Rotavirus
- SARS-Associated Coronavirus

Fungi: 1 minute contact time

- *Candida albicans*
- *Trichophyton mentagrophytes* (Athlete's Foot Fungus)

Kills *Clostridium difficile*[§] spores in 3 MINUTES in the presence of three-part soil load.

§Follow special instructions for cleaning prior to disinfection. †Bloodborne pathogens

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS.

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. **FIRST AID: If in Eyes:** Hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. **PHYSICAL OR CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners. This product not intended for use on fabrics.

STORAGE AND DISPOSAL: Store this product in a cool, dry area, away from direct sunlight and heat. Refillable container. Refill empty Clorox Healthcare® Bleach Germicidal Cleaner spray container only from Clorox Healthcare® Bleach Germicidal Cleaner refill bottle. Cleaning before refilling is the responsibility of the refiller. Take care to update expiration date on refillable bottle and replace bottle if it becomes worn or damaged. Do not reuse this container for any other purpose. Rinse container with water prior to disposal. Cleaning the container before final disposal is the responsibility of the person disposing of the container. When discarding container, place in trash collection or take to recycling center. Do not contaminate water, food or feed by storage and disposal.





8013.68



Kills
C. difficile^s
 spores in
 3 minutes

HEALTHCARE®

Bleach Germicidal Cleaner

Limpiador Germicida con Blanqueador

- Bactericidal, Virucidal[‡], Fungicidal, Sporicidal*
- Ready-to-Use 1:10 Dilution of Bleach and Detergent
- Meets OSHA Bloodborne Pathogens Standard
- Aids in Reducing Cross-Contamination Between Treated Surfaces

KEEP OUT OF REACH OF CHILDREN
CAUTION: READ COMPLETE DIRECTIONS AND PRECAUTIONS ON REAR PANELS.

MANTÉNGASE FUERA DEL ALCANCE DE LOS NIÑOS
PRECAUCIÓN: LEA LAS INSTRUCCIONES Y PRECAUCIONES COMPLETAS EN LOS PANELES DE ATRAS.
 Las instrucciones y la información de precaución se pueden encontrar en español en el panel posterior.

ACTIVE INGREDIENT:
 Sodium Hypochlorite.....0.65%
 OTHER INGREDIENTS.....99.35%
 TOTAL.....100.00%
 Contains no phosphorous.

1 QT (32 FL. OZ) 946 mL

NI-33254

Peel for
 Additional Directions

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Always refer to manufacturer's care instructions before using on equipment or devices.
To Refill: Remove trigger sprayer from empty bottle. **Pour:** Unscrew cap on refill and pour contents directly into empty bottle. **Use:** Replace trigger sprayer and use as you normally would. **To Clean, Disinfect and Deodorize Hard, Nonporous Surfaces:** Spray 6"-8" from surface until surface is thoroughly wet. Allow this product to remain wet for 1 minute. To kill *Clostridium difficile*^s spores: let stand for 3 minutes. Then wipe with a clean, damp cloth. For heavily soiled areas, a precleaning is required. Rinse with potable water for food-contact surfaces. **To Clean:** Follow directions above or apply with a clean cloth, mop, or sponge saturated with Clorox Healthcare® Bleach Germicidal Cleaner until surface is completely wet. (When applying with a saturated cloth or sponge, gloves should be worn.) No standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel.

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS.

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. **FIRST AID: If in Eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. If present, remove contact lenses after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor or going for treatment. **PHYSICAL OR CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners. This product not intended for use on fabrics.

QUESTIONS OR COMMENTS?


Call (800) 234-7700. For SDS information, please go to www.cloroxhealthcare.com. A list of this product's ingredients is available at www.IngredientsInside.com.

Mfd. for Clorox Professional Products Company,
 1221 Broadway, Oakland, CA 94612.

Clorox Healthcare is a registered trademark of Clorox Professional Products Company.

© 2016 Made in USA of global ingredients.

EPA Reg. No. 56392-7, EPA Est. No. 8251-OH-1 (Q6);
 5813-GA-1 (A4), GA-2 (VG). Beginning of batch code indicates Est. No.



NI-33255

STORAGE AND DISPOSAL:

Store this product in a cool, dry area, away from direct sunlight and heat.

Refillable container. Refill empty Clorox Healthcare® Bleach Germicidal Cleaner spray container only from Clorox Healthcare® Bleach Germicidal Cleaner refill bottle. Cleaning before refilling is the responsibility of the refiller. Take care to update expiration date on refillable bottle and replace bottle if it becomes worn or damaged. Do not reuse this container for any other purpose. Rinse container with water prior to disposal. Cleaning the container before final disposal is the responsibility of the person disposing of the container. When discarding container, place in trash collection or take to recycling center. Do not contaminate water, food or feed by storage and disposal.

Para Limpiar, Desinfectar y Desodorizar, superficies duras y sin poros: Rocíe de 6 a 8 pulgadas de la superficie hasta que esté bien mojada. Deje que este producto permanezca mojado durante 1 minuto. Para eliminar las esporas de *Clostridium difficile*®: deje reposar 3 minutos. Luego limpie con un paño limpio y húmedo. Para áreas muy sucias, se requiere una limpieza previa. Para superficies que entren en contacto con los alimentos, enjuague con agua potable. **Para Limpiar:** Siga las instrucciones de arriba o aplique con un paño, esponja o trapeador limpio y saturado con el Limpiador Germicida con Blanqueador Clorox Healthcare® hasta que la superficie esté completamente mojada. (Se deben usar guantes cuando se aplique con un paño o esponja saturada.) No requiere tiempo de reposo. Si se observan rayas, limpie con un paño o con una toalla de papel limpia y húmeda.

DECLARACIONES PREVENTIVAS: RIESGOS PARA LOS SERES HUMANOS Y ANIMALES DOMÉSTICOS.

PRECAUCIÓN: Causa irritación moderada a los ojos. Evite el contacto con los ojos o la ropa. Lávese bien con agua y jabón después de usar y antes de comer, beber, mascar goma, usar tabaco o usar el inodoro. **PRIMEROS AUXILIOS: Contacto con los Ojos:** Mantenga los ojos abiertos y enjuague con agua lentamente, con cuidado, durante 15–20 minutos. Si hay lentes de contacto, retírelos después de los primeros 5 minutos, luego continúe enjuagando los ojos. Llame a un médico o a un centro de control de envenenamiento para obtener más asesoría sobre el tratamiento. Cuando llame a un centro de control de envenenamiento o a un médico, o intente obtener tratamiento, tenga a la mano el envase o la etiqueta del producto. **RIESGOS FÍSICOS O QUÍMICOS:** No use este producto con amoníaco ni ácidos como vinagre, removedores de óxido o productos para limpiar el inodoro. Este producto no está hecho para usarse en telas.

HINGED AREA, NO PRINT ALLOWED

DIRECTIONS FOR USE continued: Special Instructions for Cleaning Prior to Disinfection Against *Clostridium difficile* spores. Personal Protection:

Wear appropriate barrier protection such as gloves, gowns, masks, and eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Spray 6"–8" from surface until surface is completely wet. Let stand for 3 minutes. Wipe with a clean, damp cloth or paper towel or allow to air dry. **KILLS HIV, HBV AND HCV ON PRE-CLEANED, NONPOROUS SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV). **Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV and HCV on Surfaces/Objects Soiled with Blood/Body Fluids. Personal Protection:** When handling items soiled with blood or body fluids, use disposable gloves, gowns, masks, and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying this product. **Contact Time:** Let stand for one minute. For all other organisms, see directions for contact times. **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

PEEL

HINGED AREA, NO PRINT ALLOWED

HINGED AREA, NO PRINT ALLOWED

Multipurpose, broad spectrum cleaner disinfectant for use in healthcare and other areas where control of cross-contamination on surfaces is required. Meets OSHA Bloodborne Pathogens Standard.

Efficacy tests have demonstrated that this product is an effective bactericide, virucide* and fungicide in the presence of organic soil (5% blood serum).

ORGANISMS

Bacteria: 1 minute contact time

- | | | |
|---|---|--|
| <ul style="list-style-type: none"> • <i>Acinetobacter baumannii</i> • <i>Bordetella pertussis</i> • <i>Campylobacter jejuni</i> • Carbapenem-resistant <i>Klebsiella pneumoniae</i> (CRKP) • <i>Enterobacter aerogenes</i> • <i>Enterobacter cloacae</i> – New Delhi Metallo-Beta Lactamase-1 (NDM-1) • <i>Enterococcus hirae</i> • ESBL (Extended Spectrum Beta Lactamase) producing <i>Escherichia coli</i> • ESBL (Extended Spectrum Beta Lactamase) producing <i>Klebsiella pneumoniae</i> | <ul style="list-style-type: none"> • <i>Escherichia coli</i> – New Delhi Metallo-Beta Lactamase-1 (NDM-1) • <i>Klebsiella pneumoniae</i> • <i>Klebsiella pneumoniae</i> – New Delhi Metallo-Beta Lactamase-1 (NDM-1) • Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) • Multi-drug-resistant <i>Enterococcus faecium</i> • <i>Proteus mirabilis</i> • <i>Pseudomonas aeruginosa</i> • <i>Salmonella enterica</i> • <i>Serratia marcescens</i> | <ul style="list-style-type: none"> • <i>Staphylococcus aureus</i> • <i>Staphylococcus epidermidis</i> • <i>Streptococcus pyogenes</i> • <i>Stenotrophomonas maltophilia</i> • <i>Streptococcus pneumoniae</i> • Vancomycin-intermediate-resistant <i>Staphylococcus aureus</i> (VISA) • Vancomycin-resistant <i>Staphylococcus aureus</i> (VRSA) • Vancomycin-resistant <i>Enterococcus faecalis</i> (VRE) • Vancomycin-resistant <i>Enterococcus faecium</i> (VRE) |
|---|---|--|

Bacteria: 3 minute contact time

- * *Clostridium difficile* spores

Viruses*: 1 minute contact time

- | | | |
|---|--|---|
| <ul style="list-style-type: none"> • *Adenovirus Type 2 • *Avian Influenza A (H3N2) Virus • *Canine Parvovirus • *Enterovirus Type D68 • *Feline Panleukopenia Virus • *Hepatitis A Virus (HAV) • *Hepatitis B Virus (HBV) (Duck Hep B Virus [DHBV] as surrogate)† | <ul style="list-style-type: none"> • *Hepatitis C Virus (HCV) (Bovine Viral Diarrhea Virus [BVDV] as surrogate)† • *Herpes Simplex Virus Type 2 (HSV-2) • *Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus)† • *Influenza A H1N1 Virus • *Influenza A Virus | <ul style="list-style-type: none"> • *Influenza B Virus • *Measles Virus • *MERS – Coronavirus • *Norovirus (Feline Calicivirus as surrogate) • *Poliovirus • *Rhinovirus Type 37 • *Rotavirus • *SARS-Associated Coronavirus |
|---|--|---|

Fungi: 1 minute contact time

- *Candida albicans*
- *Trichophyton mentagrophytes* (Athlete's Foot Fungus)

Kills *Clostridium difficile spores in 3 MINUTES in the presence of 3-part soil load.**

- †Follow special instructions for cleaning prior to disinfection.
- ‡Bloodborne pathogens

Compatible with a wide variety of surfaces.
 Check cloroxhealthcare.com/compatible.

SURFACES

May be used on hard, nonporous surfaces such as:

- | | | | |
|--|--|--|---|
| <ul style="list-style-type: none"> • Chrome • Enamel | <ul style="list-style-type: none"> • Fiberglass (sealed) surfaces • Formica® | <ul style="list-style-type: none"> • Glass • Glazed ceramic tile • Laminated surfaces | <ul style="list-style-type: none"> • Plastic†† • Stainless steel • Vinyl |
|--|--|--|---|

Avoid contact with:

- | | | | |
|---|---|---|--|
| <ul style="list-style-type: none"> • Aluminum • Brass • Copper | <ul style="list-style-type: none"> • Corian® • Fabric • Galvanized steel | <ul style="list-style-type: none"> • Rubber • Sealed granite • Sealed marble | <ul style="list-style-type: none"> • Silver • Wood |
|---|---|---|--|

Medical hard, nonporous surfaces:

- | | | | |
|--|---|--|--|
| <ul style="list-style-type: none"> • Autoclaves†† • Bed pans • Bed railings | <ul style="list-style-type: none"> • Blood glucose monitors†† • Exam tables | <ul style="list-style-type: none"> • High touch surfaces • Mattress covers | <ul style="list-style-type: none"> • Medical equipment surfaces†† |
|--|---|--|--|

Dental hard, nonporous surfaces:

- Counters
- Operator hard surfaces

General hard, nonporous surfaces:

- | | | | |
|--|---|---|---|
| <ul style="list-style-type: none"> • Bathroom fixtures • Countertops | <ul style="list-style-type: none"> • Doorknobs • Exterior toilet surfaces | <ul style="list-style-type: none"> • Garbage cans • Showers | <ul style="list-style-type: none"> • Sinks • Tables |
|--|---|---|---|

USE SITES

Medical:

- | | | | |
|--|--|--|--|
| <ul style="list-style-type: none"> • Dialysis clinics • Emergency rooms • Home healthcare | <ul style="list-style-type: none"> • Hospitals • Intensive care units • Isolation areas | <ul style="list-style-type: none"> • Nursing homes • Outpatient clinics • Operating rooms | <ul style="list-style-type: none"> • Patient restrooms • Patient rooms |
|--|--|--|--|

Dental:

- Dental clinics
- Labs

††For any use site that includes electronics, refer to manufacturers' instructions before use.

‡‡ For unfamiliar surfaces, test in a small, inconspicuous area first. Consult medical equipment warranty information and equipment care and cleaning guides prior to selecting and using any surface disinfectant on plastic surfaces.

HINGED AREA, NO PRINT ALLOWED



Shown on page 2

Shown this page

Shown on page 3



8013.68



HEALTHCARE®

Bleach Germicidal Cleaner

Limpiador Germicida con Blanqueador

**Kills
C. difficile^S
 spores in
 3 minutes**

- Bactericidal, Virucidal[‡], Fungicidal, Sporicidal*
- Ready-to-Use 1:10 Dilution of Bleach and Detergent

- Meets OSHA Bloodborne Pathogens Standard
- Aids in Reducing Cross-Contamination Between Treated Surfaces

KEEP OUT OF REACH OF CHILDREN
CAUTION. READ COMPLETE DIRECTIONS AND PRECAUTIONS ON REAR PANELS.

MANTÉNGASE FUERA DEL ALCANCE DE LOS NIÑOS
PRECAUCIÓN. LEA LAS INSTRUCCIONES Y PRECAUCIONES COMPLETAS EN LOS PANELES DE ATRAS.
 Las instrucciones y la información de precaución se pueden encontrar en español en el panel posterior.

ACTIVE INGREDIENT:
 Sodium Hypochlorite.....0.65%
 OTHER INGREDIENTS.....99.35%
 TOTAL.....100.00%

Contains no phosphorous.

2 QT (64 FL OZ) 1.89 L

NI-33504

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Always refer to manufacturer's care instructions before using on equipment or devices.

To Refill Spray Bottles: Remove trigger sprayer from empty bottle. Unscrew cap on refill and pour contents directly into empty bottle. Replace trigger sprayer and use as you normally would. **To Clean, Disinfect and Deodorize Hard, Nonporous Surfaces:** Spray 6"-8" from surface until surface is thoroughly wet. Allow this product to remain wet for 1 minute. To kill *Clostridium difficile*⁹ spores: let stand for 3 minutes. Then wipe with a clean, damp cloth. For heavily soiled areas, a precleaning is required. Rinse with potable water for food-contact surfaces. **To Clean:** Follow directions above or apply with a clean cloth, mop, or sponge saturated with Clorox Healthcare[®] Bleach Germicidal Cleaner until surface is completely wet. (When applying with a saturated cloth or sponge, gloves should be worn.) No standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel. **Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores.**

Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks, and eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Materials Disposal: Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Spray 6"-8" from surface until surface is completely wet. Let stand for 3 minutes. Wipe with a clean, damp cloth or paper towel or allow to air dry. **KILLS HIV, HBV, AND HCV ON PRE-CLEANED, NONPOROUS SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV). **Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV and HCV on Surfaces/Objects Soiled with Blood/Body Fluids.** **Personal Protection:** When handling items soiled with blood or body fluids, use disposable gloves, gowns, masks, and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying this product. **Contact Time:** Let stand for 1 minute. For all other organisms, see directions for contact times. **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Multipurpose, broad spectrum cleaner disinfectant for use in healthcare and other areas where control of cross-contamination on surfaces is required. Meets OSHA Bloodborne Pathogens Standard.

QUESTIONS OR COMMENTS? Call (800) 234-7700.

For SDS information, please go to www.cloroxhealthcare.com. A list of this product's ingredients is available at www.IngredientsInside.com.

Mfd. for Clorox Professional Products Company, 1221 Broadway, Oakland, CA 94612.

Clorox Healthcare is a registered trademark of Clorox Professional Products Company.

© 2016 Made in USA of global ingredients.

EPA Reg. No. 56392-7. EPA Est. No. 8251-OH-1 (06); 5813-GA-1 (A4); GA-2 (VG). Beginning of batch code indicates Est. No.

Compatible with a wide variety of surfaces. Check cloroxhealthcare.com/compatible.

SURFACES: May be used on hard, nonporous surfaces such as chrome, enamel, fiberglass (sealed) surfaces, Formica®, glass, glazed ceramic tile, laminated surfaces, plastic†; stainless steel, vinyl. **Avoid contact with:** aluminum, brass, copper, Corian® fabric, galvanized steel, rubber, sealed granite, sealed marble, silver, wood. **Medical hard, nonporous surfaces:** autoclaves†, bed pans, bed railings, blood glucose monitors†, exam tables, high touch surfaces, mattress covers, medical equipment surfaces†. **Dental hard, nonporous surfaces:** counters, operatory hard surfaces. **General hard, nonporous surfaces:** bathroom fixtures, countertops, doorknobs, exterior toilet surfaces, garbage cans, showers, sinks, tables. **USE SITES: Medical** – dialysis clinics, emergency rooms, hospitals, home healthcare, intensive care units, isolation areas, nursing homes, outpatient clinics, operating rooms, patient restrooms, patient rooms. **Dental** – dental clinics, labs.

†† For any use site that includes electronics, refer to manufacturers' instructions before use.

‡‡ For unfamiliar surfaces, test in a small, inconspicuous area first. Consult medical equipment warranty information and equipment care and cleaning guides prior to selecting and using any surface disinfectant on plastic surfaces.

Efficacy tests have demonstrated that this product is an effective bactericide, virucide† and fungicide in the presence of organic soil (5% blood serum).

ORGANISMS:

Bacteria: 1 minute contact time

- *Acinetobacter baumannii*
- *Bordetella pertussis*
- *Campylobacter jejuni*
- Carbapenem-resistant *Klebsiella pneumoniae* (CRKP)
- *Enterobacter aerogenes*
- *Enterobacter cloacae* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- *Enterococcus hirae*

- ESBL (Extended Spectrum Beta Lactamase) producing *Escherichia coli*
- ESBL (Extended Spectrum Beta Lactamase) producing *Klebsiella pneumoniae*
- *Escherichia coli* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- *Klebsiella pneumoniae*

- *Klebsiella pneumoniae* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Multi-drug-resistant *Enterococcus faecium*
- *Proteus mirabilis*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*

- *Serratia marcescens*
- *Staphylococcus aureus*
- *Staphylococcus epidermidis*
- *Streptococcus pyogenes*
- *Stenotrophomonas maltophilia*
- *Streptococcus pneumoniae*
- Vancomycin-intermediate-resistant *Staphylococcus aureus* (VISA)

- Vancomycin-resistant *Staphylococcus aureus* (VRSA)
- Vancomycin-resistant *Enterococcus faecalis* (VRE)
- Vancomycin-resistant *Enterococcus faecium* (VRE)

Bacteria: 3 minute contact time

- *Clostridium difficile*† spores

Viruses†: 1 minute contact time

- Adenovirus Type 2
- Avian Influenza A (H3N2) Virus
- Canine Parvovirus
- Enterovirus Type D68
- Feline Panleukopenia Virus
- Hepatitis A Virus (HAV)

- Hepatitis B Virus (HBV) (Duck Hep B Virus [DHBV] as surrogate)†
- Hepatitis C Virus (HCV) (Bovine Viral Diarrhea Virus [BVDV] as surrogate)†
- Herpes Simplex Virus Type 2 (HSV-2)
- Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus)†

- Influenza A H1N1 Virus
- Influenza A Virus
- Influenza B Virus
- Measles Virus
- MERS – Coronavirus
- Norovirus (Feline Calicivirus as surrogate)

- Poliovirus
- Rhinovirus Type 37
- Rotavirus
- SARS-Associated Coronavirus

Fungi: 1 minute contact time

- *Candida albicans*
- *Trichophyton mentagrophytes* (Athlete's Foot Fungus)

Kills *Clostridium difficile*† spores in 3 MINUTES in the presence of three-part soil load.


† Follow special instructions for cleaning prior to disinfection.

† Bloodborne pathogens

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. **FIRST AID: If in Eyes:** Hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. **PHYSICAL OR CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners. This product not intended for use on fabrics.

STORAGE AND DISPOSAL: Store this product in a cool, dry area, away from direct sunlight and heat. Nonrefillable container. Do not reuse or refill this container. When discarding container, place in trash collection or take to recycling center. Do not contaminate water, food, or feed by storage and disposal.





STATE OF HAWAII
Department of Agriculture
ACCEPTED
LICENSE NO.

8013.68



HEALTHCARE®

**Kills
C. difficile^s
spores in
3 minutes**

Bleach Germicidal Cleaner

Limpiador Germicida con Blanqueador

- Bactericidal, Virucidal[‡], Fungicidal, Sporicidal*
- Ready-to-Use 1:10 Dilution of Bleach and Detergent
- Meets OSHA Bloodborne Pathogens Standard
- Aids in Reducing Cross-Contamination Between Treated Surfaces

KEEP OUT OF REACH OF CHILDREN

CAUTION: READ COMPLETE DIRECTIONS AND PRECAUTIONS ON REAR PANELS.

MANTÉNGASE FUERA DEL ALCANCE DE LOS NIÑOS

PRECAUCIÓN: LEA LAS INSTRUCCIONES Y PRECAUCIONES COMPLETAS EN LOS PANELES DE ATRÁS. Las instrucciones y la información de precaución se pueden encontrar en español en el panel posterior.

ACTIVE INGREDIENT:
Sodium Hypochlorite.....0.65%
OTHER INGREDIENTS.....99.35%
TOTAL.....100.00%

Contains no phosphorous.

1 GALLON (128 FL OZ) 3.78 L

NI-33505

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Always refer to manufacturer's care instructions before using on equipment or devices.

To Refill Spray Bottles: Remove trigger sprayer from empty bottle. Unscrew cap on refill and pour contents directly into empty bottle. Replace trigger sprayer and use as you normally would. **To Clean, Disinfect and Deodorize Hard, Nonporous Surfaces:** Spray 6"–8" from surface until surface is thoroughly wet. Allow this product to remain wet for 1 minute. To kill *Clostridium difficile*[§] spores: let stand for 3 minutes. Then wipe with a clean, damp cloth. For heavily soiled areas, a precleaning is required. Rinse with potable water for food-contact surfaces. **To Clean:** Follow directions above or apply with a clean cloth, mop, or sponge saturated with Clorox Healthcare® Bleach Germicidal Cleaner until surface is completely wet. (When applying with a saturated cloth or sponge, gloves should be worn.) No standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel. **Special Instructions for Cleaning Prior to Disinfection Against *Clostridium difficile* spores. Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks, and eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Spray 6"–8" from surface until surface is completely wet. Let stand for 3 minutes. Wipe with a clean, damp cloth or paper towel or allow to air dry. **KILLS HIV, HBV, AND HCV ON PRE-CLEANED, NONPOROUS SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV). **Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV and HCV on Surfaces/Objects Soiled with Blood/Body Fluids. Personal Protection:** When handling items soiled with blood or body fluids, use disposable gloves, gowns, masks, and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying this product. **Contact Time:** Let stand for 1 minute. For all other organisms, see directions for contact times. **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal. This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Compatible with a wide variety of surfaces. Check cloroxhealthcare.com/compatible.

SURFACES: May be used on hard, nonporous surfaces such as chrome, enamel, fiberglass (sealed) surfaces, Formica®, glass, glazed ceramic tile, laminated surfaces, plastic††, stainless steel, vinyl. **Avoid contact with:** aluminum, brass, copper, Corian®, fabric, galvanized steel, rubber, sealed granite, sealed marble, silver, wood. **Medical hard, nonporous surfaces:** autoclaves††, bed pans, bed railings, blood glucose monitors††, exam tables, high touch surfaces, mattress covers, medical equipment surfaces††.

Dental hard, nonporous surfaces: counters, operatory hard surfaces. **General hard, nonporous surfaces:** bathroom fixtures, countertops, doorknobs, exterior toilet surfaces, garbage cans, showers, sinks, tables. **USE SITES: Medical** – dialysis clinics, emergency rooms, home healthcare, hospitals, intensive care units, isolation areas, nursing homes, outpatient clinics, operating rooms, patient restrooms, patient rooms. **Dental** – dental clinics, labs.

†† For any use site that includes electronics, refer to manufacturers' instructions before use.

‡‡ For unfamiliar surfaces, test in a small, inconspicuous area first. Consult medical equipment warranty information and equipment care and cleaning guides prior to selecting and using any surface disinfectant on plastic surfaces.

EPA Reg. No. 56392-7. EPA Est. No. 8251-OH-1 (Q6); 5813-GA-1 (A4); GA-2 (V6).
Beginning of batch code indicates Est. No.

Multipurpose, broad spectrum cleaner disinfectant for use in healthcare and other areas where control of cross-contamination on surfaces is required. Meets OSHA Bloodborne Pathogens Standard.

Efficacy tests have demonstrated that this product is an effective bactericide, virucide[§] and fungicide in the presence of organic soil (5% blood serum).

ORGANISMS:

Bacteria: 1 minute contact time

- *Acinetobacter baumannii*
- *Bordetella pertussis*
- *Campylobacter jejuni*
- Carbapenem-resistant *Klebsiella pneumoniae* (CRKP)
- *Enterobacter aerogenes*
- *Enterobacter cloacae* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- *Enterococcus hirae*
- ESBL (Extended Spectrum Beta Lactamase) producing *Escherichia coli*
- ESBL (Extended Spectrum Beta Lactamase) producing *Klebsiella pneumoniae*
- *Escherichia coli* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- *Klebsiella pneumoniae*
- *Klebsiella pneumoniae* – New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Multi-drug-resistant *Enterococcus faecium*
- *Proteus mirabilis*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Serratia marcescens*
- *Staphylococcus aureus*
- *Staphylococcus epidermidis*
- *Streptococcus pyogenes*
- *Stenotrophomonas maltophilia*
- *Streptococcus pneumoniae*
- Vancomycin-intermediate-resistant *Staphylococcus aureus* (VISA)
- Vancomycin-resistant *Staphylococcus aureus* (VRSA)
- Vancomycin-resistant *Enterococcus faecalis* (VRE)
- Vancomycin-resistant *Enterococcus faecium* (VRE)

Bacteria: 3 minute contact time

- *Clostridium difficile*[§] spores

Viruses: 1 minute contact time

- Adenovirus Type 2
- Avian Influenza A (H3N2) Virus
- Canine Parvovirus
- Enterovirus Type D68
- Feline Panleukopenia Virus
- Hepatitis A Virus (HAV)
- Hepatitis B Virus (HBV) (Duck Hep B Virus [DHBV] as surrogate)[†]
- Hepatitis C Virus (HCV) (Bovine Viral Diarrhea Virus [BVDV] as surrogate)[†]
- Herpes Simplex Virus Type 2 (HSV-2)
- Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus)[†]
- Influenza A H1N1 Virus
- Influenza A Virus
- Influenza B Virus
- Measles Virus
- MERS – Coronavirus
- Norovirus (Feline Calicivirus as surrogate)
- Poliovirus
- Rhinovirus Type 37
- Rotavirus
- SARS-Associated Coronavirus

Fungi: 1 minute contact time

- *Candida albicans*
- *Trichophyton mentagrophytes* (Athlete's Foot Fungus)

Kills *Clostridium difficile*[§] spores in 3 MINUTES in the presence of three-part soil load.

[§]Follow special instructions for cleaning prior to disinfection. [†]Bloodborne pathogens

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS.

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. **FIRST AID: If in Eyes:** Hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. **PHYSICAL OR CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners. This product not intended for use on fabrics.

STORAGE AND DISPOSAL: Store this product in a cool, dry area, away from direct sunlight and heat. Nonrefillable container. Do not reuse or refill this container. When discarding container, place in trash collection or take to recycling center. Do not contaminate water, food, or feed by storage and disposal.

QUESTIONS OR COMMENTS?

Call (800) 234-7700. For SDS information, please go to www.cloroxhealthcare.com.

A list of this product's ingredients is available at www.IngredientsInside.com.

Mfd. for Clorox Professional Products Company, 1221 Broadway, Oakland, CA 94612.

Clorox Healthcare is a registered trademark of Clorox Professional Products Company.
© 2016 Made in USA of global ingredients.

