



Kills *C. difficile*^{§§} spores in 3 minutes

DISPATCH[®]

HOSPITAL CLEANER DISINFECTANT TOWELS WITH BLEACH



- ▶ Ready-to-Use Hospital Disinfectant
- ▶ Pre-moistened Towels with 1:10 Dilution of Bleach and Detergent
- ▶ Kills Bacteria, Viruses* and Fungi in 1 minute on Hard, Nonporous Surfaces*
- ▶ Kills *Mycobacterium bovis* (TB) in 2 minutes

Bacteria:

- *Acinetobacter baumannii**
- *Clostridium difficile* spores §§***
- *Enterobacter aerogenes**
- *Enterococcus faecalis*,
Vancomycin Resistant*
- Multi-Drug Resistant
*Enterococcus faecium**
- *Enterococcus hirae**
- ESBL (Extended Spectrum Beta
Lactamase) producing *Escherichia coli**
- *Klebsiella pneumoniae**
- Methicillin Resistant
*Staphylococcus aureus**
- *Pseudomonas aeruginosa**
- *Salmonella enterica**
- *Serratia marcescens**

Organisms

- *Staphylococcus aureus**
- *Streptococcus pyogenes**
- *Mycobacterium bovis* (TB) 20° C**

Fungi

- *Trichophyton mentagrophytes**

Viruses*

- † Adenovirus type 2*
- † Avian Influenza A Virus*, §
- † Canine parvovirus*
- † Feline parvovirus*
- † Hepatitis A Virus*
- †, †† Human Immunodeficiency Virus
type 1 (HIV-1)*
- †, †† Hepatitis B Virus (HBV)
(as duck HBV)*
- †, †† Hepatitis C Virus (HCV)
(as bovine Diarrhea Virus)*

- Herpes simplex virus type 2*
- † Influenza A Virus*
- † Influenza A Virus H1N1*
- † Influenza B Virus*
- † Norovirus (as Feline Calicivirus)*
- † Poliovirus*
- † Rhinovirus*
- † Rotavirus*

*1-minute contact time
 ** 2-minute contact time
 ***3-minute contact time
 † Kills TB in 2 minutes
 †† Bloodborne pathogens
 §Kills Avian Influenza A virus on
 precleaned hard, nonporous surfaces
 §§Follow special instructions for
 cleaning prior to disinfection

ACTIVE INGREDIENT:	
Sodium Hypochlorite.....	0.65%
OTHER INGREDIENTS:	99.35%
TOTAL:	100.00%
(yields available chlorine0.62%)	

KEEP OUT OF REACH OF CHILDREN

CAUTION:

READ COMPLETE DIRECTIONS AND PRECAUTIONS ON REAR PANELS.
 NOT TO BE USED AS BABY WIPES.
 KEEP CONTAINER TIGHTLY SEALED WHEN NOT IN USE.



NET CONTENTS: 150 – 6.75" x 8" TOWELS; 2 LB (32 OZ) 907 g



DISPATCH® HOSPITAL CLEANER DISINFECTANT
TOWELS WITH BLEACH

Pre-moistened disposable towels for use in healthcare and other areas where control of cross-contamination is required. Meets OSHA Bloodborne Pathogens Standard. Efficacy tests have demonstrated that this product is an effective bactericide, virucide†, tuberculocide and fungicide in the presence of organic soil (5% blood serum). Kills 99.9999% *C. difficile*^{SS} spores in 3 minutes, in the presence of three-part soil load.

USE SITES:

Medical: Blood glucose monitors, carts, central service areas, clinics, dialysis clinics, dialysis machines, emergency rooms (ERs), home healthcare settings, hospitals, intensive care units, isolation areas, laboratories, nursing homes, operating rooms, outpatient clinics, patient restrooms, patient rooms, pharmacies, physician offices and ultrasound transducers and probes. **Dental:** Dental offices.

MEDICAL SURFACES:

Autoclaves, bed railings, cabinets, chairs, counters, exam tables, external surfaces of medical equipment, gurneys, high-touch surfaces, IV poles, mattress covers, mammography equipment, stretchers, showers, sinks, and wheelchairs.

SURFACES: For use on hard, nonporous surfaces such as chrome, fiberglass sealed surfaces, Formica, glass, glazed ceramic, laminated surfaces, plastic††, stainless steel, vinyl.

Avoid contact with: aluminum, brass, copper, Corian, fabric, galvanized steel, natural marble, rubber, sealed granite, silver, wood.

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet. **FIRST AID:** If in eyes, hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes of rinsing, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. **PHYSICAL OR CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers or toilet bowl cleaners.

STORAGE AND DISPOSAL: Store this product in a cool, dry area away from direct sunlight and heat. When not in use, keep center cap of lid closed to prevent moisture loss. Nonrefillable container. Do not reuse or refill this container. Discard used towels and container in trash or take container to recycling center. Do not flush towels down toilet.

†† For unfamiliar surfaces, test in a small, inconspicuous area first. Consult medical equipment warranty information and equipment care and cleaning guides prior to selecting and using any surface disinfectant on plastic surfaces.

EPA Reg. No. 56392-8

EPA Est. No. 8251-OH-1 (Q6); 5813-GA-2 (VG)

Beginning of batch code indicates Est. No.

REORDER NO. 69150

See expiration date on package.

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleansing. For surfaces that may come in contact with food, a potable water rinse is required. Avoid contact with clothes, fabric, wood, rubber, and painted and paper surfaces. Do not use on silver, copper, or chipped enamel. **APPLY CAREFULLY** – This product contains bleach. Always refer to manufacturer's care instructions before using on equipment. Product expiration date on back of canister or pack. Discard wipe after use. To prevent cross-contamination, do not reuse wipe.

Dispenser. To Open: Remove lid. Discard inner seal from canister. **To Thread:** Find center of wipe roll, twist corner of first wipe into a point and thread through the hole in the canister lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out and tearing off. When not in use, keep center cap of lid closed to prevent moisture loss.

To Clean and Disinfect and Deodorize Hard, Nonporous Surfaces: Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for the contact time listed on label. Let air dry. Gross filth must be removed prior to disinfecting.

For Killing *Clostridium difficile* Spores^{§5}: Pre-clean hard, nonporous surfaces by removing gross filth. Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for 3 minutes. Let air dry.

To Clean: Follow directions above, but no standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel.

Special Instructions for Cleaning Prior to Disinfection Against *Clostridium difficile* Spores. Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks and eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Wipe surface with towel until completely wet. Let stand for 3 minutes. Wipe dry or allow to air dry. Gross soil must be removed prior to disinfecting. Gloves must be worn.

Ultrasound Transducers and Probes: For use on external surfaces of ultrasound transducers and probes only.

To Disinfect: Pre-clean gel from external surface of ultrasound transducer or probe to be disinfected with a paper towel or clean cloth. Then thoroughly wet the probe with a Dispatch® Hospital Cleaner Disinfectant Towels with Bleach towel. Let stand for one minute. To kill TB, let stand for two minutes at room temperature 20° C.

This product kills HIV-1, HBV and HCV on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Human Hepatitis B Virus (HBV) and Human Hepatitis C Virus (HCV).

§5Special instructions for cleaning and decontamination against HIV-1 (AIDS Virus), HBV and HCV on surfaces/objects soiled with blood/body fluids:

Personal Protection: When handling items soiled with blood or body fluids, personal protective equipment such as disposable gloves, gowns, masks and eye coverings must be used. **Cleaning Procedure:** Blood/body fluids must be thoroughly cleaned from surfaces and objects before disinfecting. **Contact Time:** Allow surfaces to remain wet for one minute, let air dry. For all other organisms, see directions for contact times. **Disposal of Infectious Materials:** Blood, body fluids, and cleaning materials must be autoclaved and disposed according to federal, state, and local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semicritical medical devices prior to sterilization or high-level disinfection.

QUESTIONS OR COMMENTS? Call (800) 234-7700.
For SDS information, go to www.cloroxhealthcare.com.
For more product ingredient information,
visit www.IngredientsInside.com.
Mfd. for Clorox Professional Products Company,
1221 Broadway, Oakland, CA 94612.
© 2016 Dispatch is a registered trademark of
Clorox Professional Products Company.
Made in USA of global ingredients.



NI-34661



STATE OF HAWAII
 Department of Agriculture

ACCEPTED

LICENSE NO.

8013.82

Kills *C. difficile*^{§§} spores in 3 minutes

DISPATCH[®]



HOSPITAL CLEANER DISINFECTANT TOWELS WITH BLEACH

- ▶ Ready-to-Use Hospital Disinfectant
- ▶ Pre-Moistened Towels with 1:10 Dilution of Bleach and Detergent
- ▶ Kills Bacteria, Viruses[†] and Fungi in 1 Minute on Hard, Nonporous Surfaces[†]
- ▶ Kills *Mycobacterium bovis* (TB) in 2 Minutes

Organisms

Bacteria

- *Acinetobacter baumannii*^{*}
- *Clostridium difficile* spores^{§§,***}
- *Enterobacter aerogenes*^{*}
- *Enterococcus faecalis*, Vancomycin Resistant^{*}
- Multi-Drug Resistant *Enterococcus faecium*^{*}
- *Enterococcus hirae*^{*}
- ESBL (Extended Spectrum Beta Lactamase) Producing *Escherichia coli*^{*}
- *Klebsiella pneumoniae*^{*}
- Methicillin Resistant *Staphylococcus aureus*^{*}

- *Pseudomonas aeruginosa*^{*}
- *Salmonella enterica*^{*}
- *Serratia marcescens*^{*}
- *Staphylococcus aureus*^{*}
- *Streptococcus pyogenes*^{*}
- *Mycobacterium bovis* (TB) 20°C^{**}

Fungi

- *Trichophyton mentagrophytes*^{*}

Viruses[†]

- †Adenovirus Type 2^{*}
- †Avian Influenza A Virus^{*,§}
- †Canine Parvovirus^{*}
- †Feline Parvovirus^{*}
- †Hepatitis A Virus^{*}

- †,##Human Immunodeficiency Virus Type 1 (HIV-1)^{*}
- †,##Hepatitis B Virus (HBV) (as Duck HBV)^{*}
- †,##Hepatitis C Virus (HCV) (as Bovine Diarrhea Virus)^{*}
- †Herpes Simplex Virus Type 2^{*}
- †Influenza A Virus^{*}
- †Influenza A Virus H1N1^{*}
- †Influenza B Virus^{*}
- †Norovirus (as Feline Calicivirus)^{*}
- †Poliovirus^{*}
- †Rhinovirus^{*}
- †Rotavirus^{*}

* 1 minute contact time
 ** 2 minute contact time
 *** 3 minute contact time
 † Kills TB in 2 minutes

‡ Bloodborne pathogens
 § Kills Avian Influenza A Virus on pre-cleaned hard, nonporous surfaces
 §§ Follow special instructions for cleaning prior to disinfection

ACTIVE INGREDIENT:
 Sodium Hypochlorite.....0.65%
 OTHER INGREDIENTS:99.35%
 TOTAL:100.00%
 (Yields available chlorine...0.62%)

KEEP OUT OF REACH OF CHILDREN
CAUTION:
 READ COMPLETE DIRECTIONS AND PRECAUTIONS ON REAR PANELS
 NOT TO BE USED AS BABY WIPES
 KEEP PACKAGE TIGHTLY SEALED WHEN NOT IN USE



NET CONTENTS: 60 – 9" x 10" TOWELS; 1 LB 7 OZ

Shown this page

Shown on page 2

Shown on page 2



DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleansing. For surfaces that may come in contact with food, a potable water rinse is required. Expiration date on back of canister or pack. Discard wipe after use. To prevent cross-contamination, do not reuse wipe.

Dispenser: Open package. Pull out wipes as needed and reseal. When not in use, keep sealed to prevent moisture loss.

To Clean and Disinfect and Deodorize Hard, Nonporous Surfaces: Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for the contact time listed on label. Let air dry. Gross filth must be removed prior to disinfecting.

For Killing *Clostridium difficile* Spores^{††}: Pre-clean hard, nonporous surfaces by removing gross filth. Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for 3 minutes. Let air dry.

To Clean: Follow directions above, but no standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel.

Special Instructions for Cleaning Prior to Disinfection Against *Clostridium difficile* Spores. Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks and eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Materials Disposal: Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Wipe surface with towel until completely wet. Let stand for 3 minutes. Wipe dry or allow to air dry. Gross soil must be removed prior to disinfecting. Gloves must be worn.

Ultrasound transducers and probes: For use on external surfaces of ultrasound transducers and probes only. **To Disinfect:** Pre-clean gel from external surface of ultrasound transducer or probe to be disinfected with a paper towel or clean cloth. Then thoroughly wet the probe with a Dispatch[®] Hospital Cleaner Disinfectant Towels with Bleach towel. Let stand for one minute. To kill TB, let stand for two minutes at room temperature 20° C.

This product kills HIV-1, HBV and HCV on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Human Hepatitis B Virus (HBV) and Human Hepatitis C Virus (HCV).

††Special instructions for cleaning and decontamination against HIV-1 (AIDS Virus), HBV and HCV on surfaces/objects soiled with blood/body fluids: Personal Protection: When handling items soiled with blood or body fluids, personal protective equipment such as disposable gloves, gowns, masks and eye coverings must be used. **Cleaning Procedure:** Blood/body fluids must be thoroughly cleaned from surfaces and objects before disinfecting. **Contact Time:** Allow surfaces to remain wet for one minute, let air dry. For all other organisms, see directions for contact times. **Disposal of Infectious Materials:** Blood, body fluids, and cleaning materials must be autoclaved and disposed according to federal, state, and local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Pre-moistened disposable towels for use in healthcare and other areas where control of cross-contamination is required. Meets OSHA Bloodborne Pathogens Standard. Efficacy tests have demonstrated that this product is an effective bactericide, virucide[†], tuberculocide, and fungicide in the presence of organic soil (5% blood serum). Kills 99.9999% *C. difficile*^{§§} spores in 3 minutes, in the presence of three-part soil load.

USE SITES: Medical: Blood glucose monitors, carts, central service areas, clinics, dialysis clinics, dialysis machines, emergency rooms (ERs), home healthcare settings, hospitals, intensive care units, isolation areas, laboratories, nursing homes, operating rooms, outpatient clinics, patient restrooms, patient rooms, pharmacies, physician offices and ultrasound transducers and probes. **Dental:** Dental offices. **Medical Surfaces:** autoclaves, bed railings, cabinets, chairs, counters, exam tables, external surfaces of medical equipment, gurneys, high touch surfaces, IV poles, mattress covers, mammography equipment, stretchers, showers, sinks, and wheelchairs. **Surfaces:** For use on hard, nonporous surfaces such as chrome, fiberglass sealed surfaces, Formica, glass, glazed ceramic, laminated surfaces, plastic^{††}, stainless steel and vinyl. **Avoid Contact With:** aluminum, brass, copper, Corian, fabric, galvanized steel, natural marble, rubber, sealed granite, silver and wood. ††For unfamiliar surfaces, test in a small, inconspicuous area first. Consult medical equipment warranty information and equipment care and cleaning guides prior to selecting and using any surface disinfectant on plastic surfaces.

PRECAUTIONARY STATEMENTS:**HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet. **FIRST AID:** If in eyes, hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes of rinsing, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

PHYSICAL OR CHEMICAL HAZARDS: Do not use this product with ammonia or acids such as vinegar, rust removers or toilet bowl cleaners.

STORAGE AND DISPOSAL: Store this product in a cool, dry area away from direct sunlight and heat. Nonrefillable container. Do not reuse or refill this container. Discard used towel and empty package in trash. Do not flush towels down toilet.

**Manufactured for:**

**Clorox Professional Products Company,
1221 Broadway, Oakland, CA 94612.**

QUESTIONS OR COMMENTS? Call (800) 234-7700.

For SDS information, go to www.cloroxhealthcare.com.

For more product ingredient information, visit www.IngredientsInside.com.

Made in USA of global ingredients.

© 2016 Dispatch is a registered trademark of Clorox Professional Products Company.

EPA Reg. No. 56392-8

EPA Est. No. 8251-OH-1 (Q6); 5813-GA-2 (VG)

Beginning of batch code indicates Est. No.

REORDER NO. 69260


See expiration date on package.

NI-34664

**KEEP OUT OF REACH
OF CHILDREN
CAUTION:**

READ COMPLETE
DIRECTIONS AND
PRECAUTIONS ON
REAR PANELS

NOT TO BE USED
AS BABY WIPES
KEEP PACKAGE
TIGHTLY SEALED
WHEN NOT IN USE



DIRECTIONS FOR USE: It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleansing. For surfaces that may come in contact with food, a potable water rinse is required. Product expiration date on back of canister or pack. To prevent cross-contamination do not reuse wipe.

Dispenser: Open package. Pull out wipes as needed and reseal. When not in use, keep sealed to prevent moisture loss.

To Clean and Disinfect and Deodorize Hard, Nonporous Surfaces: Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for the contact time listed on label. Let air dry. Gross filth must be removed prior to disinfecting.

For Killing *Clostridium difficile* Spores⁶⁶: Pre-clean hard, nonporous surfaces by removing gross filth. Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for 3 minutes. Let air dry.

To Clean: Follow directions above, but no standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel.

**NET CONTENTS:
40 – 8" x 10" TOWELS 13.4 OZ**

Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores.

Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks and eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Wipe surface with towel until completely wet. Let stand for 3 minutes. Wipe dry or allow to air dry. Gross soil must be removed prior to disinfecting. Gloves must be worn.

Ultrasound transducers and probes: For use on external surfaces of ultrasound transducers and probes only. **To Disinfect:** Pre-clean gel from external surface of ultrasound transducer or probe to be disinfected with a paper towel or clean cloth. Then thoroughly wet the probe with a Dispatch® Hospital Cleaner Disinfectant Towels with Bleach towel. Let stand for one minute. To kill TB, let stand for two minutes at room temperature 20°C.

This product kills HIV-1, HBV and HCV on precleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Human Hepatitis B Virus (HBV) and Human Hepatitis C Virus (HCV). **Special instructions for cleaning and decontamination against HIV-1 (AIDS Virus), HBV and HCV on surfaces/objects soiled with blood/body fluids:** **Personal Protection:** When handling items soiled with blood or body fluids, personal protective equipment such as disposable gloves, gowns, masks and eye coverings must be used. **Cleaning Procedure:** Blood/body fluids must be thoroughly cleaned from surfaces and objects before disinfecting. **Contact Time:** Allow surfaces to remain wet for one minute, let air dry. For all other organisms, see directions for contact times. **Disposal of Infectious Materials:** Blood, body fluids, and cleaning materials must be autoclaved and disposed according to federal, state, and local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

STORAGE AND DISPOSAL:

Store this product in a cool, dry area away from direct sunlight and heat. Nonrefillable container. Do not reuse or refill this container. Discard used towel and empty package in trash. Do not flush towels down toilet.



Manufactured for:
Clorox Professional Products Company,
 1221 Broadway, Oakland, CA 94612.
 QUESTIONS OR COMMENTS? Call (800) 234-7700.
 For SDS information, go to www.cloroxhealthcare.com.
 For more product ingredient information, visit www.ingredientsinside.com

Made in USA of global ingredients.

© 2016 Dispatch is a registered trademark of Clorox Professional Products Company.

REORDER NO. 69240

See expiration date on package.

NI-3463

Pre-moistened disposable towels for use in healthcare and other areas where control of cross-contamination is required. Meets OSHA Bloodborne Pathogens Standard. Efficacy tests have demonstrated that this product is an effective bactericide, virucide[‡], tuberculocide and fungicide in the presence of organic soil (5% blood serum). Kills 99.9999% *C. difficile*^{§§} spores in 3 minutes, in the presence of three-part soil load.

Organisms

- | | |
|--|---|
| Bacteria: | Fungi: |
| ■ <i>Acinetobacter baumannii</i> [*] | ■ <i>Trichophyton mentagrophytes</i> [*] |
| ■ <i>Clostridium difficile</i> spores ^{§§, ***} | Viruses: |
| ■ <i>Enterobacter aerogenes</i> [*] | ■ ‡Adenovirus Type 2 [*] |
| ■ <i>Enterococcus faecalis</i> ,
Vancomycin Resistant [*] | ■ ‡Avian Influenza A Virus [§] |
| ■ Multi-Drug Resistant
<i>Enterococcus faecium</i> [*] | ■ ‡Canine Parvovirus [*] |
| ■ <i>Enterococcus hirae</i> [*] | ■ ‡Feline Parvovirus [*] |
| ■ ESBL (Extended Spectrum
Beta Lactamase) Producing
<i>Escherichia coli</i> [*] | ■ ‡Hepatitis A Virus [*] |
| ■ <i>Klebsiella pneumoniae</i> [*] | ■ ‡, ‡, ‡Human Immunodeficiency
Virus Type 1 (HIV-1) [*] |
| ■ Methicillin resistant
<i>Staphylococcus aureus</i> [*] | ■ ‡, ‡, ‡Hepatitis B Virus (HBV)
(as Duck HBV) [*] |
| ■ <i>Pseudomonas aeruginosa</i> [*] | ■ ‡, ‡, ‡Hepatitis C Virus (HCV)
(as Bovine Diarrhea Virus) [*] |
| ■ <i>Salmonella enterica</i> [*] | ■ ‡Herpes Simplex Virus Type 2 [*] |
| ■ <i>Serratia marcescens</i> [*] | ■ ‡Influenza A Virus [*] |
| ■ <i>Staphylococcus aureus</i> [*] | ■ ‡Influenza A Virus H1N1 [*] |
| ■ <i>Streptococcus pyogenes</i> [*] | ■ ‡Influenza B Virus [*] |
| ■ <i>Mycobacterium bovis</i> (TB) 20°C ^{**} | ■ ‡Norovirus (as Feline Calicivirus) [*] |
| | ■ ‡Poliovirus [*] |
| | ■ ‡Rhinovirus [*] |
| | ■ ‡Rotavirus [*] |

*1-minute contact time
 **2-minute contact time
 ***3-minute contact time
 †Kills TB in 2 minutes
 ‡Bloodborne pathogens
 §Kills Avian Influenza A virus on precleaned hard, nonporous surfaces
 §§Follow special instructions for cleaning prior to disinfection

USE SITES:

Medical: Blood glucose monitors, carts, central service areas, clinics, dialysis clinics, dialysis machines, emergency rooms (ERs), home healthcare settings, hospitals, intensive care units, isolation areas, laboratories, nursing homes, operating rooms, outpatient clinics, patient restrooms, patient rooms, pharmacies, physician offices, and ultrasound transducers and probes. **Dental:** Dental offices. **Medical Surfaces:** autoclaves, bed railings, cabinets, chairs, counters, exam tables, external surfaces of medical equipment, gurneys, high touch surfaces, IV poles, mattress covers, mammography equipment, stretchers, showers, sinks and wheelchairs. **Surfaces:** For use on hard nonporous surfaces such as chrome, fiberglass sealed surfaces, Formica, glass, glazed ceramic, laminated surfaces, plastic,^{††} stainless steel and vinyl. **Avoid contact with:** aluminum, brass, copper, Corian, fabric, galvanized steel, natural marble, rubber, sealed granite, silver and wood. ††For unfamiliar surfaces, test in a small, inconspicuous area first. Consult medical equipment warranty information and equipment care and cleaning guides prior to selecting and using any surface disinfectant on plastic surfaces.

PRECAUTIONARY STATEMENTS:

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. **FIRST AID:** If in eyes, hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes of rinsing, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

PHYSICAL OR CHEMICAL HAZARDS: Do not use this product with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners.

EPA Reg. No. 56392-8
 EPA Est. No. 8251-OH-1 (Q6); 5813-GA-2 (VG).
 Beginning of batch code indicates Est. No.



STATE OF HAWAII
Department of Agriculture

ACCEPTED

LICENSE NO.

8013.82

Single Wrap

Kills *C. difficile*^{SS} Spores in 3 Minutes

DISPATCH[®]

HOSPITAL CLEANER DISINFECTANT TOWELS WITH BLEACH

- ▶ Ready-to-Use Hospital Disinfectant
- ▶ Pre-Moistened Towels with 1:10 Dilution of Bleach and Detergent
- ▶ Kills Bacteria, Viruses[†] and Fungi in 1 minute on Hard, Nonporous Surfaces[‡]
- ▶ Kills *Mycobacterium bovis* (TB) in 2 minutes

NI-34660

†Kills TB in 2 minutes

KEEP OUT OF REACH OF CHILDREN

CAUTION:

READ COMPLETE DIRECTIONS AND PRECAUTIONS ON DISPENSER CARTON
NOT TO BE USED AS BABY WIPES

ACTIVE INGREDIENT:	
Sodium Hypochlorite.....	0.65%
OTHER INGREDIENTS.....	99.35%
TOTAL.....	100.00%

(yields available chlorine...0.62%)

1 - 7" x 8" TOWEL; 0.22 OZ

REORDER NO. 689101
See Expiration Date

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Dispenser: Tear open packet. Unfold. Discard wipe after use.

To Clean and Disinfect and Deodorize Hard, Nonporous Surfaces: Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for the contact time listed on label. Let air dry. Gross filth must be removed prior to disinfecting.

For Killing *Clostridium difficile* Spores[†]: Pre-clean hard, nonporous surfaces by removing gross filth. Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for 3 minutes. Let air dry.

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS
CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. **FIRST AID:** If in eyes, hold eye open and rinse slowly and gently with water for 15-20 minutes. If present, remove contact lenses after the first 5 minutes of rinsing, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. **PHYSICAL OR CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners.

STORAGE AND DISPOSAL: Store this product in a cool, dry area away from direct sunlight and heat. Nonrefillable container. Do not reuse or refill this container. Discard used towel and empty package in trash. Do not flush towels down toilet.

Manufactured for Clorox Professional Products

Company, 1221 Broadway, Oakland, CA 94612.

QUESTIONS OR COMMENTS? Call (800) 234-7700. For SDS information, go to www.cloroxhealthcare.com. For more product ingredient information, visit www.ingredientinside.com.

© 2016 Dispatch is a registered trademark of Clorox Professional Products Company.

EPA Reg. No. 56392-8

EPA Est. No. 8251-OH-1 (Q6);

5813-GA-2 (VG)


Beginning of batch code indicates

Est. No.

Shown on page 2



Shown on page 3



Shown on page 4



Remove lid to create dispenser.

Kills *C. difficile*^{SS} Spores in 3 Minutes

DISPATCH[®]



HOSPITAL CLEANER DISINFECTANT TOWELS WITH BLEACH

- ▶ Ready-to-Use Hospital Disinfectant
- ▶ Pre-Moistened Towels with 1:10 Dilution of Bleach and Detergent
- ▶ Kills Bacteria, Viruses[†] and Fungi in 1 Minute on Hard, Nonporous Surfaces[†]
- ▶ Kills *Mycobacterium bovis* (TB) in 2 Minutes

ACTIVE INGREDIENT:
 Sodium Hypochlorite 0.65%
 OTHER INGREDIENTS:.....99.35%
 TOTAL:.....100.00%
 (yields available chlorine...0.62%)

KEEP OUT OF REACH OF CHILDREN
CAUTION:
 READ COMPLETE DIRECTIONS AND PRECAUTIONS
 ON SIDE PANELS
 NOT TO BE USED AS BABY WIPES



CONTENTS: 50 – 7" x 8" TOWELS (0.22 OZ EACH) TOTAL 11 OZ

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleansing. For surfaces that may come in contact with food, a potable water rinse is required. Expiration date on back of canister or pack. Discard wipe after use. To prevent cross-contamination, do not reuse wipe.

Dispenser: Tear open packet. Unfold.

To Clean and Disinfect and Deodorize Hard, Nonporous Surfaces: Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for the contact time listed on label. Let air dry. Gross filth must be removed prior to disinfecting. **For Killing *Clostridium difficile* Spores⁹⁵:** Pre-clean hard, nonporous surfaces by removing gross filth. Wipe surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for 3 minutes. Let air dry. **To Clean:** Follow directions above, but no standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel.

Special Instructions for Cleaning Prior to Disinfection Against *Clostridium difficile* Spores.

Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks and eye covering.

Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Wipe surface with towel until completely wet. Let stand for 3 minutes. Wipe dry or allow to air dry. Gross soil must be removed prior to disinfecting. Gloves must be worn.

Ultrasound transducers and probes: For use on external surfaces of ultrasound transducers and probes only. **To Disinfect:** Pre-clean gel from external surface of ultrasound transducer or probe to be disinfected with a paper towel or clean cloth. Then thoroughly wet the probe with a Dispatch[®] Hospital Cleaner Disinfectant Towels with Bleach towel. Let stand for one minute. To kill TB, let stand for two minutes at room temperature 20°C.

This product kills HIV-1, HBV and HCV on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Human Hepatitis B Virus (HBV) and Human Hepatitis C Virus (HCV).

##Special instructions for cleaning and decontamination against HIV-1 (AIDS Virus), HBV and HCV on surfaces/objects soiled with blood/body fluids: Personal Protection: When handling items soiled with blood or body fluids, personal protective equipment such as disposable gloves, gowns, masks and eye coverings must be used. **Cleaning Procedure:** Blood/body fluids must be thoroughly cleaned from surfaces and objects before disinfecting. **Contact Time:** Allow surfaces to remain wet for one minute, let air dry. For all other organisms, see directions for contact times. **Disposal of Infectious Materials:** Blood, body fluids, and cleaning materials must be autoclaved and disposed according to federal, state, and local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Remove flaps to create dispenser.

PRECAUTIONARY STATEMENTS:**HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. **FIRST AID:** If in eyes, hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes of rinsing, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. **PHYSICAL OR CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers or toilet bowl cleaners.

STORAGE AND DISPOSAL: Store this product in a cool, dry area away from direct sunlight and heat. Nonrefillable container. Do not reuse or refill this container. Discard used towel and empty package in trash. Do not flush towels down toilet.

REORDER NO. 69101

See expiration date on package.

**Manufactured for:**

**Clorox Professional Products Company,
1221 Broadway, Oakland, CA 94612.**

QUESTIONS OR COMMENTS? Call (800) 234-7700.
or SDS information, go to www.cloroxhealthcare.com. For more
product ingredient information, visit www.IngredientsInside.com.

Made in USA of global ingredients.

© 2016 Dispatch is a registered trademark of
Clorox Professional Products Company.

EPA Reg. No. 56392-8

EPA Est. No. 8251-OH-1 (Q6); 5813-GA-2 (VG)

Beginning of batch code indicates Est. No.

Kills *C. difficile*^{SS} Spores in 3 Minutes

DISPATCH[®]



HOSPITAL CLEANER DISINFECTANT
 TOWELS WITH BLEACH

ACTIVE INGREDIENT:
 Sodium Hypochlorite..... 0.65%
 OTHER INGREDIENTS:..... 99.35%
 TOTAL:..... 100.00%
 (yields available chlorine...0.62%)

KEEP OUT OF REACH OF CHILDREN

CAUTION:

READ COMPLETE DIRECTIONS
 AND PRECAUTIONS ON SIDE PANELS
 NOT TO BE USED AS BABY WIPES

**CONTENTS: 50 – 7" x 8" TOWELS
 (0.22 OZ EACH) TOTAL 11 OZ**



Remove flaps to create dispenser.

Pre-moistened disposable towels for use in healthcare and other areas where control of cross-contamination is required. Meets OSHA Bloodborne Pathogens Standard. Efficacy tests have demonstrated that this product is an effective bactericide, virucide[†], tuberculocide and fungicide in the presence of organic soil (5% blood serum). Kills 99.9999% of *C. difficile*^{SS} spores in 3 minutes, in the presence of three-part soil load.

USE SITES: Medical: Blood glucose monitors, carts, central service areas, clinics, dialysis clinics, dialysis machines, emergency rooms (ERs), home healthcare settings, hospitals, intensive care units, isolation areas, laboratories, nursing homes, operating rooms, outpatient clinics, patient restrooms, patient rooms, pharmacies, physician offices and ultrasound transducers and probes. **Dental:** dental offices.

Medical Surfaces: autoclaves, bed railings, cabinets, chairs, counters, exam tables, external surfaces of medical equipment, gurneys, high touch surfaces, IV poles, mattress covers, mammography equipment, stretchers, showers, sinks, and wheelchairs. **Surfaces:** For use on hard, nonporous surfaces such as chrome, fiberglass sealed surfaces, Formica, glass, glazed ceramic, laminated surfaces, plastic,^{††} stainless steel and vinyl. Avoid contact with: aluminum, brass, copper, Corian, fabric, galvanized steel, natural marble, rubber, sealed granite, silver and wood.

†† For unfamiliar surfaces, test in a small, inconspicuous area first. Consult medical equipment warranty information and equipment care and cleaning guides prior to selecting and using any surface disinfectant on plastic surfaces.

Organisms

Bacteria:

- *Acinetobacter baumannii*^{*}
- *Clostridium difficile* spores^{SS,***}
- *Enterobacter aerogenes*^{*}
- *Enterococcus faecalis*,
Vancomycin Resistant^{*}
- Multi-Drug Resistant
Enterococcus faecium^{*}
- *Enterococcus hirae*^{*}
- ESBL (Extended Spectrum
Beta Lactamase) Producing
Escherichia coli^{*}
- *Klebsiella pneumoniae*^{*}
- Methicillin Resistant
Staphylococcus aureus^{*}
- *Pseudomonas aeruginosa*^{*}
- *Salmonella enterica*^{*}
- *Serratia marcescens*^{*}
- *Staphylococcus aureus*^{*}
- *Streptococcus pyogenes*^{*}
- *Mycobacterium bovis* (TB) 20°C^{**}

*1-minute contact time
 ** 2-minute contact time
 *** 3-minute contact time
 †Kills TB in 2 minutes
 †† Bloodborne pathogens

Fungi:

- *Trichophyton mentagrophytes*^{*}

Viruses[†]:

- †Adenovirus Type 2^{*}
- †Avian Influenza A Virus^{*,§}
- †Canine Parvovirus^{*}
- †Feline Parvovirus^{*}
- †Hepatitis A Virus^{*}
- †,††Human Immunodeficiency
Virus Type 1 (HIV-1)^{*}
- †,††Hepatitis B Virus (HBV)
(as Duck HBV)^{*}
- †,††Hepatitis C Virus (HCV)
(as Bovine Diarrhea Virus)^{*}
- †Herpes Simplex Virus Type 2^{*}
- †Influenza A Virus^{*}
- †Influenza A Virus H1N1^{*}
- †Influenza B Virus^{*}
- †Norovirus (as Feline Calicivirus)^{*}
- †Poliovirus^{*}
- †Rhinovirus^{*}
- †Rotavirus^{*}

§Kills Avian Influenza A virus on
 precleaned hard, nonporous surfaces
 §§Follow special instructions for
 cleaning prior to disinfection