

Canister artwork shown at approximately 110% of actual size

LICENSED

PERIOD 2014-2016 LIC. NO.

8013.82

LICENSED

PERIOD 2011-2013 LIC. NO.



Kills *C. difficile* (*C. diff*) Spores

- Kills *Clostridium difficile* (*C. diff*) Spores
- Kills *C. diff* Spores

DISPATCH®

HOSPITAL CLEANER DISINFECTANT TOWELS WITH BLEACH



- ▶ Contains 6500 ppm Sodium Hypochlorite
- ▶ Pre-Moistened Towels With 1:10 Dilution of Bleach (6.0%) and Detergent
- ▶ Cleans, Disinfects, & Deodorizes In One Step
- ▶ Kills Pandemic 2009 H1N1 Influenza A Virus
- ▶ Kills *Clostridium difficile* Spores in 5 Minutes on Precleaned, Hard, Non-Porous Surfaces
- ▶ Ready-to-use disinfection for *C. difficile* spores*
- ▶ For daily use against *C. difficile* spores*
- ▶ To reduce *C. difficile* spores on surfaces in acute care facilities and nursing homes*

Kills 99.99% bacteria, viruses† & fungi in 1 MINUTE:

Acinetobacter baumannii, †Avian Influenza A Virus, †Canine Parvovirus, *Enterobacter aerogenes*, *Escherichia coli* ESBL, †Feline Panleukopenia Virus, †Hepatitis A Virus (HAV), †Hepatitis B Virus (HBV), †Hepatitis C Virus (HCV), †Herpes Simplex Virus (HSV-2), †HIV-1 (AIDS Virus), †Influenza A Virus, †Influenza A H1N1, *Klebsiella pneumoniae*, Methicillin resistant *Staphylococcus aureus* (MRSA), †Norovirus, †Poliovirus Type 1 (Mahoney), *Pseudomonas aeruginosa*, †Rhinovirus, †Rotavirus, *Salmonella enterica*, *Staphylococcus aureus*, *Streptococcus pyogenes*, *Trichophyton mentagrophytes* (Athlete's Foot Fungus), Vancomycin resistant *Enterococcus faecalis* (VRE)

Kills *Mycobacterium bovis* - BCG, (TB) in 2 MINUTES

Kills *Clostridium difficile* Spores in 5 MINUTES*

ACTIVE INGREDIENT:	
Sodium Hypochlorite	0.65%
OTHER INGREDIENTS:	99.35%
TOTAL:	100.00%
	(yields 0.62% available chlorine)

KEEP OUT OF REACH OF CHILDREN CAUTION

NOT TO BE USED AS BABY WIPES
KEEP CONTAINER TIGHTLY SEALED WHEN NOT IN USE
READ COMPLETE DIRECTIONS AND PRECAUTIONS ON
REAR PANELS

NET CONTENTS: 150 – 6.75" x 8" TOWELS; 2 lbs



R0917010



Shown on page 2

Shown on page 1

Shown on page 3

Canister artwork shown at approximately 110% of actual size



**DISPATCH® HOSPITAL CLEANER DISINFECTANT
TOWELS WITH BLEACH**

THIS IS NOT A BABY WIPE

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

Dispenser Directions: To Open: Remove lid. Discard inner seal from canister. **To Thread:** Find center of wipe roll, twist corner of first wipe into a point and thread through the hole in the canister lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out and tearing off. When not in use, keep center cap of lid closed to prevent moisture loss. Gloves should be worn. **††Select-A-Size:** For a 6.75" x 8" towel — pull one towel gently through the hole, angle sharply and detach at the first perforation. For a 6.75" x 16" towel — pull two towels gently straight up through the hole. Hold second towel, angle sharply and detach at the second perforation. NOTE: Select-A-Size feature for the larger size towel works best after several smaller towels have been used.

To Clean, Disinfect, & Deodorize in One Step: Wipe surface with towel until completely wet. Let stand for one minute. To kill TB, let stand for two minutes at room temperature (20 – 25°C). Wipe dry or allow to air dry. Gross soil must be removed prior to disinfecting. Gloves should be worn.

To Clean: Follow directions above, but no standing time is required. If streaking is observed, wipe with a clean, damp cloth or paper towel.

***Special Instructions for Cleaning Prior to Disinfection against Clostridium difficile spores. Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks, or eye covering. **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Wipe surface with towel until completely wet. Let stand for 5 minutes. Wipe dry or allow to air dry. Gross soil must be removed prior to disinfecting. Gloves should be worn.

Ultrasound Transducers and Probes:** For use on external surfaces of ultrasound transducers and probes only.

To Disinfect: Pre-clean gel from external surface of ultrasound transducer or probe to be disinfected with a paper towel or clean cloth. Then thoroughly wet the probe with a Dispatch Hospital Cleaner Disinfectant Towels with Bleach towel. Let stand for one minute. To kill TB, let stand for two minutes at room temperature. Efficacy tests have demonstrated that this product is an effective bactericide, virucide, tuberculocide, and fungicide in the presence of organic soil (5% blood serum).

†Kills HIV, HBV, and HCV on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV).

Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV, and HCV on Surfaces/Objects Soiled with Blood/Body Fluids. Personal Protection: When handling items soiled with blood/body fluids, personal protective equipment such as disposable gloves, gowns, masks, and eye coverings should be used. **Cleaning Procedure:** Blood/body fluids must be thoroughly cleaned from surfaces and objects before disinfecting. **Contact Time:** Let stand for one minute. **Disposal of Infectious Materials:** Blood, body fluids, and cleaning materials must be autoclaved and disposed of according to federal, state, and local regulations for infectious waste disposal.

**This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

Manufactured for:

© 2011 Clorox Professional Products Company
1221 Broadway, Oakland, CA 94612 800-234-7700
www.cloroxprofessional.com

A list of this product's ingredients is available at www.CloroxCSR.com
Clorox and Dispatch are registered trademarks of the Clorox Company



R0917010



Shown on page 2

Shown on page 1

Shown on page 3

Canister artwork shown at approximately 110% of actual size



DISPATCH® HOSPITAL CLEANER DISINFECTANT
TOWELS WITH BLEACH

THIS IS NOT A BABY WIPE

Pre-moistened disposable towels for use in healthcare and other areas where control of cross contamination is required. Meets OSHA bloodborne pathogen standards. **Surfaces:** May be used on hard non-porous surfaces such as stainless steel, glazed ceramic tile, hard plastics, non-porous vinyl, painted surfaces, glazed porcelain, sealed fiberglass surfaces, Plexiglas, laminated plastic countertops, and glass.

USE SITES:

Hospital and Healthcare: clinics, dental offices, emergency rooms (ER), hospitals, hospital central service areas, intensive care units (ICU), laboratory, nursing homes, physical therapy, physician offices, radiology, recovery rooms, and surgery. **Critical Care Areas:** ER, ICU, neonatal intensive care units (NICU), pediatric intensive care units (PICU), and surgery.

APPLICATIONS:

Hard non-porous surfaces such as autoclaves, bed railings, blood glucose monitors, cabinets, carts, chairs, counters, exam tables, gurneys, infant incubators and care cribs, IV poles, phlebotomy trays, PVC tubing, stretchers, tables, ambulance equipment, diagnostic equipment, dialysis machines, mammography equipment, mattress covers, patient monitoring equipment, respiratory equipment, ultrasound transducers and probes, and wheelchairs. **Nursery or Daycare Centers:** Changing tables, cribs, trash cans, toys, and other hard, non-porous surfaces listed on this label.

PRECAUTIONARY STATEMENTS: HAZARDS TO HUMANS AND DOMESTIC ANIMALS


CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling. **FIRST AID:** If in eyes, hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes of rinsing, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. **PHYSICAL AND CHEMICAL HAZARDS:** Do not use this product with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners. This product not intended for use on fabrics.

STORAGE AND DISPOSAL: Store this product in a cool, dry area away from direct sunlight and heat. When not in use, keep center cap of lid closed to prevent moisture loss. Nonrefillable container. Do not reuse or refill this container. Discard used towels and container in trash or take container to recycling center. Do not flush towels down toilet.

EPA Reg. No. 56392-8
EPA Est. No. 8251-OH-1 (Q6); 5813-GA-2 (VG)
Beginning of batch code indicates Est. No.

REORDER NO. 69150

See Expiration Date on Package



R0917010



Shown on page 2

Shown on page 1

Shown on page 3