



Department of Agriculture
STATE OF HAWAII

LICENSED

PERIOD 2014-2016 LIC. NO.

8153.5

Safetec

Sanizide^{PRO}

Ready-to-use Liquid Disinfectant

**Hospital / Medical / Surgical
Surface and Equipment Disinfectant / Cleaner**

BACTERICIDAL • FUNGICIDAL • VIRUCIDAL • TUBERCULOCIDAL
Kills 99.9% of bacteria - Fragrance free - Deodorizes

ACTIVE INGREDIENTS

n-Alkyl (60% C₁₄, 30% C₁₆, 5% C₁₂, 5% C₁₈) dimethyl benzyl ammonium chloride0.154%
n-Alkyl (68% C₁₂, 32% C₁₄) dimethyl ethylbenzyl ammonium chloride0.154%
Isopropanol21.000%

OTHER INGREDIENTS

.....78.692%
Total100.000%

KEEP OUT OF REACH OF CHILDREN

CAUTION

See side panel for additional Precautionary Statements

NET CONTENTS: 1 Gallon (3.785 L)

Reorder No. XXXXX Made in the U.S.A. EPA Reg. No. 70144-1-67161 EPA Est. No. 67161-NY-001

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

CAUTION: May be harmful if inhaled or absorbed through skin. May cause moderate eye irritation. Avoid breathing vapors and contact with eyes, skin or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove contaminated clothing and wash clothing before reuse.

FIRST AID

If inhaled	<ul style="list-style-type: none"> • Move person to fresh air. • If person is not breathing, call 911 or an ambulance then give artificial respiration, preferably by mouth-to-mouth, if possible.
If on skin or clothing	<ul style="list-style-type: none"> • Take off contaminated clothing. • Rinse skin immediately with plenty of water for 15-20 minutes.
If in eyes	<ul style="list-style-type: none"> • Hold eye open and rinse slowly and gently with water for 15-20 minutes. • Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.

Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. You may also contact 1-800-255-3924 for emergency medical treatment information.

Physical or Chemical Hazards

Do not use or store near heat or open flame.

STORAGE AND DISPOSAL

DO NOT CONTAMINATE WATER, FOOD, OR FEED BY STORAGE OR DISPOSAL.

PESTICIDE STORAGE- Store in original container in a cool, well-ventilated area inaccessible to children and pets.
PESTICIDE DISPOSAL- Wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility. **CONTAINER DISPOSAL-** Non-refillable container. Do not reuse or refill this container. Triple rinse as follows: Fill container 1/4 full with water and reclose the container. Agitate vigorously and dispose of rinsate in accordance with Pesticide Disposal instructions. Repeat two more times. Offer for recycling if available. Puncture and dispose of in a sanitary landfill, or by other procedures approved by state and local authorities.

Manufactured in the U.S.A. for: **SAFETEC OF AMERICA, Inc.**

887 Kensington Ave., Buffalo, NY 14215 Telephone: 1-800-456-7077

Hospital, Surgical & Home Healthcare Use

Use this product to clean, disinfect and sanitize** soft and hard non-food contact surfaces, non-porous inanimate surfaces such as those made of acrylic, metals (aluminum, brass, stainless steel, chrome), glass, plastics (polycarbonate, polypropylene, polystyrene and polyvinylchloride), clear plastics, Plexiglas®, vinyl and leather surfaces as well as laminated and painted surfaces.

May be used on items such as anesthesia machines, monitors, beds, bedrails, bedside devices, biohazard team equipment, birthing room fixtures, blood bank equipment surfaces, blood pressure equipment, chest pumps, coated mattresses, computer keyboards and screens, colon equipment surfaces, drug preparation equipment, empty whirlpool tanks***, emergency room stretchers, exam room surfaces, floors, glucometers, gonial lenses, empty hydrotherapy units, infant care equipment, isolettes®, infusion (IV) pump surfaces, oxygen hoods, patient beds, physician office surfaces, respiratory therapy equipment, skin care equipment surfaces, telemetry equipment, tonometer equipment and other similar items.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Apply this undiluted product by sponge, rag, cloth, mop, spray or immersion to wet all surfaces thoroughly.

When used as directed this product effectively kills the following microorganisms at room temperature with 2-minute disinfectant contact time against:

Acinetobacter baumannii, Avian Influenza A (H3N2) Reassortant Virus*†, Escherichia coli - with extended beta- lactamase resistance (ESBL) (E.coli), Hepatitis B Virus*† (HBV), Hepatitis C Virus*† (HCV), Herpes Simplex Virus*† (HSV-2), Human Immunodeficiency Virus† (HIV-1), Influenza A Virus*†, Klebsiella pneumoniae- Carbapenem-resistant (KPC), Listeria monocytogenes Pseudomonas aeruginosa, Rotavirus*†, Serratia marcescens, Salmonella enterica, Staphylococcus aureus, Methicillin-Drug Resistant Staphylococcus aureus (MRSA), Mycobacterium bovis BCG (tuberculosis) (TB), Trichophyton mentagrophytes (Ringworm), Vancomycin- Resistant Enterococcus faecalis (VRE), Swine Influenza A (H1N1) Virus (formerly known as "swine flu")*†
* Enveloped/lipophilic type virus ** Non-enveloped/hydrophilic type virus

Non-Food Contact Surfaces

For surface cleaning: Prior to surface disinfection, pre-clean surfaces, making sure to remove all gross filth and heavy soil, by saturating surfaces with this product by towel or paper towel then immediately wipe surface clean. Repeat, if necessary, until surfaces are visibly clean. Discard paper towels or launder cloth towels before reusing.

For cleaning walls and floors: In a clean bucket dilute 2 to 4 oz of this product per gallon of warm, clean potable water. Thoroughly damp-mop using a clean microfibre mop, cloth towel or sponge. Allow to air dry. **For use as a pre-cleaning spot or stain removal spray treatment:** Thoroughly spray soiled area, allow to remain wet for at least 10 seconds before laundering as usual.

For surface disinfection: Apply this product to pre-cleaned surfaces by sponge, cloth, mop, spray or immersion to wet all surfaces thoroughly. Allow to remain wet for 2minutes, then wipe dry using a clean mop, paper or cloth towel, or rinse surface using potable water and wipe surface dry or allow to air dry. Discard or launder spent materials before reuse.

Instruments/Devices

For use as a pre-cleaning instrument/device soak: Fill appropriately sized container with a sufficient amount of this product to allow for complete submersion of items. Place items into solution and allow to soak for 2 minutes. Be certain to flush lumens of hollow devices with this product. Remove items and rinse. For critical and semi-critical items, follow the appropriate cleaning and disinfection or sterilization process. Change solution as needed, i.e., when the solution becomes diluted or visibly soiled.

For use as a pre-cleaning spray: Place item into a suitable container. Thoroughly spray item with product until drenched, thoroughly rinse. Follow with an appropriate terminal cleaning and sterilization/disinfection process.

For use as an ultrasonic cleaning solution: Prepare the ultrasonic bath using a dilution ratio of 1 fl. oz. of this product with 1 quart of clean, potable water. Place items into the ultrasonic bath, activate unit and allow unit to operate until items are visibly clean. Remove items and thoroughly rinse under running cool tap water. Proceed with appropriate disinfection or sterilization process.

For use as a non-critical medical item disinfection immersion solution: Thoroughly pre-clean items being sure to remove all gross filth and heavy soil, then rinse and dry. Allowing for complete submersion of item(s), place pre-cleaned item(s) into an appropriate container that has been filled with this product. Flush lumens of hollow devices with this product. Allow item(s) to remain submerged for 2 minutes. Remove, rinse with clean potable water and dry. Change solution daily or when it becomes visibly soiled or discolored.

For use as an instrument spray: Place item into a suitable container. Thoroughly spray item with this product until drenched, thoroughly rinse. Follow with an appropriate terminal cleaning and sterilization/disinfection process.

HBV, HCV and HIV-1 Precautions: This product kills HBV, HCV, and HIV-1 on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of HBV, HCV and/or the human immunodeficiency virus Type 1 (HIV-1) associated with AIDS. **SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HBV, HCV and HIV-1 OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS: PERSONAL PROTECTION:** Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks, and eye coverings. **CLEANING PROCEDURE:** Blood and body fluids must be thoroughly cleaned from surfaces/objects before disinfection. **CONTACT TIME:** Allow surface to remain wet for 2 minutes. **INFECTIOUS MATERIALS DISPOSAL:** Cleaning materials used that may contain blood/body fluids must be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

This product is not to be used as a terminal sterilant or high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter sterile areas of the body. This product may be used to pre-clean or decontaminate critical and semi-critical medical devices prior to sterilization or high-level disinfection.