

A BISYNTHETIC PHENOL SOLUTION IMPREGNATED IN A UNIFORM RELEASE TOWELETTE

DIRECTIONS FOR USE:

Birex Disinfectant Wipes are intended for use in hospitals, medical, dental and veterinary offices, laboratories, industrial clean rooms, nursing homes, ambulances, police vehicles, hotels, schools, homes, health spas, grocery stores, shopping centers, boats, autos, and toilet seats. Use **Birex Disinfectant Wipes** on inanimate surfaces and equipment such as stainless steel, formica counters, glass tables, carts, baskets, counters, cabinets, telephones, grocery cart handles, grocery cart child seats, shopping cart handles, shopping cart child seats, and other hard non-porous surfaces.

Dispensing Directions: Open cap lid and pull up firmly on lid to remove the entire cap from the canister. Remove protective seal. Locate wipe at the center of the roll and insert corner of wipe through opening in lid. Replace lid. Pull sheet up and slightly to the side. The next wipe pops up automatically. Always snap lid securely shut between uses to prevent moisture loss.

IT IS A VIOLATION OF FEDERAL LAW TO USE THIS PRODUCT IN A MANNER INCONSISTENT WITH ITS LABELING. READ ENTIRE LABEL BEFORE USING THIS PRODUCT. DO NOT USE AS A BABY WIPE.

Cleaning and Disinfecting Directions: Remove one **Birex Disinfectant Wipe**, unfold, wipe surface to remove excess organic debris and thoroughly clean surface. Use a second **Birex Disinfectant Wipe** to thoroughly wet surface. Ensure surface remains wet for 10 minutes to destroy or inactivate bacteria including: Mycobacterium tuberculosis var. bovis, Staphylococcus aureus, Salmonella typhosa, Pseudomonas aeruginosa, Escherichia Coli 0113; fungi including: Trichophyton mentagrophytes, and the following viruses*: Herpes simplex virus 1 & 2, *Vaccinia virus, *Influenza virus types A1 (New Jersey), *Influenza B (Maryland) and *HIV-1 (AIDS virus). Use on non-food contact surfaces only. When using in grocery stores or shopping centers and similar locations: wipe thoroughly in order to clean non-food contact surfaces on carts, such as grocery cart handles and grocery cart child seats. Must not be used in hospital nurseries.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

CAUTION: NOT A BABY WIPE



+H810BI2401

Lemon
Scent



STATE OF HAWAII
Department of Agriculture

ACCEPTED
LICENSE NO. 8262.7
Disc.

disinfectant wipes

Tuberculocidal • Virucidal* • Bactericidal • Fungicidal
Kills Pandemic 2009 H1N1 influenza A virus (formerly called swine flu)

EPA Reg. No. 46851-10-51003, EPA Est. No. 46851-CO-003

Active Ingredients:

o-phenylphenol.....	0.28%
o-benzyl-p-chlorophenol	0.03%
Inert Ingredients.....	99.69%
Total.....	100.00%

CAUTION: KEEP OUT OF REACH OF CHILDREN

Net Contents: 2 lbs. 4 oz.
240 towelettes
pre-moistened, 6" x 6.75"

REORDER NO. BI240

Distributed By:



BIOTROLTM
www.biotrol.com
Earth City, MO 63045
(800) 822-8550



SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

BIREX DISINFECTANT WIPES KILL HIV ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings or other setting in which there is an expected likelihood of soiling inanimate surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

WHEN HANDLING BLOOD OR BODILY FLUIDS

PERSONAL PROTECTION: When handling items soiled with blood or body fluids or when conducting procedures which might lead to contact with blood or body fluids, always wear disposable latex gloves, gowns, masks, and eye coverings.

CLEANING PROCEDURE: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this disinfectant.

DISPOSAL OF INFECTIOUS MATERIAL: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

CONTACT TIME: Use a 10 minute contact time for disinfection against viruses, bacteria, and fungi claimed.

PRECAUTIONARY STATEMENT

HAZARD TO HUMANS AND DOMESTIC ANIMALS

CAUTION: Harmful if swallowed. Wash thoroughly with soap and water before eating, drinking, chewing gum or using tobacco. Harmful if inhaled. Avoid breathing vapors.

FIRST AID: *If swallowed:* Call a Poison Control Center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by Poison Control Center or doctor. Do not give anything by mouth to an unconscious person.

If inhaled: Move person to fresh air. If person is not breathing, call 911 or ambulance, then give artificial respiration, preferably mouth-to-mouth, if possible. Call a Poison Control Center or doctor for further treatment advice.

If in eyes: Hold eye open and rinse slowly and gently with water for 15 - 20 minutes. Remove contact lenses, if present, after the first 5 minutes, and then continue rinsing. Call a Poison Control Center or doctor for treatment advice.

HOTLINE NUMBER: In case of emergency, contact INFOTRAC at 1-800-535-5053. Have the product container or label with you when calling a Poison Control Center or doctor or going for treatment.

STORAGE AND DISPOSAL: Do not contaminate water, food or feed by storage or disposal.

PESTICIDE STORAGE: Store at room temperature.

PESTICIDE DISPOSAL: Do not reuse towelettes. Dispose of used towelette and empty container in regular trash.

CONTAINER DISPOSAL: Rinse and discard empty container. Do not reuse empty container.