



STATE OF HAWAII
Department of Agriculture
ACCEPTED
LICENSE NO. **8262.9**

16.6875" W

3.75" 0.75"

4.25" H

0.5"

These wipes are disposable, single-use towelette wipes pre-saturated with disinfectant solution. This product is intended for use on pre-cleaned, hard, non-porous inanimate surfaces in areas where the disinfection of surfaces is required.

USE SITES:

Birex® Quat™ Wipes may be used on dental equipment surfaces/devices, counters, curing lights, finished surfaces, eye protection devices, instrument stands, keyboards, laboratory equipment/surfaces, medical equipment surfaces/devices, office fixtures, OR lights, painted surfaces, telephones, work stations, X-ray equipment/surfaces and devices/surfaces potentially contaminated with disease-causing microorganisms; and use sites made of glass, metal, Plexiglas®, plastics (i.e. polycarbonate, polypropylene, polyvinylchloride and polystyrene), clear plastic, stainless steel, brass, chrome, copper, acrylic, and non-porous vinyl, Formica®, glazed tile, laminated surfaces and other similar surfaces.

DIRECTIONS FOR USE:

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

DO NOT USE THIS PRODUCT AS A BABY OR SKIN WIPE, OR ON CONTACT LENSES!

CLEANING AND SANITIZING/DISINFECTION OF SURFACES:

Use a fresh wipe to pre-clean surfaces of all gross filth and heavy soil. Repeat as necessary until surfaces are visibly clean. To effectively sanitize or disinfect the pre-cleaned surfaces, use a fresh new wipe or turn the wipe over to the clean side to thoroughly wet the surfaces and allow surfaces to remain wet for the appropriate time indicated for the purpose intended.

295056C Rev 04/2016 M60017-83788B-0416

PRECAUTIONARY STATEMENTS:

Hazards to Humans and Domestic Animals

CAUTION: Harmful if inhaled or absorbed through skin. Causes moderate eye irritation. Avoid breathing vapors and contact with eyes, skin or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove contaminated clothing and wash clothing before reuse.

FIRST AID: If inhaled: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible. **If on skin or clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. **If in eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a Poison Control Center or doctor for treatment advice. Have the product container or label with you when calling a Poison Control Center or doctor, or going for treatment. You may also contact 1-800-535-5053 for emergency medical treatment information.

Physical or Chemical Hazards

Do not use or store near heat or open flame.

STORAGE AND DISPOSAL:

Do not contaminate food or feed by storage or disposal. **Storage:** Store in a cool, well-ventilated area. **Product and Container Disposal:** Dispose of wipe in trash after use. Do not flush. Tightly close lid between uses to retain moisture. **Non-refillable Container:** Do not reuse or refill this empty container. Offer empty container for recycling. If recycling is not available, discard container in trash.

BIREX®
Quat™ WIPES
Dental Office Disinfectant Cleaner

Fragrance Free

Tuberculocidal & 10 Second Hard, Non-Porous Surface Sanitizer with Biodegradable Detergent

160
6 in. x 6.75 in. Wipes
Net Weight: 1 lb. 9.8 oz (730.39g)

REORDER NO. 295507

CAUTION
See side panel for additional Precautionary Statements

KEEP OUT OF REACH OF CHILDREN

ACTIVE INGREDIENTS:	
n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18)	
dimethyl benzyl ammonium chloride	0.154%
n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride	0.154%
Isopropanol	21.000%
OTHER INGREDIENTS	78.692%
TOTAL	100.000%



This product disinfects the following microorganisms at room temperature with a 2-minute contact time when used as directed: *Mycobacterium bovis BCG* (tuberculosis) (TB), *Acinetobacter baumannii*, *Escherichia coli*-with extended beta-lactamase resistance (ESBL) (E. coli), Hepatitis B Virus* (HBV), Hepatitis C Virus* (HCV), Herpes Simplex II Virus (HSV 2) G Strain, Human Immunodeficiency Virus* (HIV-1) (associated with AIDS), Influenza A Virus H1N1, Influenza Virus Strain A2/Hong Kong, *Klebsiella pneumoniae*-Carbapenem-resistant, *Listeria monocytogenes*, *Pseudomonas aeruginosa*, Rotavirus Strain WA, *Salmonella enterica*, *Serratia marcescens*, *Staphylococcus aureus*, Methicillin-Resistant *Staphylococcus aureus* (MRSA), *Trichophyton mentagrophytes* (ringworm) and Vancomycin-Resistant *Enterococcus faecium* (VRE).

SANITIZATION: This product sanitizes pre-cleaned non-food contact hard surfaces against the following microorganisms at room temperature with a 10-second contact when used as directed: *Staphylococcus aureus* and *Enterobacter aerogenes*.

***HBV, HCV & HIV-1 PRECAUTIONS:** This product effectively inactivates HBV, HCV and HIV-1 on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of HBV, HCV and/or the human immunodeficiency virus Type 1 (HIV-1) associated with AIDS. **Special instructions for cleaning and decontamination against HBV, HCV and HIV-1 of surfaces/objects soiled with blood/body fluids: Personal Protection:** Wear appropriate barrier protection including gloves, gowns, masks, and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces/objects before disinfection. **Contact Time:** Allow the surface to remain wet for 2 minutes prior to wiping. **Infectious Materials Disposal:** Cleaning materials used that may contain blood/body fluids should be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

EPA Reg. No. 70144-2-51003, EPA Est. No. 56952-WI-001

This product is not to be used as a terminal sterilant or high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter sterile areas of the body. This product may be used to pre-clean or decontaminate critical and semi-critical medical devices prior to sterilization or high-level disinfection.

DISPENSING DIRECTIONS: To start feed, remove cover and discard seal. From the center of the wipe roll, pull up a wipe corner, twist it into a point and thread it through the hole located on the container cover. Pull through about one inch. Replace cover. Pull out first wipe and tear off at an angle. Remaining wipes feed automatically, ready for the next use. When not in use, keep center cap closed to prevent evaporation.

Manufactured in the U.S.A. for:
BIOTROL™
13705 Shoreline Crt. E.
Earth City, MO 63045
(800) 822-8550



No Print Area • Glue Adhesive