



# BetONE™ Disinfectant Wipes

One Step Cleaning, Disinfecting  
and Deodorizing

Fragrance Free • Bleach Free

**KILLS PATHOGENS/  
BACTERIA/FUNGI/VIRUSES\*  
In One Minute!**

**KEEP OUT OF REACH  
OF CHILDREN  
CAUTION**

See Side Panel for Additional  
Precautionary Statements  
and First Aid.

EPA EST. NO. 10088-WI-2  
EPA REG. NO. 88494-2-4170

**NET CONTENTS:**  
1 lb. 6 oz.  
100 Pre-Moistened Wipes  
(7 in. by 9 in.)

**ACTIVE INGREDIENTS**

Ethyl Alcohol .....72.50%  
Didecyl Dimethyl Ammonium Chloride.....0.33%  
OTHER INGREDIENTS .....27.17%  
TOTAL ..... 100.00%

**DISINFECTANT  
3111B1-00**

This product will clean and disinfect hard, non-porous, surfaces of bed railings, cabinets, carts, chairs, counters, exam tables, gurneys, isolettes, IV poles, stethoscopes, stretchers, tables, telephones, toilet seats, toilets, bathtubs, bathroom fixtures, shower stalls, sinks, tiles, vanity tops, hospital beds, traction devices, MRI, CAT, scales, paddles, wheelchairs, grocery cart handles, grocery carts, shopping carts, grocery cart child seats, changing tables, cribs, diaper pails, trash-cans, garbage cans, cabinets, chairs, computers, countertops, desks, desktops, doorknobs, fax machines, garbage cans, headsets, keyboards, railings, seats, sinks, tables, tabletops, telephones, trashcans, appliance exteriors, faucets, bathtubs, shower stalls, and stove tops, gym equipment, exercise machines, non-porous wrestling/gymnastic/exercise/yoga mats, and athletic equipment, personal protective safety equipment, protective headgear, hard hats, half mask respirators, full face breathing apparatus, gas masks, goggles, spectacles, face shields, silicone rubber or PVC hearing protectors and vinyl covered ear muffs. Rinse all equipment that comes in prolonged contact with skin before reuse with clean warm water, and allow to air dry. This product is formulated to disinfect hard, non-porous, inanimate environmental surfaces such as floors, walls, metal surfaces, stainless steel surfaces, glazed porcelain, glazed ceramic tile, plastic surfaces, bathrooms, shower stalls, bathtubs, and cabinets.

**DISPENSING INSTRUCTIONS:** To Dispense Wipes: Open lid, remove entire cover, remove and discard seal. From the center of wipes roll pull up a corner of the wipe, twist it and thread it through the hole in the container cover. Pull through about one inch. Replace cover tightly. Pull out first wipe and tear off at 45° angle. Remaining wipes feed automatically, ready for next use. When not in use, keep the tip of the exposed wipe inside the well of the cover, and the lid snapped closed.

**DIRECTIONS FOR USE:** It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. **CLEANING:** Use one of more wipes, as necessary, to wet surface sufficiently and to thoroughly clean the surface. **DISINFECTING:** To disinfect hard, non-porous surfaces, use one or more wipes, as necessary, to thoroughly wet the surface to be treated. Treated surface must remain visibly wet for one minute to achieve complete disinfection of all pathogens listed on this label. For visibly soiled surfaces, clean before following disinfecting instructions. Allow surfaces to air dry. For surfaces that may come in contact with food, a potable water rinse is required. **Virucidal:** BetONE™ Disinfectant Wipes is an effective virucide on hard, non-porous surfaces, in the presence of an organic load, against Poliovirus serotype 1 [PV-1] [ATCC VR-1562], Hepatitis B Virus, Hepatitis C Virus, Rhinovirus [ATCC VR-340], Norovirus, Human Coronavirus [ATCC VR-1558], Rotavirus Strain WA, Respiratory Syncytial Virus [ATCC VR-26], Human Immunodeficiency Virus [HIV-1], Herpes Simplex Virus Type 1 [ATCC VR-539], and Influenza Virus Type A [ATCC VR-544], when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry. This product has demonstrated effectiveness against Influenza A virus and is expected to inactivate all Influenza A viruses including Pandemic 2009 H1N1 Influenza A [formerly called swine flu].

**KILLS HIV-1, HBV and HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 [HIV-1] [associated with AIDS], Hepatitis B Virus [HBV], Hepatitis C Virus [HCV]. **SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV, and HCV FOR SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS: PERSONAL PROTECTION:** When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks and eye coverings. **CLEANING PROCEDURES:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of BetONE™ Disinfectant Wipes. **CONTACT TIME:** Allow surface to remain wet for 1 full minute. **DISPOSAL OF INFECTIOUS MATERIALS:** Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

**BACTERICIDAL:** BetONE™ Disinfectant Wipes is an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Escherichia coli* 0157:H7 [ATCC 35150], *Pseudomonas aeruginosa* [ATCC 15442], *Salmonella enterica* [ATCC 10708], *Staphylococcus aureus* [ATCC 6538], *Streptococcus pyogenes* [ATCC 19615], Vancomycin resistant *Enterococcus faecium* [VRE] [ATCC 51559], Methicillin resistant *Staphylococcus aureus* [MRSA] [ATCC 33591], *Klebsiella pneumoniae* [ATCC 4352], *Listeria monocytogenes* [ATCC 984], *Acinetobacter baumannii* [ATCC 19606], *Campylobacter jejuni* [ATCC 29428], and *Burkholderia cepacia* [ATCC 25416], when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

(M)SDS No. 3111



RLB1632

©2013 Betco Corporation  
1001 Brown Avenue | Toledo, Ohio 43607  
Made In U.S.A. | All Rights Reserved.

888-GO BETCO (888-462-3826) | Betco.com



11311113

**TUBERCULOCIDAL:** BetONE™ Disinfectant Wipes is an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Mycobacterium bovis* [ATCC 35743], when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry. **FUNGICIDAL:** This product is an effective fungicide against *Trichophyton mentagrophytes* [ATCC 9533], causative agent of Athlete's Foot Fungus, when used as directed on hard surfaces found in bathrooms, shower stalls, locker rooms, exercise facilities or other clean, hard, non-porous, surfaces commonly contacted by bare feet. The surface must remain visibly wet for 1 minute. Allow surfaces to air dry. **EFFICACY TESTS HAVE DEMONSTRATED THAT BetONE™ DISINFECTANT WIPES IS AN EFFECTIVE, VIRUCIDE, BACTERICIDE, TUBERCULOCIDE AND FUNGICIDE IN THE PRESENCE OF ORGANIC SOIL.** **Food Processing and Food Packaging Facilities** Before use in federally inspected meat and poultry food processing plants and dairies, food products and packaging materials must be removed from the room or carefully protected. This product will not harm most surfaces: including stainless steel, aluminum, plastics, fiberglass, vinyl, and most painted surfaces. For plastic and painted surfaces, spot test on an inconspicuous area before use. Do not use on unpainted wood. Not for use on people or pets.

**PRECAUTIONARY STATEMENTS  
Hazards to Humans and Domestic Animals  
CAUTION:**

Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet.

**FIRST AID:**

**IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15 - 20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center, doctor, or going for treatment. For emergency information concerning this product, call the National Pesticides Information Center at 1-800-858-7378 seven days a week, 6:30 am to 4:30 pm Pacific Time. NPIC Web site: [www.npic.orst.edu](http://www.npic.orst.edu).

**PHYSICAL OR CHEMICAL HAZARDS:**

Flammable. Keep away from heat or open flame.

**STORAGE AND DISPOSAL**

Do not contaminate water, food or feed by storage or disposal.

**PESTICIDE STORAGE:**

Always store this product and all pesticides in the original container. If a leaky container must be contained within another, mark the outer container to identify the contents. Protect from extreme heat or cold.

**PESTICIDE DISPOSAL:**

Waste resulting from the use of the product must be disposed of on site or at an approved waste disposal facility.

**CONTAINER DISPOSAL:**

Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available, or discard empty container in the trash.

**LIMITED WARRANTY AND DISCLAIMER NOTICE:**

No warranty, either expressed or implied, is given to the usage contrary to the instructions given herein.

5.75"

**ATHEA**  
LABORATORIES | PACKAGING

1900 WEST CORNELL STREET  
MILWAUKEE, WI 53209  
T: 800.743.6417 | F: 414.354.9219

**WORK HISTORY**

DESCRIPTION: Dieline - 105MM standard	Vers.	Date
PRODUCT NUMBER: CAN-105-MC	2	11/19/12
CUSTOMER: Specimen		
SIZE: 5.75" h x 12.875" w - bleed 1/4"		
APPLICATION: Pressure Sensitive - auto		
COLOR(S): up to 7 color		
SIGNATURE: SP		



DIE LINE  
DOES NOT PRINT

Final Unwind →

REWIND #: 4

3 inch core  
12 inch diameter maximum