



Manufactured by:
Virox Technologies Inc.
 2770 Coventry Rd.
 Oakville, ON Canada L6H 6R1
 1-800-387-7578 | www.virox.com

**NO VARNISH
AREA**

* See reference sheet for additional features, claims and directions for use.

**FOR USE AS A ONE-STEP CLEANER/
DISINFECTANT**

1. Pre-clean heavily soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. All surfaces must remain visibly wet for 1 minute. Use a 5 minute contact time for TB and a 10 minute contact time for fungi.
4. Allow to air dry.

FOR USE AS NON-FOOD CONTACT SANITIZER

1. Pre-clean heavily soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow surface to remain wet for 30 seconds.
4. Wipe surfaces dry or allow to air dry.
5. Change towelette after each use.

PART NO. ACCWP20-TBUS
 EPA REG. NO. 74559-3
 EPA EST. NO. 74559-CAN-1

**KEEP OUT
OF REACH OF
CHILDREN**

READ M.S.D.S. BEFORE
USING PRODUCT
 M.S.D.S. is available
 online at www.virox.com

Net Contents:
4 ounces (113 grams)
20 Sheets: 6 in. x 6.75 in.
 (15.24 cm x 17.14 cm)



One-Step Disinfectant Cleaner Wipes

Bactericidal, *Virucidal, Tuberculocidal,
Fungicidal, Sanitizer

ACTIVE INGREDIENT:
 Hydrogen Peroxide0.5%
 Inert Ingredients99.5%
 Total100.0%



USA PAT# 7,632,523

FOR USE IN HEALTH CARE, LABORATORY, MEDICAL FACILITIES AND HOSPITALS.

USE ON NON-FOOD CONTACT SURFACES.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Not for personal cleansing. **THIS IS NOT A BABY WIPE.** Do not flush in toilet. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

STORAGE AND DISPOSAL:

Do not contaminate water, food or feed by storage or disposal.
For Canisters that are not intended to be refilled by consumer.

CONTAINER HANDLING: Dispose of wipe in trash after use. Do NOT flush. Tightly close lid between uses to retain moisture. Nonrefillable container. Do not reuse or refill this container. Offer empty container for recycling. If recycling is not available, discard container in trash.

AL136 (13/206)

BROAD-SPECTRUM NON-FOOD CONTACT SANITIZER

BACTERICIDAL

WIPES

***VIRUCIDAL**

TUBERCULOCIDAL

FUNGICIDAL

Department of Agriculture
STATE OF HAWAII

LICENSED

PERIOD 2013-2015 LIC. NO.

8477.2

DIELINE DOES NOT PRINT

Part No: AL136

Rev No: 13/206

Replaces: AL136 (13/138)

DECO 28 Greensboro Dr.
Etobicoke, ON
M9W 1E1

LABELS & TAGS

Proof Date: JULY 29, 2013
CLIENT: VIROX TECHNOLOGIES
CLIENT P/N: AL136 ACCEL TB
 WIPES ONE STEP DISINFECT
 CLEANER WIPES
DECO P/N: VIR003-DPN022930A
Size: 2" x 10.118"
Production Artist: CR

ARTWORK SPECIFICATIONS

PANTONE® 185	PANTONE® REF BLUE	PANTONE® Cool Gray 2	PANTONE® BLACK	SPOT UV

This proof is for CONTENT ONLY. It is not to be used for color match. Colors will be matched using the Pantone Matching System (PMS).

COPY POSITION

COPY POSITION: 1-4 LABELS WOUND OUT			
COPY POSITION: 5-8 LABELS WOUND IN			

9 ANY COPY POSITION

CLIENT ARTWORK APPROVAL AREA

PLEASE COMPLETE APPROVAL AREA IN FULL
Client assumes sole responsibility for printed artwork and copy accuracy.

<input checked="" type="checkbox"/> All Copy Correct <input checked="" type="checkbox"/> Layout Correct <input checked="" type="checkbox"/> Approved as shown	<input type="checkbox"/> Changes required
	<input type="checkbox"/> Not approved
	<p><i>All 3 must be checked to proceed with production!</i></p> <p>APPROVED <i>By Bhargavi Patel at 1:54 pm, Jul 29, 2013</i></p>

NOTE TO CLIENT: This block with signature must be faxed or returned to Deco Labels & Tags before production can begin. **Client assumes sole responsibility for printed artwork and copy accuracy.**

PROOF REV. 002 AUG 2011 / Approved by General Manager

8477.2



LICENSED PERIOD 2013-2015 LIC. NO.

Accel TB WIPES

One-Step Disinfectant Cleaner Wipes Bactericidal, *Virucidal, Tuberculocidal, Fungicidal, Sanitizer

ACTIVE INGREDIENT: Hydrogen Peroxide 0.5%, Inert Ingredients 99.5%, Total 100.0%



KEEP OUT OF REACH OF CHILDREN

See reference sheet for additional features, claims and directions for use.

WIPES

FOR USE IN HEALTH CARE, LABORATORY, MEDICAL FACILITIES AND HOSPITALS

This product is a non-woven, disposable, pre-saturated wipe with Hydrogen Peroxide disinfectant cleaner.

Accel TB Wipe is a one-step disinfectant cleaner and odor neutralizer designed for general cleaning, and disinfecting, of hard, non-porous environmental surfaces.

Accel TB Wipe's non-dulling formula eliminates the time and labor normally required for rinsing. It is designed for use on hard, non-porous environmental surfaces: tables, chairs, countertops, telephones, glazed tile and sinks.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Not for personal cleansing. This is not a baby wipe. Do not flush in toilet.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body.

FOR USE AS A CLEANER AND/OR DEODORIZER

Pull towlette from canister and wipe hard, non-porous environmental surfaces. Allow to air dry.

READ M.S.D.S. BEFORE USING PRODUCT MSDS available online at www.virox.com

FOR USE AS ONE-STEP CLEANER/ DISINFECTANT:

- 1. Pre-clean heavily soiled areas. 2. Pull towlette from canister and wipe hard, non-porous environmental surfaces. 3. All surfaces must remain visibly wet for 1 minute. Use a 5 minute contact time for Tuberculosis (TB) and a 10 minute contact time for fungi. 4. Allow to air dry.

*KILLS HIV AND HCV AND HBV ON PRE-CLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings and other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) Hepatitis B Virus (HBV), or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV OR HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Disposable latex or vinyl gloves, gowns, face masks, and eye coverings as appropriate, must be worn during all cleaning of body fluids, blood, and decontamination procedures.

Cleaning Procedures: Blood/body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Allow surface to remain wet for 1 minute to kill HIV-1, HBV and HCV. Use a 5 minute contact time for Tuberculosis (TB) and a 10 minute contact time for fungi.

Disposal of Infectious Material: Blood/body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

PART NO. ACCWIPB-TBUS EPA REG. NO. 74559-3 EPA EST. NO. 74559-CAN-1

Accel

Net Contents: 1.1 pounds (493 grams) 60 Sheets: 7 in. x 8 in. (17.8 cm x 20.3 cm)

BROAD-SPECTRUM NON-FOOD CONTACT SANITIZER

FUNGICIDAL • TUBERCULOCIDAL • *VIRUCIDAL • BACTERICIDAL

BACTERICIDAL • *VIRUCIDAL • TUBERCULOCIDAL • FUNGICIDAL

BROAD-SPECTRUM NON-FOOD CONTACT SANITIZER



Proof Date: AUG 9, 2013 CLIENT: VIROX TECHNOLOGIES CLIENT P/N: AL 127 ACCEL RTU US WIPES DECO P/N: VIR003-DPN023264 Size: 5.709" x 10.118" Production Artist: CR

ARTWORK SPECIFICATIONS: Pantone color swatches (185, REF BLUE, Cool Gray 2, BLACK, SPOT UV), copy position diagrams (1-4 labels wound out, 5-8 labels wound in), and a '9 ANY COPY POSITION' label.

CLIENT ARTWORK APPROVAL AREA: PLEASE COMPLETE APPROVAL AREA IN FULL. Client assumes sole responsibility for printed artwork and copy accuracy. Includes checkboxes for 'All Copy Correct', 'Layout Correct', 'Approved as shown', 'Changes required', and 'Not approved'. Signature: APPROVED By Bhargavi Patel at 9:20 am, Aug 12, 2013. NOTE TO CLIENT: This block with signature must be faxed or returned to Deco Labels & Tags before production can begin.

*This product has demonstrated effectiveness against Influenza A virus and is expected to inactivate all Influenza A viruses including Pandemic 2009 H1N1 Influenza A virus.

TUBERCULOCIDAL: 5 Min. In the presence of 5% serum load and 5 minute contact time at 20°C on hard, non-porous environmental surfaces. Mycobacterium bovis

FUNGICIDAL: 10 Min. Trichophyton mentagrophytes

NON-FOOD CONTACT SANITIZER FOR HARD NON-POROUS SURFACES : 30 Sec. In the presence of 5% serum load and 30 second contact time at 20°C on hard, non-porous environmental surfaces.

Escherichia coli O157:H7 Klebsiella pneumoniae Pseudomonas aeruginosa Salmonella enterica, (formerly known as Salmonella choleraesuis) Staphylococcus aureus

Notice to User: Accel TB may not be suitable for all items or surfaces. Prolonged exposure to copper, brass, anodized aluminum and some coatings should be avoided.

Manufactured by: Virox Technologies Inc. 2770 Coventry Rd. Oakville, ON Canada L6H 6R1 1-800-387-7578 www.virox.com



NON VARNISH AREA

AL127 (13/221)

Part No: AL127 Rev: 13/221 Replaces: AL130 Latest revision (Replace only after instructed by Virox)

DIELINE DOES NOT PRINT



FOR USE AS A ONE-STEP CLEANER/DISINFECTANT

1. Pre-clean heavily soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. All surfaces must remain visibly wet for 1 minute. Use a 5 minute contact time for TB and a 10 minute contact time for fungi.
4. Allow to air dry.

FOR USE AS NON-FOOD CONTACT SANITIZER

1. Pre-clean heavily soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow surface to remain wet for 30 seconds.
4. Wipe surfaces dry or allow to air dry.
5. Change towelette after each use.

READ M.S.D.S. BEFORE USING PRODUCT
M.S.D.S. is available online at www.virox.com



Manufactured by:
Virox Technologies Inc.
2770 Coventry Rd.
Oakville, ON Canada L6H 6R1
1-800-387-7578
www.virox.com

PART NO.
ACCWIPTBUS-12X12RP
EPA REG. NO.
74559-3
EPA EST. NO.
74559-CAN-1

Accel

Net Contents:
4 lbs. 2.6 oz.
(1.89 kg)
110 sheet
12 in. x 12 in.
(30.48 cm x 30.48 cm)

Accel[®] TB
WIPES

LICENSED
PERIOD 2013-2015 LIC. NO.

One-Step Disinfectant Cleaner Wipes

Bactericidal, *Virucidal, Tuberculocidal,
Fungicidal, Sanitizer

ACTIVE INGREDIENT:

Hydrogen Peroxide	0.5%
Inert Ingredients	99.5%
Total	100.0%



USA PAT#: 7,632,523

KEEP OUT OF REACH OF CHILDREN

WIPES

FOR USE IN HEALTH CARE, LABORATORY, MEDICAL FACILITIES, AND HOSPITALS

Use on non-food contact surfaces.

DIRECTIONS FOR USE: STORAGE AND DISPOSAL:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Not for personal cleansing. **THIS IS NOT A BABY WIPE.** Do not flush in toilet. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

REFILL POUCH: Look for product expiration date on back of refill pouch. Remove lid from empty bucket. Place refill pouch into bucket. Tear open refill pouch. **DO NOT** remove wipe roll from refill pouch. Pull out corner of first wipe from centre of roll and push 1/2" to 3/4" (1 to 2 cm) through "X" opening on lid. **DO NOT** push finger through opening. Replace lid and pull wipe at an angle. The next wipe pops up automatically. When finished, close lid to retain moisture. If bucket becomes contaminated or damaged replace immediately. Replace bucket immediately if tabs below lid become damaged.

For buckets, lids, or pouches intended to be refilled by the consumer.

CONTAINER HANDLING: Dispose of wipe in trash after use. Do not flush in toilet. Tightly reseal bucket flap between uses to retain moisture. Do not reuse or refill except as described in the directions for use. Refill only with this product. Offer empty container for recycling. If recycling is not available, wrap and discard in trash.

* See reference sheet for additional features, claims and directions for use.



AL151(14/010)

BACTERICIDAL

*VIRUCIDAL

TUBERCULOCIDAL

FUNGICIDAL

BROAD-SPECTRUM NON-FOOD CONTACT SANITIZER



Accel[®] TB WIPES

PART NO.
ACCWIP11x12-TBUS
15/075

EPA REG. NO.
74559-3
EPA EST. NO.
74559-CAN-1



USA PAT#: 7,632,523

KEEP OUT OF REACH OF CHILDREN

Net Contents: 6.2 lbs. (2.8 kg)
160 Wipes, 11 in. x 12 in. (28 cm x 30 cm)

BACTERICIDAL
*VIRUCIDAL
TUBERCULOCIDAL

FUNGICIDAL
BROAD-SPECTRUM NON-FOOD
CONTACT SANITIZER

BACTERICIDAL: 1 MIN.
*VIRUCIDAL: 1 MIN.
TUBERCULOCIDAL: 5 MIN.

FUNGICIDAL: 10 MIN.
BROAD-SPECTRUM NON-FOOD
CONTACT SANITIZER: 30 SEC.

BACTERICIDAL: 1 MIN.
*VIRUCIDAL: 1 MIN.
TUBERCULOCIDAL: 5 MIN.

FUNGICIDAL: 10 MIN.
BROAD-SPECTRUM NON-FOOD
CONTACT SANITIZER: 30 SEC.



LICENSED

PERIOD 2013-2015 L.I.C. NO.

8477.2

FOR USE AS ONE-STEP CLEANER/DISINFECTANT

1. Pre-clean heavily soiled areas.
2. Pull towelette from dispenser and wipe hard, non-porous environmental surfaces.
3. All surfaces must remain visibly wet for 1 minute. Use a 5 minute contact time for Tuberculosis (Tb) and a 10 minute contact time for fungi.
4. Allow to air dry.

FOR USE AS A NON-FOOD CONTACT SANITIZER

1. Pre-clean heavily soiled areas.
2. Pull towelette from dispenser and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain wet for 30 seconds.
4. Wipe surfaces dry or allow to air dry.
5. Change towelette after each use.

READ S.D.S. BEFORE USING PRODUCT
S.D.S. is available online at www.virox.com



Manufactured by:
Virox Technologies Inc.
2770 Coventry Rd.
Oakville, ON Canada L6H 6R1
1-800-387-7578 | www.virox.com

BACTERICIDAL: 1 MIN.
*VIRUCIDAL: 1 MIN.
TUBERCULOCIDAL: 5 MIN.

FUNGICIDAL: 10 MIN.
BROAD-SPECTRUM NON-FOOD
CONTACT SANITIZER: 30 SEC.

BACTERICIDAL
*VIRUCIDAL
TUBERCULOCIDAL

FUNGICIDAL
BROAD-SPECTRUM NON-FOOD
CONTACT SANITIZER

BACTERICIDAL: 1 MIN.
*VIRUCIDAL: 1 MIN.
TUBERCULOCIDAL: 5 MIN.

FUNGICIDAL: 10 MIN.
BROAD-SPECTRUM NON-FOOD
CONTACT SANITIZER: 30 SEC.

FOR USE IN HEALTH CARE, LABORATORY, MEDICAL FACILITIES, AND HOSPITALS

Use on non-food contact surfaces.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Not for personal cleansing. **THIS IS NOT A BABY WIPE.** Do not flush in toilet. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

STORAGE AND DISPOSAL:

Store this product in a cool dry area.

For buckets intended to be refilled by the consumer.

CONTAINER HANDLING: Dispose of wipe in trash after use. Do not flush in toilet. Tightly reseal bucket flap between uses to retain moisture. Do not reuse or refill except as described in the directions for use. Refill only with this product. Offer empty container for recycling. If recycling is not available, wrap and discard in trash.

See reference sheet for additional features,
claims and directions for use.



**FOR USE AS ONE-STEP
CLEANER/DISINFECTANT**

1. Pre-clean heavily soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. All surfaces must remain visibly wet for 1 minute. Use a 5 minute contact time for Tuberculosis (Tb) and a 10 minute contact time for fungi.
4. Allow to air dry.

**FOR USE AS A NON-FOOD
CONTACT SANITIZER**

1. Pre-clean heavily soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain wet for 30 seconds.
4. Wipe surfaces dry or allow to air dry.
5. Change towelette after each use.

READ S.D.S. BEFORE USING PRODUCT
S.D.S. is available online at www.virox.com



Manufactured by:
Virox Technologies Inc.
2770 Coventry Rd.
Oakville, ON Canada L6H 6R1
1-800-387-7578
www.virox.com

PART NO.
ACCWIP11x12RP-TBUS
15/075

EPA REG. NO.
74559-3

EPA EST. NO.
74559-CAN-1

Accel

Net Contents:
6.2 lbs.
(2.8 kg)

160 Wipes,
11 in. x 12 in.
(28 cm x 30 cm)



One-Step Disinfectant Cleaner Wipes

Bactericidal, *Virucidal, Tuberculocidal,
Fungicidal, Sanitizer

ACTIVE INGREDIENT:

Hydrogen Peroxide	0.5%
Inert Ingredients.....	99.5%
Total	100.0%



USA PAT#: 7,632,523

KEEP OUT OF REACH OF CHILDREN

WIPES

**FOR USE IN HEALTH CARE, LABORATORY,
MEDICAL FACILITIES, AND HOSPITALS**

Use on non-food contact surfaces.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Not for personal cleansing. **THIS IS NOT A BABY WIPE.** Do not flush in toilet. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semicritical medical devices prior to sterilization or high level disinfection.

REFILL POUCH: Look for product expiration date on back of refill pouch. Remove lid from empty bucket. Place refill pouch into bucket. Tear open refill pouch. DO NOT remove wipe roll from refill pouch. Pull out corner of first wipe from centre of roll and push 1/2" to 3/4" (1 to 2 cm) through "X" opening on lid. DO NOT push finger through opening. Replace lid and pull wipe at an angle. The next wipe pops up automatically. When finished, close lid to retain moisture. If bucket becomes contaminated or damaged replace immediately. Replace bucket immediately if tabs under lid become damaged.

STORAGE AND DISPOSAL:

For buckets, lids, or pouches intended to be refilled by the consumer.

CONTAINER HANDLING: Dispose of wipe in trash after use. Do not flush in toilet. Tightly reseal bucket flap between uses to retain moisture. Do not reuse or refill except as described in the directions for use. Refill only with this product. Offer empty container for recycling. If recycling is not available, wrap and discard in trash.

See reference sheet for additional features, claims and directions for use.



LICENSED

PERIOD 2013-2015 LIC. NO.

8477.2



8 07174 53440 5

BACTERICIDAL

*VIRUCIDAL

TUBERCULOCIDAL

FUNGICIDAL

BROAD-SPECTRUM NON-FOOD CONTACT SANITIZER

AL174 (15/079)

Accel[®] TB WIPES

Reference Sheet



Manufactured by:
Virox Technologies Inc.
2770 Coventry Rd., Oakville, ON L6H 6R1
1-800-387-7578 | www.viroxaccel.com



FOR USE IN HEALTH CARE, LABORATORY, MEDICAL FACILITIES AND HOSPITALS

Accel TB is a one-step disinfectant cleaner and odor neutralizer designed for general cleaning, and disinfecting, of hard, non-porous environmental surfaces. It cleans quickly by removing dirt, grime, mold, mildew, and other common soils found in hospitals, nursing homes, and medical offices.

Accel TB's non-dulling formula eliminates the time and labor normally required for rinsing. It is designed for use on hard, non-porous environmental surfaces listed on this label made from: plastic, chrome, stainless steel, aluminum associated with walls, tables, chairs, countertops, telephones, glazed tile, sinks.

BACTERICIDAL: 1 MINUTE

HOSPITAL DISINFECTANT: BACTERICIDAL in the presence of 5% serum load and 1 minute contact time on hard, non-porous environmental surfaces.

BACTERICIDAL ACTIVITY:

Acinetobacter baumannii [ATCC 19606]

Escherichia coli O157:H7 [ATCC 35150]

Pseudomonas aeruginosa [ATCC 15442]

Salmonella enterica [ATCC 10708] formerly known as *Salmonella choleraesuis*

Staphylococcus aureus [ATCC 6538]

ANTIBIOTIC-RESISTANT BACTERIA:

Enterococcus faecalis VRE [ATCC 51575]

Staphylococcus aureus MRSA [ATCC 33592]

*VIRUCIDAL: 1 MINUTE

***VIRUCIDAL** in the presence of 5% serum load and 1 minute contact time on hard, non-porous environmental surfaces.

*VIRUSES:

*Herpes Simplex Virus, Type 1 [HSV-1] [ATCC VR-733]

*Herpes Simplex Virus, Type 2 [HSV-2] [ATCC VR-734]

*Human Coronavirus [ATCC VR-740]

*Influenza A/Hong Kong [ATCC VR-544]

*Poliovirus Type 1, Strain Brunhilde [ATCC VR-1000]

*Rhinovirus Type 14 [ATCC VR-284]

*Rotavirus WA

*Kills HIV-1 [AIDS virus] [HTLV-III_g] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

*Kills HBV & HCV when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

*Kills Norovirus [Feline Calicivirus as the surrogate] [VR-782] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

*This product has demonstrated effectiveness against Influenza A virus and is expected to inactivate all Influenza A Viruses including Pandemic 2009 H1N1 Influenza A virus.

VETERINARY *VIRUSES:

*Avian Influenza A, H3N2 [ATCC VR-2072]

*Feline Calicivirus, Strain F9 [ATCC VR-782]

TUBERCULOCIDAL: 5 MINUTES

TUBERCULOCIDAL in the presence of 5% serum load and 5 minute contact time on hard, non-porous environmental surfaces.

Mycobacterium bovis [Tb] [OT 451C150]

FUNGICIDAL: 10 MINUTES

FUNGICIDAL This product is fungicidal against the pathogenic fungi, *Trichophyton mentagrophytes* [ATCC 9533] [the fungus which causes Athlete's Foot] in the presence of 5% serum load and 10 minute contact time when used as directed on hard surfaces found in bathroom, shower stalls, locker rooms, or other clean, non-porous, hard surfaces commonly contacted by bare feet.

Trichophyton mentagrophytes [ATCC 9533]

BROAD-SPECTRUM NON-FOOD CONTACT SANITIZER: 30 SECONDS

NON-FOOD CONTACT SURFACE SANITIZER FOR HARD NON-POROUS SURFACES – in the presence of 5% serum load and 30 second contact time on hard, non-porous environmental surfaces.

Enterococcus faecalis VRE [ATCC 51575]

Escherichia coli O157:H7 [ATCC 35150]

Klebsiella pneumoniae [ATCC 4352]

Pseudomonas aeruginosa [ATCC 15442]

Staphylococcus aureus [ATCC 6538]

Staphylococcus aureus MRSA [ATCC 33592]

Salmonella enterica [ATCC 10708] formerly known as *Salmonella choleraesuis*

***KILLS HIV AND HCV AND HBV ON pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids** in health care settings [Hospitals, Nursing Homes] and other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/ body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) Hepatitis B Virus (HBV), or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV OR HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Disposable latex or vinyl gloves, gowns, face masks, and eye coverings as appropriate, must be worn during all cleaning of body fluids, blood, and decontamination procedures.

Cleaning Procedures: Blood/body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Allow surface to remain wet for 1 minute to kill HIV-1, HBV and HCV. Use a 5 minute contact time for Tuberculosis [Tb] and 10 minute contact for fungi.

Disposal of Infectious Material: Blood/body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

The logo for Accel TB Wipes features a grid of blue dots above the text "Accel TB Wipes" in red and blue, with "One-Step Disinfectant Wipes" in blue below it.

Accel[®] TB Wipes

One-Step Disinfectant Wipes

EPA Reg. No. 74559-3

Animal Housing Facilities

See Reference Sheet for complete list of organisms

This product cleans by removing dirt, grime, blood, urine, fecal matter and other common soils found in animal housing facilities, grooming facilities, farms, kennels, pet stores, veterinary clinics, laboratories or other small animal facilities. It also eliminates odor leaving surfaces smelling clean and fresh.

Directions for Use

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

For cleaning and disinfecting the following hard nonporous surfaces: feeding and watering equipment, utensils, instruments, cages, kennels, stables, catteries, etc.

1. Remove all animals and feeds from premises, animal transportation vehicles, crates, etc.
2. Remove all litter, droppings and manure from floors, walls and surfaces of facilities (including pens, chutes, and barns) occupied or traversed by animals.
3. Empty all troughs, racks, and other feeding and watering appliances.
4. Thoroughly clean all surfaces with soap and detergent and rinse with potable water. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
5. All surfaces must remain visibly wet for 1 minute.
6. Allow to air dry
7. Ventilate buildings and other closed spaces.
8. Do not house animals or employ equipment until treatment has been absorbed, set or dried.
9. Thoroughly scrub all treated feed racks, troughs, mangers, automatic feeders, fountains, and waterers with soap or detergent, and rinse with potable water before reuse.

NOTE: Not for use on handling and restraining equipment such as leashes, muzzles, halters or ropes, as well as forks, shovels, and scrapers used for removing litter and manure.

Environmental Hazards

This product is toxic to birds, fish and aquatic invertebrates. Caution should be used when applying indoors because pets may be at risk.



BACTERICIDAL • *VIRUCIDAL • TUBERCULOCIDAL • FUNGICIDAL • SANITIZER

Accel

Accel[®] TB
WIPES

One-Step Disinfectant Cleaner Wipes

Bactericidal, *Virucidal, Tuberculocidal, Fungicidal, Sanitizer

ACTIVE INGREDIENT:

Hydrogen Peroxide 0.5%
Inert Ingredients 99.5%
Total 100.0%

Net Contents:

2.6 ounces (75 grams)

10 Sheets:

7.48 in. x 10 in. (19.0 cm x 25.4 cm)

KEEP OUT OF REACH OF CHILDREN.

See reference sheet for additional features, claims, and directions for use.

Manufactured For:
Virox Technologies Inc.
1-800-387-7578 | www.virox.com

PEEL & SEAL / PULL HERE



Department of Agriculture
STATE OF HAWAII

EPA REG. NO. 74559-3

EPA EST. NO. 075272-WI-001

AL156 (14/112)

LICENSED 8477.2

PERIOD 2013-2015 LIC. NO.

Accel[®] TB WIPES

SOFT PACK

ACCWIPSP-TBUS
Reference Sheet

KEEP OUT OF REACH OF CHILDREN
EPA EST. NO. 075272-WI-001

FOR USE IN HEALTH CARE, LABORATORY, MEDICAL FACILITIES AND HOSPITALS

Accel TB is a one-step disinfectant cleaner and odor neutralizer designed for general cleaning, and disinfecting, of hard, non-porous environmental surfaces. It cleans quickly by removing dirt, grime, mold, mildew, and other common soils found in hospitals, nursing homes, and medical offices.

Accel TB's non-dulling formula eliminates the time and labor normally required for rinsing. It is designed for use on hard, non-porous environmental surfaces listed on this label made from: plastic, chrome, stainless steel, aluminum associated with walls, tables, chairs, countertops, telephones, glazed tile, sinks - any washable non-food contact surface where disinfection is required.

BACTERICIDAL: 1 MINUTE

HOSPITAL DISINFECTANT: BACTERICIDAL in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

BACTERICIDAL ACTIVITY:

Acinetobacter baumannii [ATCC 19606]
Escherichia coli O157:H7 [ATCC 35150]
Klebsiella pneumoniae (ATCC 4352)
Pseudomonas aeruginosa [ATCC 15442]
Salmonella enterica [ATCC 10708] formerly known as *Salmonella choleraesuis*
Shigella dysenteriae (ATCC 11835)
Staphylococcus aureus [ATCC 6538]

ANTIBIOTIC-RESISTANT BACTERIA:

Community Associated Methicillin Resistant *Staphylococcus aureus* – CA-MRSA (NARSA NRS 384) (Genotype USA 300)
Community Associated Methicillin Resistant *Staphylococcus aureus* – CA-MRSA (NARSA NRS 123) (Genotype USA 400)
Enterococcus faecalis VRE [ATCC 51575]
Escherichia coli with extended beta-lactamase resistance (ATCC BAA-196) (ESBL)
Staphylococcus aureus MRSA [ATCC 33592]

*VIRUCIDAL: 1 MINUTE

***VIRUCIDAL** in the presence of 5% serum load and 1 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

*VIRUSES:

*Adenovirus type 8 (ATCC VR-1368, Strain Trim)
*Herpes Simplex Virus, Type 1 [HSV-1] [ATCC VR-733]
*Herpes Simplex Virus, Type 2 [HSV-2] [ATCC VR-734]
*Human Coronavirus [ATCC VR-740]
*Influenza A/Hong Kong [ATCC VR-544]



Manufactured for:
Virox Technologies Inc.
2770 Coventry Rd., Oakville, ON L6H 6R1
1-800-387-7578 | www.virox.com



*Poliovirus Type 1, Strain Brunhilde [ATCC VR-1000]

*Rhinovirus Type 14 [ATCC VR-284]

*Rotavirus WA

*Kills HIV-1 [AIDS virus] [HTLV-III_g] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

*Kills HBV & HCV when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

*Kills Norovirus [Feline Calicivirus as the surrogate] [VR-782] when used as directed on hard, non-porous inanimate surfaces with a 1 minute contact time.

*This product has demonstrated effectiveness against Influenza A virus and is expected to inactivate all Influenza A Viruses including Pandemic 2009 H1N1 Influenza A virus.

VETERINARY *VIRUSES:

*Avian Influenza A, H3N2 [ATCC VR-2072]

*Feline Calicivirus, Strain F9 [ATCC VR-782]

TUBERCULOCIDAL: 5 MINUTES

TUBERCULOCIDAL in the presence of 5% serum load and 5 minute contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Mycobacterium bovis [Tb] [OT 451C150]

FUNGICIDAL: 10 MINUTES

FUNGICIDAL This product is fungicidal against the pathogenic fungi, *Trichophyton mentagrophytes* [ATCC 9533] [the fungus which causes Athlete's Foot] in the presence of 5% serum load and 10 minute contact time at 20°C (68°F) when used as directed on hard surfaces found in bathroom, shower stalls, locker rooms, or other clean, non-porous, hard surfaces commonly contacted by bare feet.

Trichophyton mentagrophytes [ATCC 9533]

BROAD-SPECTRUM NON-FOOD CONTACT SANITIZER: 30 SECONDS

NON-FOOD CONTACT SURFACE SANITIZER FOR HARD NON-POROUS SURFACES – in the presence of 5% serum load and 30 second contact time at 20°C (68°F) on hard, non-porous environmental surfaces.

Enterococcus faecalis VRE [ATCC 51575]

Escherichia coli O157:H7 [ATCC 35150]

Klebsiella pneumoniae [ATCC 4352]

Pseudomonas aeruginosa [ATCC 15442]

Staphylococcus aureus [ATCC 6538]

Staphylococcus aureus MRSA [ATCC 33592]

Salmonella enterica [ATCC 10708] formerly known as *Salmonella choleraesuis*

Use on non-food contact surfaces

DIRECTIONS FOR USE:

It is a violation of Federal Laws to use this product in a manner inconsistent with its labeling.

Not for personal cleansing. This is not a baby wipe. Do not flush in toilet.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body.

For Use as One-Step Cleaner/Disinfectant:

1. Pre-clean heavily soiled areas
2. Pull towelette from dispenser and wipe hard, non-porous environmental surfaces.
3. All surfaces must remain visibly wet for 1 minute. Use a 5 minute contact time for TB and a 10 minute contact time for fungi.
4. Allow to air dry.

For Use as a Non-Food Contact Sanitizer:

1. Pre-clean heavily soiled areas.
2. Pull towelette from dispenser and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain wet for 30 seconds.
4. Wipe surfaces dry or allow to air dry.
5. Change towelette after each use.

STORAGE AND DISPOSAL:

Store this product in a cool area.

For soft-packs that are **not intended to be refilled** by consumer.

CONTAINER HANDLING: Dispose of wipe in trash after use. Do not flush. Tightly close flap between uses to retain moisture. Nonrefillable container. Do not reuse or refill this container. Offer empty container for recycling. If recycling is not available, discard container in trash.

READ M.S.D.S BEFORE USING PRODUCT

M.S.D.S available on line at: www.virox.com

***KILLS HIV AND HCV AND HBV ON pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids** in health care settings and other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) Hepatitis B Virus (HBV), or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV OR HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Disposable latex or vinyl gloves, gowns, face masks, and eye coverings as appropriate, must be worn during all cleaning of body fluids, blood, and decontamination procedures.

Cleaning Procedures: Blood/body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Allow surface to remain wet for 1 minute to kill HIV-1, HBV and HCV. Use a 5 minute contact time for Tuberculosis [TB] and 10 minute contact for fungi.

Disposal of Infectious Material: Blood/body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.