



CAT. # GC110

**KILLS C DIFF SPORES IN 1 STEP IN 3 MINS.
100 GLUCO-CHLOR® 3" X 3" BLEACH WIPES**

SPORICIDAL ♦ BACTERICIDAL ♦ FUNGICIDAL ♦ +VIRUCIDAL

Contains 1:10 dilution (6500 ppm) of Sodium hypochlorite

One step cleaner and disinfectant that kills *Clostridium difficile* spores, MRSA, VRE, *Staphylococcus aureus*, *Salmonella enterica*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes*, +HIV-1 and +HBV. Pre-clean prior to disinfecting surfaces for +HCV.

EPA Reg. No. 69687-1
EPA Est. No. 69687-IL-001

800-647-7180

sales@medtrol.com
www.medtrol.com



STATE OF HAWAII
Department of Agriculture
ACCEPTED
LICENSE NO. **8538.1**
Disc.

100 Pre-saturated Wipes
3" x 3" (7.62cm x 7.62cm)

CAT. # GC110

KILLS C DIFF SPORES IN 1 STEP IN 3 MINS.

GLUCO-CHLOR®

SPORICIDAL ♦ BACTERICIDAL ♦ FUNGICIDAL ♦ +VIRUCIDAL

Disinfects hard, non porous healthcare surfaces, blood glucose meters†, patient care equipment and point-of-care equipment

One step cleaner and disinfectant that kills *Clostridium difficile* spores, MRSA, VRE, *Staphylococcus aureus*, *Salmonella enterica*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes*, +HIV-1 and +HBV. Pre-clean prior to disinfecting surfaces for +HCV.

Kills
C DIFF spores
in ONE STEP
in 3 mins.

Active Ingredient: Sodium hypochlorite.....0.65%
Other Ingredients*.....99.35%
Total.....100.00%

*Other ingredients do not include the weight of wipe.

Kills
+HBV
+HCV
in 30 sec.



Effective against
Blood Borne Pathogens
Kills +HIV-1, +HBV, +HCV

CAT. # GC110

- Gluco-Chlor® meets OSHA blood borne pathogen standards in 29 CFR 1910-1030
- Gluco-Chlor® meets the CDC guidelines for disinfection of blood glucose meters
- Gluco-Chlor® meets the CMS guidelines for Tag F441 in 42 CFR 483.65 for disinfection of blood glucose meters.
- Gluco-Chlor® is a one-step cleaner and disinfectant that is sporicidal and meets the CDC, CMS Tag F441, SHEA-IDSA, APIC and Joint Commission guidelines for killing *Clostridium difficile* spores.



Medtrol, Inc.
7157 North Austin Ave., Niles, IL 60714

800-647-7180
sales@medtrol.com
www.medtrol.com

EPA Reg. No. 69687-1
EPA Est. No. 69687-IL-001
Made in U.S.A.

**KEEP OUT OF REACH OF CHILDREN
CAUTION**

(See side panels for Suggested Areas of Usage, First Aid Statement, Special Instructions for Cleaning and Decontamination against +HIV-1, +HBV and +HCV, Precautionary Statements, and Storage and Disposal)

Net Contents 10 oz.
(100 X 0.1 oz.)
U.S. Patent No. 6,530,477
©2013 Medtrol, Inc.



8 93081 00002 7

38"

SUGGESTED AREAS OF USAGE: Medical, Dental and Laboratory Counters, Exam Tables, Blood Glucose Meters†, Carts, Patient Care Equipment, Point-of-Care Equipment, Digital Screens, on Plastics, Glazed Porcelain, Stainless Steel, Formica, Vinyl, and Rubber Surfaces.

FIRST AID: If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have the product container with you when calling a poison control center or doctor, or going for treatment.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

DISINFECTION:

1. Always use personal protective equipment.+
2. Open Gluco-Chlor® packet.
3. Remove pre-saturated 3" x 3" wipe.
4. Apply wipe and wipe desired surface to be disinfected.
5. A 30 second contact time is required to kill HBV and HCV. A 3 minute contact time is required to kill *Clostridium difficile* spores. A 5 minute contact time is required to kill HIV and other organisms listed on the label. Reapply as necessary to ensure that the surface remains wet for entire contact time.
6. Allow surface to air dry and discard used wipe and empty packet (see storage and disposal).

***KILLS HIV-1, HBV AND HCV ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of HIV-1(associated with AIDS), HBV and HCV.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV AND HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

- **PERSONAL PROTECTION:** When handling items soiled with blood or body fluids, use disposable gloves, gown, masks and eye protection.
- **CLEANING PROCEDURES:** Blood/body fluids must be thoroughly cleaned from surfaces/objects before application of Gluco-Chlor®.
- **CONTACT TIME:** Allow surface to remain wet for 30 seconds to kill HBV and HCV, for 3 minute to kill *Clostridium difficile* spores, and 5 minutes to kill HIV and other organisms listed on the label.
- **DISPOSAL OF INFECTIOUS MATERIAL:** Blood/body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

Complete product line information available by calling:
 Medtrol, Inc. 7157 North Austin Avenue, Niles, Illinois 60714
 Ph: 800-647-7180 Fax: 800-255-3027
 sales@medtrol.com www.medtrol.com



PRECAUTIONARY STATEMENTS:

Hazards to humans and domestic animals

CAUTION: Causes moderate eye irritation. Avoid contact with eyes and clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove and wash contaminated clothing before reuse.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barriers or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

†This product is for the exterior surfaces of Glucometers and is not for use on the needle and/or the monitors.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal

Pesticide Storage: Store in a cool, dry, and well-ventilated area away from heat, open flame, direct sunlight and inaccessible to children.

Pesticide (wipe) Disposal: Wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

Container Handling: Dispose of wipe in trash after use. Do not flush. Offer empty packet for recycling if available or discard in trash.

Rotate dated products while in inventory.