



STATE OF HAWAII  
Department of Agriculture

**ACCEPTED**

LICENSE NO. **8538.2**

**Disc.**

BLEACH WIPES



CAT. # SC160



# SUPER◆CHLOR®

Ready to use EPA registered bleach wipes  
DISINFECTS • CLEANS • DEODORIZES



SPORICIDAL ◆ BACTERICIDAL ◆ FUNGICIDAL ◆ +VIRUCIDAL

Contains:1:10 dilution (6500ppm) of Sodium hypochlorite

One step cleaner and disinfectant that kills *Clostridium difficile* spores, MRSA, VRE, *Staphylococcus aureus*, *Salmonella enterica*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes*, +HIV-1 and +HBV. Preclean prior to disinfecting surfaces for +HCV.

- Kills 99.99% of germs
- These wipes are highly absorbent and tear resistant
- For use to disinfect hard, non-porous healthcare and environmental surfaces
  - Patient Care Areas
  - Patient Care Equipment
  - Point-of-Care Equipment
  - Glucometers†
  - Doctor's Offices • Dental Offices
  - Emergency Medical Services
  - Nursing Homes • Schools
  - Laboratories

Meets OSHA blood borne pathogen standards in 29 CFR 1910.1030 and guidelines of the CDC and CMS Tag F441. Aids in reducing cross contamination between treated surfaces.

A one step cleaning and disinfecting sporicidal bleach wipe that meets the CDC, CMS Tag F441, SHEA-IDS, APIC and Joint Commission guidelines for killing *Clostridium difficile* spores.

Active Ingredient: Sodium hypochlorite .....	0.65%
Other Ingredients*.....	99.35%
Total.....	100.00%

\*Other ingredients do not include the weight of wipes.

## KEEP OUT OF REACH OF CHILDREN CAUTION

(See side panels for directions for use and additional precautionary statements.)

EPA Reg. No. 69687-1  
EPA Est. No. 69687-IL-001  
©2013 Medtrol, Inc.

MADE IN U.S.A.

NET CONTENTS:

1 lb. 5 oz.  
160 Pre-saturated Wipes  
6" x 6" (15.2cm x 15.2cm)



4.500"

Disinf  
SUGGES  
Equipment  
FIRST A  
if presen  
Have the

Directi  
It is a vi

Disinf

1. A
2. C
3. R
4. A
5. A
- to
- li
6. A

To open

1. A
2. R
3. R
4. L
5. P
6. D
7. R
8. S
9. D

Pesticid  
away fro

Pesticid  
disposal

Contain  
after emp  
allowed

K

287

185

278

14.875"

### Disinfecting wipe for hard, non-porous healthcare surfaces and environmental surfaces.

**SUGGESTED AREAS OF USAGE:** On Medical, Dental and Laboratory Counters, Exam Tables, Carts, Patient Care Equipment, Point-of-Care Equipment, Digital Screens, on Plastics, Glazed Porcelain, Stainless Steel, Formica, Vinyl, Rubber Surfaces, and Blood Glucose Meters†.

**FIRST AID:** If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have the product container with you when calling a poison control center or doctor, or going for treatment.

### Directions for use:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

### Disinfection:

1. Always use personal protective equipment.\*
2. Open Super-Chlor® canister.
3. Remove pre-saturated 6"x 6" wipe.
4. Apply wipe and wipe desired surface to be disinfected.
5. A 30 second contact time is required to kill HBV and HCV. A 3 minute contact time is required to kill *Clostridium difficile* spores. A 5 minute contact time is required to kill HIV and other organisms listed on the label. Reapply as necessary to ensure that the surface remains wet for entire contact time.
6. Allow surface to air dry and discard used wipe and empty canister (see storage and disposal).

### To open canister:

1. Always use personal protective equipment.\*
2. Remove lid from canister.
3. Remove seal if present and discard.
4. Lift flap in canister lid.
5. Pull up corner of center wipe and snap off, twist to a point and thread through dispenser hole in lid.
6. **DO NOT PUSH FINGER THROUGH OPENING.**
7. Replace lid. Pull out exposed wipe and snap off. The next wipe pops up automatically.
8. Snap center flap down when finished to retain moisture.
9. Do not flush wipe into toilet.

## STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal

**Pesticide Storage:** Keep wipes in a tightly closed canister, when not in use. Store canister in cool, dry, and well-ventilated area away from heat, open flame, direct sunlight and inaccessible to children.

**Pesticide (wipe) Disposal:** Wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

**Container Handling:** Nonrefillable rigid canister. Do not re-use or refill this container. Triple rinse (or equivalent), promptly after emptying. Offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration, or, if allowed by state and local authorities, by burning. If burned, stay out of smoke.

**\*KILLS HIV-1, HBV AND HCV ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of HIV-1 (associated with AIDS), HBV and HCV.

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV AND HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.**

- **PERSONAL PROTECTION:** When handling items soiled with blood or body fluids, use disposable gloves, gown, masks and eye protection.
- **CLEANING PROCEDURES:** Blood/body fluids must be thoroughly cleaned from surfaces/objects before application of Super-Chlor®.
- **CONTACT TIME:** Allow surface to remain wet for 30 seconds to kill HBV and HCV, for 3 minute to kill *Clostridium difficile* spores, and 5 minutes to kill HIV and other organisms listed on the label.
- **DISPOSAL OF INFECTIOUS MATERIAL:** Blood/body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

### PRECAUTIONARY STATEMENTS:

Hazards to humans and domestic animals

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes and clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove and wash contaminated clothing before reuse.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barriers or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

† This product is for the exterior surfaces of GLUcometers and is not for use on the needle and/or the monitors.

Complete product line information available by calling:

Medtrol, Inc.  
7157 North Austin Avenue  
Niles, Illinois 60714  
Ph: 800-647-7180 Fx: 800-255-3027  
sales@medtrol.com www.medtrol.com

