



MEDLINE

MICRO-KILL™ Bleach

GERMICIDAL BLEACH SOLUTION



**VIRUCIDAL

- Adenovirus Type 5
- Avian Influenza A (H3N2)
- Avian Influenza A (H5N1)
- Avian Influenza A (H7N9)
- Canine Parvovirus
- Feline Panleukopenia Virus
- Hepatitis A Virus
- +Hepatitis B Virus
- +Hepatitis C Virus
- Herpes Simplex Virus Type 2
- Human Coronavirus
- +Human Immunodeficiency Virus Type 1
- Influenza A Virus (Novel H1N1)
- Influenza A Virus (Hong Kong)
- Norovirus (Feline Calicivirus)
- Polio Virus Type 1
- Respiratory Syncytial Virus
- Rhinovirus Type 37
- Rotavirus

*SPORICIDAL

- Clostridium difficile* spores

FUNGICIDAL

- Candida albicans*
- Trichophyton mentagrophytes*

BACTERICIDAL

- Acinetobacter baumannii*-MDR
- Bordetella pertussis*
- Campylobacter jejuni*
- CA-Methicillin Resistant
- Staphylococcus aureus*
- Enterobacter aerogenes*
- Enterobacter cloacae*
- Enterococcus faecalis*-VRE
- Enterococcus faecium*-MDR
- Escherichia coli*
- Escherichia coli* O157:H7
- Escherichia coli*-ESBL
- Escherichia coli*-NDM-1

- Klebsiella pneumoniae*
- Klebsiella pneumoniae*-ESBL
- Klebsiella pneumoniae*-CRKP (CRE)
- Klebsiella pneumoniae*-NDM-1
- Legionella pneumophila*
- Linezolid Resistant *Staphylococcus aureus*
- Listeria monocytogenes*
- Methicillin Resistant *Staphylococcus aureus* (MRSA)
- Proteus mirabilis*
- Pseudomonas aeruginosa*
- Salmonella enterica*
- Serratia marcescens*
- Shigella dysenteriae*
- Staphylococcus aureus*
- Staphylococcus aureus*-VISA
- Staphylococcus aureus*-VRSA
- Streptococcus pneumoniae*
- Streptococcus pyogenes*

This solution is sporicidal*, virucidal**, bactericidal and fungicidal and can be used for disinfecting hard, non-porous environmental and healthcare surfaces.

1:10 Bleach (Sodium Hypochlorite) Solution

Active Ingredient: Sodium Hypochlorite.....0.65%

Other Ingredients.....99.35%

Total.....100.00%

EPA Reg. No. 37549-2 EPA Est. No. 10088-WI-01

KEEP OUT OF REACH OF CHILDREN CAUTION

(See side panels for directions for use, first aid, and additional precautionary statements.)

SUGGESTED AREAS OF USAGE:

On Medical, Dental, Laboratory Counters, Exam Tables, Carts, Patient Care Equipment, Point-of-Care Equipment, Telephones, Sink Tops, Toilet Seats, Digital Screens, Plastics, Glazed Porcelain, Glazed ceramic tiles, Stainless Steel, Formica, Vinyl, Hard Plastic Surfaces, Laminated Countertops, Glass, Plexiglas, Suitable For Use on Equipment Requiring Bleach Based Products, Bed Rails, IV Poles, Dialysis Machines, and Blood Glucose Meters†.

DIRECTIONS FOR USE:

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

HOSPITAL DISINFECTION:

1. Always use personal protective equipment.*
2. Open Micro-Kill Bleach Germicidal Bleach Solution container.
3. Replace the closure with the trigger sprayer.
4. Spray Micro-Kill Bleach Germicidal Bleach Solution on desired surface.
5. Scrub and wipe with towel of your choice over desired surface to be disinfected. A 30 second contact time is required to kill all of the bacteria and viruses** on the label except a 1 minute contact time is required to kill *Candida albicans* and *Trichophyton mentagrophytes* and a 5 minutes contact time is required to kill *Clostridium difficile* spores*. Reapply as necessary to ensure that the surface remains wet for the entire contact time.
6. Allow surface to air dry and discard empty container (see storage and disposal).

***KILLS HIV-1, HBV AND HCV ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of HIV-1 (associated with AIDS), HBV and HCV.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV, AND HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

- PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use disposable gloves, gown, masks and eye protection.
- CLEANING PROCEDURES: Blood/body fluids must be thoroughly cleaned from surfaces/objects before application of Micro-Kill Bleach Germicidal Bleach Solution.
- CONTACT TIME: Allow surface(s) to remain wet for 30 seconds to kill all of the bacteria and viruses** on the label except a 1 minute contact time is required to kill *Candida albicans* and *Trichophyton mentagrophytes* and a 5 minute contact time is required to kill *Clostridium difficile* spores*.
- DISPOSAL OF INFECTIOUS MATERIAL: Blood/body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

*SPECIAL INSTRUCTIONS FOR CLEANING PRIOR TO DISINFECTING AGAINST CLOSTRIDIUM DIFFICILE SPORES.

- PERSONAL PROTECTION: Wear appropriate barrier protection such as gloves, gown, mask, or eye covering.
- CLEANING PROCEDURES: Fecal matter /waste must be thoroughly cleaned from surfaces /objects before disinfection by application with Micro-Kill Bleach Solution. Cleaning is to include vigorous wiping and /or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled Solution.
- INFECTIOUS MATERIAL DISPOSAL: Material used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barriers or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

† This product is for the exterior surfaces of Blood Glucose Meters and is not for use on the needle and/or the monitors.

FIRST AID:

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have the product container with you when calling a poison control center or doctor, or going for treatment.

PRECAUTIONARY STATEMENTS:

Hazards to humans and domestic animals

CAUTION: Causes moderate eye irritation. Avoid contact with eyes and clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

Physical and Chemical Hazards

This product contains bleach. Do not use this product with other chemicals such as ammonia, toilet bowl cleaners, rust removers, or acid, as this releases hazardous gases.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal

Pesticide Storage: Store in a cool, dry area away from direct sunlight and heat to avoid deterioration. In case of spill, flood areas with large quantities of water.

Product Disposal: Products or rinsates that cannot be used must be diluted with water before disposal in a sanitary sewer.

Container Handling and Disposal:

REFILLABLE CONTAINERS. Refill this container with pesticide only. Do not reuse this container for any other purpose. To clean the container before final disposal; Triple rinse (or equivalent), promptly after emptying. Offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration, or, if allowed by state and local authorities, by burying. If burned, stay out of smoke.



Manufactured for Medline Industries Inc., One Medline Place, Mundelein, IL 60060. Micro-Kill is a trademark of Medline Industries, Inc. 1-800-MEDLINE. Made in the USA. RJ15ATH For product emergency, contact: Poison Control Center, 1-800-222-1222 For Technical Information call: 1-800-MEDLINE, For an SDS, call 1-800-MEDLINE and request document EVSCHEM910

REF EVSCHEM910



(01) 10888277362284

NET CONTENTS: 1 Quart (32 fl. oz.) 946 mL

ATHEA
LABORATORIES | PACKAGING

1900 WEST CORNELL STREET
MILWAUKEE, WI 53209
T: 800.743.6417 | F: 414.354.9219

WORK HISTORY

DESCRIPTION: Dieline - 120MM standard	Vers.	Date
PRODUCT NUMBER: BQDW	2	09.29.14
CUSTOMER: Die Line		
SIZE: 4.75"h x 10"w - 1/8" Bleed		
APPLICATION: Pressure Sensitive - Flexographic		
COLOR(S): up to 6 colors		
SIGNATURE: BL		

4.75"