

ACTIVATE™

5.25% INSTITUTIONAL

BLEACH



LICENSED

PERIOD 2012-2014 LIC. NO. For use with the Activate™ Bleach Dilution Sprayer

HOSPITAL CLEANER & DISINFECTANT

- Cleans & Disinfects in One Step
- For Hospital, Veterinary and Food-contact Use on Hard, Non-porous Surfaces
- For Instrument Pre-cleaning & Disinfection of Non-critical Instruments
- Kills HIV-1 & HBV in 30 Seconds
- Kills Bacteria in 2 Minutes

Kills Clostridium difficile Spores ✓
in 4 minutes

Active Ingredient:

Sodium Hypochlorite. 5.25%

Other Ingredients: 94.75%

TOTAL: 100.00%

Contains 5% available chlorine by weight.

PEEL OPEN
HERE

NET CONTENTS 11 OZ (325 mL)

ACTIVATE 5.25% INSTITUTIONAL BLEACH

The Activate™ Bleach Dilution Sprayer automatically dilutes and mixes this product to a ready-to-use 10% solution (1:10 dilution with 1 part bleach to 9 parts water) equal to approximately 5,000 ppm available chlorine for hospital, laboratory and veterinary use on hard, non-porous surfaces.

- Cleans & Disinfects in One Step**
- Kills Clostridium difficile Spores in 4 Minutes**
- Kills HIV-1 & HBV in 30 Seconds**
- Kills Bacteria in 2 Minutes:**

- *Staphylococcus aureus*
- *Methicillin resistant Staphylococcus aureus (MRSA)*
- *Pseudomonas aeruginosa*
- *Vancomycin resistant Enterococcus faecalis (VRE)*
- *Acinetobacter baumannii*
- *Salmonella enterica*
- *Streptococcus pyogenes*

FIRST AID

If in eyes: • Hold eye open and rinse slowly and gently with water for 15-20 minutes. • Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.

If on skin or clothing: • Take off contaminated clothing. • Rinse skin with plenty of water for 15-20 minutes.

If swallowed: • Have person sip a glass of water if able to swallow. • Do not induce vomiting unless told to do so by the poison control center or doctor. • Do not give anything by mouth to an unconscious person.

Call the poison control center at 1-800-222-1222 for emergency medical treatment information. Have the product container or label with you when calling the poison control center or doctor, or going for treatment.

Note to Physician: Probable mucosal damage may contraindicate the use of gastric lavage.

EPA Reg No 75266-1

EPA Est No 75266-OR-001

Manufactured by:
 Deardorff Fitzsimmons Corporation
 PO Box 539
 Merlin, OR 97532
 541-476-6065

Keep Out of Reach of Children

DANGER

Read complete directions and precautions inside before using.



Department of Agriculture
STATE OF HAWAII

LICENSED

PERIOD 2015-2017 LIC. NO.

8727.1

ACTIVATE™ 5.25% INSTITUTIONAL BLEACH

For use in hospitals, laboratories, nursing homes, clinics, physicians offices, dental operatories, veterinary hospitals, schools, institutions, reception areas, kitchen areas, and restrooms. For use on hard, non-porous surfaces such as glazed ceramic tile, stainless steel, chrome, enameled walls & ceilings, laminated plastic countertops & cabinets, sealed fiberglass, glass, and glazed porcelain.

Kills HIV-1 & HBV in 30 Seconds

Kills Bacteria in 2 Minutes:

- *Staphylococcus aureus*
- *Methicillin resistant Staphylococcus aureus (MRSA)*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Streptococcus pyogenes*
- *Vancomycin resistant Enterococcus faecalis (VRE)*
- *Acinetobacter baumannii*

Kills Clostridium difficile Spores in 4 minutes

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Install Activate™ Bleach cartridge onto Activate™ Bleach Dilution Sprayer by following Sprayer Instructions below.

To Clean and Disinfect in One Step: Spray Activate™ 6"-8" from surface until surface is thoroughly wet, or apply with a clean cloth, mop or sponge saturated with Activate™; let stand for 2 minutes. To kill HIV-1 and HBV, allow 30 seconds contact time at room temperature. Wipe with a clean cloth or paper towel and allow to air dry. For food contact surfaces, rinse with potable water. Always remove gross dirt and debris from surface prior to disinfection.

Special Instructions for Cleaning Prior to Disinfecting against Clostridium difficile spores: Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks, or eye covering.

Cleaning Procedure: Thoroughly clean fecal matter/waste from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with Activate™. Vigorously wipe and/or scrub surfaces/objects until all visible soil is removed. Pay special attention to high-touch surfaces. Clean surfaces in patient rooms in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Clean restrooms last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Immediately dispose of materials used in the cleaning process that may contain feces/waste in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Spray Activate™ 6"-8" from surface until surface is thoroughly wet, or apply with a clean cloth, mop or sponge saturated with Activate™; let stand for 4 minutes.

Wipe with a clean, damp cloth or paper towel, or allow to air dry.

Kills HIV-1 and HBV on pre-cleaned environmental surfaces and objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) and Hepatitis B Virus (HBV).

Special Instructions for Cleaning and Decontaminating Against HIV-1 and HBV on Surfaces/Objects Soiled with Blood/Body Fluids: When handling items soiled with blood/body fluids, personal protective equipment such as disposable gloves, gowns, masks, and eye coverings must be used. Blood/body fluids must be thoroughly cleaned from surfaces/objects prior to disinfecting. Spray Activate™ on surface/object until all surfaces are thoroughly wet, or apply with a clean cloth, mop or sponge saturated with Activate™; let stand for 30 seconds. Wipe with a clean cloth or paper towel, and allow to air dry. Blood/body fluids and cleaning materials must be autoclaved and disposed of according to federal, state, and local regulations for infectious waste disposal.

To Pre-clean Instruments Prior to Terminal Sterilization:

Place instruments in a suitable container. Spray Activate™ on instruments until all surfaces are thoroughly wet. Allow to stand for up to 10 minutes. Rinse instruments with water, and follow with an appropriate sterilization/high-level disinfecting process. **Cleaning of critical and semi-critical devices must be followed by an appropriate terminal sterilization/high-level disinfection process.**

Use Activate™ on stainless steel instruments only.

To disinfect non-critical, pre-cleaned instruments:

Instruments first must be thoroughly cleaned to remove gross dirt and debris, rinsed, and dried. Place instruments in a suitable container. Spray Activate™ on instruments until all surfaces are thoroughly wet; let stand for 2 minutes. To kill HIV-1 and HBV, allow 30 seconds contact time at room temperature. To kill Clostridium difficile spores, let stand for 4 minutes. Wipe with a clean, damp cloth or paper towel, and allow to air dry.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

To Clean and Disinfect in a Veterinary Application: Use to clean and disinfect hard, non-porous surfaces such as feeding and watering equipment, cages, utensils, instruments, kennels, stables, catteries, etc. Remove all animals and feed from premises, animal transportation vehicles, crates, etc. Remove all litter, droppings, and manure from walls, floors, and surfaces of facilities occupied or traversed by animals. Empty all feeding and watering equipment. Pre-clean all surfaces with soap or detergent and rinse with water. Spray Activate™ on surface/object until all surfaces are saturated; let stand for 2 minutes. To kill HIV-1 and HBV, allow 30 seconds at room temperature. Rinse all surfaces with potable water. Ventilate buildings and other closed spaces. Do not house animals or employ equipment

until treated surfaces have been thoroughly rinsed with water and allowed to dry. Thoroughly scrub all treating, feeding, and watering appliances with soap or detergent, and rinse with potable water before reuse.

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

DANGER:

CORROSIVE: Causes irreversible eye damage and skin burn. Avoid breathing spray vapors. Vacate poorly ventilated areas as soon as possible until strong odors have dissipated. Wear safety glasses and rubber gloves when handling this product. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

Physical or Chemical Hazards

Strong oxidizer. Flush drains before and after use. Do not use or mix with other chemicals, such as toilet bowl cleaners, rust removers, acid or ammonia containing products. To do so will release hazardous gases. Prolonged contact with metal may cause pitting or discoloration.

STORAGE & DISPOSAL

Do not contaminate water, food or feed by storage and disposal.

Storage: This product degrades with age; test active level of diluted product periodically with high-level chlorine test strips. Store in a cool (50°–72°F), dry area away from direct sunlight and heat to avoid deterioration. In case of spill, flood area with large quantities of water.

Product Disposal: Products or rinsates that cannot be used must be diluted with water before disposal in a sanitary sewer.

Container Handling: Nonrefillable container. Do not reuse or refill this container. Offer container for recycling if available or place in trash.



Manufactured by:
Deardorff Fitzsimmons Corporation
PO Box 539
Merlin, OR 97532
541-476-6065



WATER BOTTLE

For use only with the Activate™ Bleach Dilution Sprayer

The Activate™ Bleach Dilution Sprayer automatically dilutes and mixes Activate™ 5.25% Institutional Bleach to a ready-to-use 10% solution (1:10 dilution with 1 part bleach to 9 parts water) equal to approximately 5,000 ppm available chlorine for hospital, laboratory and veterinary use on hard, non-porous surfaces.

Activate™ Bleach Dilution Sprayer Instructions

1. Fill the sprayer's empty water bottle with cool tap water, and then insert the sprayer's tube down into the filled water bottle. Close the sprayer tab to lock the water bottle into place on the sprayer.
2. Remove the threaded plastic cap from the Activate™ bleach cartridge. Do not puncture or tamper with the plug that is in the neck of the bleach cartridge.
3. Lock the bleach cartridge onto the other side of the sprayer.
4. Pointing the sprayer away from you into a sink or basin, pump the trigger 20 times to prime the system.
5. The sprayer is now ready to use. Hold the sprayer 6" to 8" from the surface, and spray thoroughly until surface is completely wet.
6. Refill the water bottle when empty.
7. Always re-prime the system when installing a new bleach cartridge.
Note that a safety feature ensures the sprayer will not function if either the bleach cartridge or the water bottle is empty.

The water in this bottle is not for human consumption.

Container Handling: Refillable container. Refill this container with water only. Do not re-use this container for any other purpose. Offer for recycling if available or place in trash.

Manufactured by



Deardorff Fitzsimmons Corporation
PO Box 539
Merlin, OR 97532
541-476-6065



WATER BOTTLE



ACTIVATE™ 5.25% INSTITUTIONAL BLEACH

For use in hospitals, laboratories, nursing homes, clinics, physicians offices, dental operators, veterinary hospitals, schools, institutions, reception areas, kitchen areas, and restrooms. For use on hard, non-porous surfaces such as glazed ceramic tile, stainless steel, chrome, enameled walls & ceilings, laminated plastic countertops & cabinets, sealed fiberglass, glass, and glazed porcelain.

Kills HIV-1 & HBV in 30 Seconds

Kills Bacteria in 2 Minutes:

- *Staphylococcus aureus* • *Methicillin resistant Staphylococcus aureus (MRSA)*
- *Pseudomonas aeruginosa* • *Salmonella enterica* • *Streptococcus pyogenes*
- *Vancomycin resistant Enterococcus faecalis (VRE)* • *Acinetobacter baumannii*

Kills Clostridium difficile Spores in 4 minutes

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

DANGER

CORROSIVE: Causes irreversible eye damage and skin burn. Avoid contact with eyes, skin and clothing. Avoid breathing spray vapors. Vacate poorly ventilated areas as soon as possible until strong odors have dissipated. Wear safety glasses and rubber gloves when handling this product. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse.

Physical or Chemical Hazards

Strong oxidizer. Flush drains before and after use. Do not use or mix with other chemicals, such as toilet bowl cleaners, rust removers, acid or ammonia containing products. To do so will release hazardous gasses. Prolonged contact with metal may cause pitting or discoloration.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage and disposal.

Storage: This product degrades with age; test active level of diluted product periodically with high-level chlorine test strips. Store in a cool (50° – 72°F), dry area away from direct sunlight and heat to avoid deterioration. In case of spill, flood area with large quantities of water.

Product Disposal: Products or rinsates that cannot be used must be diluted with water before disposal in a sanitary sewer.

Container Handling: Nonrefillable container. Do not reuse or refill this container. Clean container promptly after emptying. Triple rinse as follows: Fill container 1/4 full with water and recap. Shake for 10 seconds. Drain for 10 seconds after the flow begins to drip. Dispose of rinsate in sanitary sewer. Offer for recycling if available or place in trash.

Manufactured by  Deardorff Fitzsimmons Corp PO Box 539, Merlin, OR 97532 541-476-6065
EPA Reg. No. 75266-1 EPA Est. No. 75266-OR-001



LICENSED

PERIOD 2012-2014 LIC. NO.



5.25% Institutional BLEACH

Hospital Cleaner & Disinfectant
Cleans & Disinfects in One Step

Kills Clostridium difficile Spores in 4 Minutes

Active Ingredient:

Sodium Hypochlorite 5.25%

Other Ingredients: 94.75%

TOTAL 100.00%

Contains 5% available chlorine by weight

Keep Out of Reach of Children

DANGER

Read complete directions and precautions before using.

FIRST AID

If in eyes: • Hold eye open and rinse slowly and gently with water for 15-20 minutes. • Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.

If on skin or clothing: • Take off contaminated clothing. • Rinse skin with plenty of water for 15-20 minutes.

If swallowed: • Have a person sip a glass of water if able to swallow. • Do not induce vomiting unless told to do so by the poison control center or doctor. • Do not give anything by mouth to an unconscious person.

Call the poison control center at 1 800 222-1222 for emergency medical treatment information. Have the product container or label with you when calling the poison control center or doctor or going for treatment.

Note to Physician: Probable mucosal damage may contraindicate the use of gastric lavage.

NET CONTENTS ONE GALLON (128 OZ)



LICENSED

PERIOD 2015-2017 LIC. NO.

8727.1

DIRECTIONS FOR USE It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Use 1 part Activate™ 5.25% Institutional Bleach to 9 parts water to achieve a 1:10 dilution (~5000 ppm available chlorine). This can be prepared by diluting 14 oz. (1 3/4 cup) Activate™ 5.25% Institutional Bleach per gallon of water.

To Clean and Disinfect in One Step: Spray Activate™ 6-8" from the surface until surface is thoroughly wet, or apply with a clean cloth, mop or sponge saturated with Activate™; let stand for 2 minutes. To kill HIV-1 and HBV, allow 30 seconds contact time at room temperature. Wipe with a clean cloth or paper towel and allow to air dry. For food contact surfaces, rinse with potable water. Always remove gross dirt and debris from surface prior to disinfection.

Special Instructions for Cleaning Prior to Disinfecting against Clostridium difficile spores: Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks, or eye covering. **Cleaning Procedure:** Thoroughly clean fecal matter/waste from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with Activate™. Vigorously wipe and/or scrub surfaces/objects until all visible soil is removed. Pay special attention to high-touch surfaces. Clean surfaces in patient rooms in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Clean restrooms last. Do not reuse soiled cloths. **Infectious Materials Disposal:** Immediately dispose of materials used in the cleaning process that may contain feces/waste in accordance with local regulations for infectious materials disposal. **Disinfection Procedure:** Spray Activate™ 6- 8" from surface until surface is thoroughly wet, or apply with a clean cloth, mop or sponge saturated with Activate™; let stand for 4 minutes. Wipe with a clean, damp cloth or paper towel, or allow to air dry.

To Pre-clean Instruments Prior to Terminal Sterilization: Place instruments in a suitable container. Spray Activate™ on instruments until all surfaces are thoroughly wet. Allow to stand for up to 10 minutes. Rinse instruments with water, and follow with an appropriate sterilization/high-level disinfecting process.

Cleaning of critical and semi-critical devices must be followed by an appropriate terminal sterilization/high-level disinfection process.

Use Activate™ on stainless steel instruments only.

To disinfect non-critical, pre-cleaned instruments: Instruments first must be thoroughly cleaned to remove gross dirt and debris, rinsed, and dried. Place instruments in a suitable container. Spray Activate™ on instruments until all surfaces are thoroughly wet; let stand for 2 minutes. To kill HIV-1 and HBV, allow 30 seconds contact time at room temperature. To kill *Clostridium difficile* spores, let stand for 4 minutes. Wipe with a clean, damp cloth or paper towel, and allow to air dry.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Kills HIV-1 and HBV on pre-cleaned environmental surfaces and objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) and Hepatitis B Virus (HBV).

Special Instructions for Cleaning and Decontaminating Against HIV-1 and HBV on Surfaces/Objects Soiled with Blood/Body Fluids: When handling items soiled with blood/body fluids, personal protective equipment such as disposable gloves, gowns, masks, and eye coverings must be used. Blood/body fluids must be thoroughly cleaned from surfaces/objects prior to disinfecting. Spray Activate™ on surface/object until all surfaces are thoroughly wet, or apply with a clean cloth, mop or sponge saturated with Activate™; let stand for 30 seconds. Wipe with a clean cloth or paper towel, and allow to air dry. Blood/body fluids and cleaning materials must be autoclaved and disposed of according to federal, state, and local regulations for infectious waste disposal.

To Clean and Disinfect in a Veterinary Application: Use to clean and disinfect hard, non-porous surfaces such as feeding and watering equipment, cages, utensils, instruments, kennels, stables, catteries, etc. Remove all animals and feed from premises, animal transportation vehicles, crates, etc. Remove all litter, droppings, and manure from walls, floors, and surfaces of facilities occupied or traversed by animals. Empty all feeding and watering equipment. Pre-clean all surfaces with soap or detergent and rinse with water. Spray Activate™ on surface/object until all surfaces are saturated, or apply with a clean cloth, mop or sponge saturated with Activate™; let stand for 2 minutes. To kill HIV-1 and HBV, allow 30 seconds at room temperature. Rinse all surfaces with potable water. Ventilate buildings and other closed spaces. Do not house animals or employ equipment until treated surfaces have been thoroughly rinsed with water and allowed to dry. Thoroughly scrub all treating, feeding, and watering appliances with soap or detergent, and rinse with potable water before reuse.