



Run Direction # 3

30% Black

DECON-PHENE®

Sterile Pharmaceutical Clean Room Formula

Germicidal Detergent with Organic Soil Tolerance

Tuberculocidal • Fungicidal • Virucidal *

ACTIVES:

Sodium ortho-benzyl parachlorophenate	4.40%
Sodium ortho-phenylphenate	2.82%
Sodium para-terta-amyphenate	2.49%
Inert Ingredients	90.29%
	100%

EPA REG. NO.211-36-68959 EPA EST, NO. 68959-PA-1

Net Contents: 1 Gallon

DANGER
KEEP OUT OF REACH
OF CHILDREN

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals • DANGER

Corrosive. Causes irreversible eye damage or skin burns. May be fatal if swallowed. Harmful if absorbed through skin or inhaled. Do not get in eyes, on skin or on clothing. Avoid breathing spray mist. Wear goggles or face shield (safety glasses). Wear protective clothing (long sleeve shirt and long pants, socks plus shoes and chemical resistant gloves such as water proof gloves).

FIRST AID

If in eyes	<ul style="list-style-type: none"> Hold eye open and rinse slowly and gently with plenty of water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.
If inhaled	<ul style="list-style-type: none"> Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible. Call a poison control center or doctor for treatment advice.
If on skin or clothing	<ul style="list-style-type: none"> Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.
If swallowed	<ul style="list-style-type: none"> Call a poison control center or doctor for treatment advice. Do not induce vomiting unless told to do so by a poison control center or doctor. Do not give any liquid to this person Do not give anything by mouth to an unconscious person

HOTLINE NUMBER

For emergency medical advice, call your local poison control center or a doctor.
Have the product container or label with you when seeking medical advice or treatment.

NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage. Measure against circulatory shock, respiratory depression, and convulsion may be needed.
Product Questions? Call 1-888-478-3745

Container and Product Sterilized and Distributed by:

Veltek Associates, Inc. 15 Lee Blvd. Malvern, PA 19355-1234 Tel: 610-644-8335 Fax: 610-644-8336

Made in U.S.A.



EFFICACY PERFORMANCE

The broad spectrum efficacy of DECON-PHENE has been confirmed using A.O.A.C. Protocol testing at a use-dilution of 1:128 in the presence of 5% blood serum (used to simulate organic soil load) at 20 degrees Celsius with an exposure time of 10 minutes.

GERMICIDAL: At 1 ounce per gallon of water, DECON-PHENE is effective against a broad spectrum of Gram positive and Gram negative pathogens as confirmed by the A.O.A.C. Use-Dilution test in the presence of 5% organic soil load at 20° C with an exposure time of 10 minutes against the following organisms:

Staphylococcus aureus	Pseudomonas aeruginosa	Salmonella choleraesuis
Saccharomyces cerevisiae	Nesisseria meningitidis	Hemophilus influenzae
Bacillus subtilis	Staphylococcus epidermidis	Streptococcus pneumoniae
Escherichia coli		

FUNGICIDAL: DECON-PHENE is effective against the following pathogenic fungi, as confirmed by the A.O.A.C. Use-Dilution test in the presence of 5% organic soil load at 20° C with an exposure time of 10 minutes, in use areas such as aseptic filling and gowning rooms, and general manufacturing areas, or any other hard, nonporous, inanimate environmental surfaces:

Trichophyton mentagrophytes	Candida albicans
-----------------------------	------------------

TUBERCULOCIDAL: DECON-PHENE is effective against Mycobacterium bovis (BCG) as demonstrated through testing using the A.O.A.C. Method "Confirmative In Vitro Test for Determining Tuberculocidal Activity", in the presence of 5% organic soil load, at 20°C with an exposure time of 15 minutes.

***VIRUCIDAL:** DECON-PHENE is effective against the following viruses on hard, nonporous, inanimate surfaces, as confirmed by the A.O.A.C. use-dilution test, in the presence of 5% organic soil load, at 20° C with an exposure time of 10 minutes:
HIV-1 (AIDS virus)

***VIRUCIDAL:** "KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type-1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 (AIDS VIRUS) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

PERSONAL PROTECTIVE EQUIPMENT (PPE) Applicators and other handlers must wear: Long-sleeved shirt and long pants; socks plus shoes; goggles or face shield and chemical resistant gloves (such as barrier laminate, butyl rubber, nitrile rubber, neoprene rubber, polyvinyl chloride, or viton).
USER SAFETY RECOMMENDATIONS Users must wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Users must remove PPE

immediately after handling this product. Wash the outside of the gloves before removing. As soon as possible, wash thoroughly and change into clean clothing. Follow manufacturer's instructions for cleaning/maintaining PPE. If no such instructions for washables, use detergents and hot water. Keep and wash PPE separately from other laundry.

CLEANING PROCEDURE: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant. After cleaning is complete, follow disinfection directions described previously. Contact time: Leave surfaces wet for 10 minutes. Let air dry.

DISPOSAL OF INFECTIOUS MATERIALS: Contaminated waste and blood and other body fluids should be autoclaved and disposed according to local regulations for infectious waste disposal.

DECON-PHENE is a dilutable germicidal detergent formulated for use where broad spectrum disinfection is required in the presence of moderate amounts of organic soil. DECON-PHENE is effective in 10 minutes at 20° C, in the presence of 5% blood serum, used to simulate organic soil. DECON-PHENE'S multi-phenolic formula is designed to clean, disinfect and deodorize any washable, inanimate non-porous surface in one easy step. When used as directed, DECON-PHENE is strong enough to kill a broad spectrum of pathogenic bacteria, yet is mild enough to have no harmful effects on the surface being disinfected.

DECON-PHENE is acceptable for use by pharmaceutical and biotechnology industries, medical device manufacturers, hospitals and any health care institutions that are dedicated to controlling the hazards of cross contamination in their establishments. For use on hard, inanimate surfaces in aseptic filling and gowning rooms, general manufacturing areas, or on machinery, tables, counters, laminar flow benches, floors, walls, stainless steel, porcelain, glass and chrome.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with the labeling.

CLEANING AND DISINFECTION: For broad spectrum disinfection of gram negative and gram positive bacteria, including Pseudomonas aeruginosa, add 1 ounce DECON-PHENE to the gallon of water. Remove heavy soil or gross filth then apply solution with a mop, cloth, sponge or trigger sprayer, so as to wet the surface thoroughly. Allow to remain wet for 10 minutes and then let air dry. Where Tuberculocidal activity is required, clean the surface prior to application for the disinfectant solution. If higher detergency is desired, increase dilution to 2 ounces per gallon of water. Prepare a fresh solution for each use.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly to the human body, either into or in contact with blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used in to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage or disposal. **STORAGE:** Store in securely closed original container. Store in area that will prevent cross contamination of other pesticides, fertilizer, food and feed. Avoid storage at temperature extremes or in sunlight. **PESTICIDE DISPOSAL:** Pesticide wastes are acutely hazardous. Improper disposal of excess pesticides is in violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your Pesticide or Environmental Control Agency, or the hazardous waste representative at the nearest EPA Regional Office for guidance. **CONTAINER DISPOSAL:** Do not reuse container. Replace cap and discard in the trash.

MSDS #: DP-98-01

REV4/08 DP-02

P-02 Rev 4/08

Colors: Black

IMPORTANT

Please check ENTIRE proof carefully including ALL new, revised, repeat copy AND human readable barcode characters.
By approval of this proof you are accepting responsibility for correctness of ALL copy contained herein.
Refer to Pantone Color Chips for accurate PMS colors. All colors print 100% unless noted.

Proof # 1 Date: 05/01/08 JH/ak/jc

Ticket # 63949 Product # 36500

Customer: Veltek Associates, Inc.

Die Size: 5" x 11"
(.125" corner radius)

Colors: black

Alterations to supplied art and re-proofs will incur additional charges.

This proof is provided for your inspection of copy, layout, finished size and color breakdown. Please review carefully & indicate instructions below.

- Proof approved for printing as is. Proceed with order.
 Proof not approved. Make corrections as indicated & reproof.

Authorized Signature _____ Date _____

F-1042 Rev4 12/7/06

Attention: Colors may not be accurately represented due to monitor variation.

