



STATE OF HAWAII
Department of Agriculture

ACCEPTED
LICENSE NO. **8754.4**

Run Direction # 3

Dieline
Will Not Print.

DECON-SPORE 200 Plus®

Sterilization of Manufacturing,
Filling and Packaging Equipment
in Aseptic Processes.

Active Ingredients:

Hydrogen Peroxide.....	27.5%
Peroxyacetic Acid.....	5.8 %
Inert ingredients.....	66.7%
Total.....	100.0%

**KEEP OUT OF REACH
OF CHILDREN
DANGER**

FOR COMMERCIAL USE
STRONG OXIDIZING AGENT

Net Contents 12 x 13 ounces
EPA Reg. No. 1677-129-68959
EPA Est. 68959-PA-1



Container and Product Sterilized and Distributed By:
Veltek Associates, Inc.
15 Lee Blvd. Malvern, PA 19355-1234
Tel: 610-644-8335 • Fax: 610-644-8336

PRECAUTIONARY STATEMENTS - HAZARDS TO HUMANS AND DOMESTIC ANIMALS

DANGER: CORROSIVE: Causes irreversible eye damage and skin burns. May be fatal if inhaled. Harmful if swallowed. Do not get in eyes, on skin, or on clothing. Do not breathe vapor. Wash thoroughly after handling with soap and water and before eating, drinking or using tobacco. Remove contaminated clothing and wash clothing before reuse. The following Personal Protective Equipment (PPE) should be used when handling the product: coveralls over long-sleeved shirt and long pants, socks and chemical-resistant footwear, goggles or face shield, chemical-resistant gloves (such as rubber or made out of any waterproof material), chemical-resistant apron. Wear a mask or pesticide respirator jointly approved by the Mine Safety and Health Administration and the National Institute for Occupational Safety and Health.

FIRST AID

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call poison control center or doctor for treatment advice.

IF SWALLOWED: Call a poison control center or doctor immediately for treatment advice. Have a person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by the poison control center or doctor. Do not give anything by mouth to an unconscious person.

IF ON SKIN OR CLOTHING: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

IF INHALED: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably mouth to mouth, if possible. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or going for treatment.

For Spill/Exposure Emergency Response Service in English (and 26 languages) Call +44-208-762-8322, Arabic +961-3-487-287, Chinese +86-10-5100-3039 From USA and Canada prefix (+) with 011.

NOTE TO PHYSICIAN: Probable mucosal damage may contraindicate the use of gastric lavage.

PHYSICAL AND CHEMICAL HAZARDS: Strong oxidizing agent. Corrosive. Do not use in concentrated form. Mix only with water according to label instructions. Never bring concentrate in contact with other sanitizers, cleaners or organic substances.

ENVIRONMENTAL HAZARDS: This pesticide is toxic to birds, fish and aquatic invertebrates. Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance, contact your State Water Board or Regional Office of the EPA.

MSDS #: DS200-0397-01-01

REV1/10-1 DS200-03-13ZA

DIRECTIONS FOR USE:

It is a violation of Federal law to use this product in a manner inconsistent with its labeling. **STERILIZATION OF MANUFACTURING, FILLING AND PACKAGING EQUIPMENT IN ASEPTIC PROCESSES:**

Prior to the use of this product, remove gross soil particles from the processing surfaces, then wash with a recommended detergent solution, followed by a thorough potable water rinse. Prepare a sterilizing solution by diluting 6.4 ounces DECON-SPORE 200 Plus® concentrate per each gallon of water (50 ml/liter) (5.0% v/v).

Circulate, coarse spray, or flood the sterilizing solution through the system. All surfaces should be exposed to the sterilizing solution for a minimum exposure time based on the product solution temperature. The following time and temperature relationships are required:

DECON-SPORE 200 Plus® Concentration	Temperature	Time
5%	68° F (20° C)	6 hours
5%	122° F (50° C)	20 minutes
5%	176° F (80° C)	5 minutes

Thoroughly rinse food contact surfaces with either a sterile water or potable water rinse. For food contact surfaces, follow with a sanitizing solution of DECON-SPORE 200 Plus®. Allow surfaces to drain thoroughly prior to any food contact. This product is an effective sporicide against *Bacillus subtilis* and *Clostridium sporogenes* when used per the label directions. Kills Pandemic 2009 H1N1 influenza A virus (formerly called swine flu).

Note: This product in its use solutions is compatible with stainless steel and aluminum surfaces. If product is intended to be used on any other surface, it is recommended that you apply product to a smaller test area to determine compatibility before proceeding with its use.

STORAGE & DISPOSAL
DO NOT CONTAMINATE WATER, FOOD OR FEED BY STORAGE OR DISPOSAL.

PESTICIDE STORAGE: Product should be kept cool and in a vented container to avoid any explosion hazard.

PESTICIDE DISPOSAL: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.

CONTAINER DISPOSAL: Non-refillable container. Do not reuse or refill this container. Offer for recycling, if available. Triple rinse container (or equivalent) promptly after emptying.

Contents: 13 Ounces (384 ml) Quan/Case: 12x13 ounces Reorder #: DS200-03-13ZA Sterile

Description: DS200-03-13ZA Rev.1/10-1 (5" x 11")
Colors: Black, PMS 200

IMPORTANT

Please check ENTIRE proof carefully including ALL new, revised, repeat copy AND human readable barcode characters. By approval of this proof you are accepting responsibility for correctness of ALL copy contained herein. Refer to Pantone Color Chips for accurate PMS colors. All colors print 100% unless noted.

Proof # 8	Date: 3/11/10	JH/ak/jc	<p>Alterations to supplied art and re-proofs will incur additional charges.</p> <p>This proof is provided for your inspection of copy, layout, finished size and color breakdown. Please review carefully & indicate instructions below.</p> <p><input type="checkbox"/> Proof approved for printing as is. Proceed with order.</p> <p><input type="checkbox"/> Proof not approved. Make corrections as indicated & reproof.</p> <p>Authorized Signature _____ Date _____</p>
Customer: Veltek Associates, Inc.			
Ticket # 69100	Product # 38864		
Die Size: 5" x 11" (0.125" corner radius)			
Colors: Black, PMS 200			

F-1042 Rev 6/17/09

Colors may not be accurately represented due to monitor variation. Process color must be approved from a digital color proof & not from a pdf file. Files supplied as RGB format will be converted to CMYK & color shifts may occur.