

LICENSED

PERIOD **2014-2016** LIC. NO.

8822.1

LICENSED

PERIOD **2011-2013** LIC. NO.

From the makers of **ClōSYS**® Oral Care Products.

MicroCLEAR
DENTAL UNIT
**WATERLINE
CLEANER**

A CONCENTRATED SOLUTION FOR USE
IN SELF-CONTAINED WATER SYSTEMS

Cleans dental unit waterlines
Dental waterline microbial contaminant control solution
Kills odor causing bacteria
Maintains dental unit waterline effluent < 10 CFU/mL

Active Ingredients:
Chlorine dioxide.....0.1% wv
Inactive Ingredients:
Other Ingredients.....99.9% wv
Total.....100.0%

**1 Gallon
(3.785 L)**

CAUTION
KEEP OUT OF REACH
OF CHILDREN

Front

FIRST AID

IF IN EYES	<ul style="list-style-type: none"> - Hold eye open and rinse slowly and gently with water for 15 to 20 minutes. - Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. - Call a poison control center or doctor for treatment advice.
IF ON SKIN OR CLOTHING	<ul style="list-style-type: none"> - Take off contaminated clothing. - Rinse skin immediately with plenty of water for 15 to 20 minutes. - Call a poison control center or doctor for treatment advice.
IF SWALLOWED	<ul style="list-style-type: none"> - Call poison control center or doctor immediately for treatment advice. - Have person sip a glass of water if able to swallow. - Do not induce vomiting unless told to do so by the poison control center or doctor. - Do not give anything by mouth to an unconscious person.
IF INHALED	<ul style="list-style-type: none"> - Move person to fresh air. - If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible. - Call a poison control center or doctor for treatment advice.

For 24-hour emergency information on this product, call CHEMTREC 1-800-262-8200

Manufactured by:
Rowpar Pharmaceuticals, Inc. · 16100 North Greenway-Hayden Loop F-400 · Scottsdale, Arizona 85260

EPA Est. No. 65427-IL-001

EPA Reg. No. 87117-1

SEE OTHER PANEL FOR PRECAUTIONARY STATEMENTS

Side

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION: Causes moderate eye irritation. Harmful if absorbed through skin or if swallowed. Avoid contact with skin, eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

PHYSICAL OR CHEMICAL HAZARDS

Chlorine dioxide is a strong oxidizing agent. Contamination with other materials such as acids, chlorine, organic chemicals, etc. may cause a chemical reaction resulting in the evolution of gaseous chlorine dioxide and heat. Explosion and/or fire could result. Chlorine dioxide is a poisonous explosive gas. Keep all chemical and foreign materials away from this solution.

STORAGE AND DISPOSAL

Pesticide Storage: Store in a cool, dark area in original container. Do not contaminate water, food, or feed by storage or disposal. Do not store with easily oxidizable materials, acids, reducers, and combustible materials. **Pesticide Disposal:** Wastes resulting from the use of this product may be disposed of on-site or at an approved waste disposal facility. **Container Disposal:** Nonrefillable container. Do not reuse or refill this container. Clean container promptly after emptying. Offer for recycling, if available.

WARRANTY STATEMENT

The directions for use of this product are believed to be adequate and must be followed carefully, it is impossible to eliminate all risks inherently associated with the use of this product. To the extent consistent with applicable law, the products sold to you are furnished "as is" by Rowpar Pharmaceuticals, Inc. The manufacturer and the seller are subject only to the manufacturer's warranties, if any, which appear on the label of the product sold to you. Except as warranted by this label, Rowpar Pharmaceuticals, Inc. or the seller makes no warranties, guarantees, or representations of any kind to the buyer or the user, either express or implied, or by usage of trade, statutory or otherwise, with regard to the product sold or use of the product, including, but not limited to, merchantability, fitness for a particular purpose or use, or eligibility of the product for any particular trade usage. To the extent consistent with applicable law, Buyer's or user's exclusive remedy, and Rowpar Pharmaceuticals, Inc. or the seller's total liability shall be limited to damages not exceeding the cost of the product. No agent or employee of Rowpar Pharmaceuticals, Inc., or the seller is authorized to amend the terms of this warranty disclaimer or the product's label or to make a presentation or recommendation different from or inconsistent with the label of this product. To the extent consistent with applicable law, Rowpar Pharmaceuticals, Inc. or the seller shall not be liable for consequential, special, or indirect damages resulting from the use, handling, application, storage, or disposal of this product or for damages in the nature of penalties, and the buyer and the user waive any right that they may have to such damages.

MicroCLEAR®

NON-STAINING · NON-CORROSIVE · pH BALANCED · DOES NOT AFFECT BOND STRENGTH · DENTAL UNIT WATERLINE CLEANER · DENTAL WATERLINE MICROBIAL CONTAMINANT CONTROL SOLUTION KILLS ODOR CAUSING BACTERIA · MAINTAINS DENTAL UNIT WATER LINE EFFLUENT < 10 CFU/mL · USE MICROCLEAR DENTAL UNIT WATERLINE CLEANER TO CONTROL WATER QUALITY · DOES NOT DAMAGE OR HARM LINES OR EQUIPMENTS EVEN AFTER LONG-TERM USE

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Initial One-Time Treatment:

1. Empty and clean the self-contained bottle at the end of the day.
2. Fill the bottle with 100 mL (4 pumps) of MicroCLEAR®
3. Run each line for 30 seconds.
4. Leave MicroCLEAR® in the lines overnight.
5. The next morning completely flush remaining solution through waterlines and remove the self-contained bottle.
6. Fill clean bottle with the daily use formula. Run this solution through the lines for 30 seconds.

Daily Use:

1. Prepare mixture by adding the required volume of concentrate (see table below) directly to the system's self-contained water bottle every time it needs to be refilled.
2. Fill the remainder of the bottle with water.

Volume of Bottle (mL)	Amount of MicroCLEAR		Number of Pumps
	(mL)	(oz)	
*750 (A-dec)	75	2.5	3
1000	100	3.5	4
*2000	200	7	8

*Most common volumes on self-contained units.

Note:

One full pump will dispense ¾ ounces (25 mL).
Use a 1:10 dilution if the volume needed is not shown.

For further information, call Rowpar Pharmaceuticals at 800-643-3337.

Back

Side