



Department of Agriculture
STATE OF HAWAII

LICENSED

PERIOD **2015-2017** LIC. NO.

8884.1



peridox[®]

Concentrate



Department of Agriculture
STATE OF HAWAII

LICENSED

PERIOD **2012-2014** LIC. NO.

Disinfectant • Virucide • Fungicide • Tuberculocide • Bactericidal • Sporicidal Disinfectant

Kills 99.9999% *Clostridium difficile* spores in 3 minutes
Disinfects against Acinetobacter baumannii in 2 minutes
Disinfects against MRSA in 2 minutes
Disinfects against VRE in 2 minutes

- Concentrate, Dilute Before Use**
- Cleans, disinfects, and deodorizes**
- Disinfects hard, nonporous surfaces in Hospitals and Health Care Facilities**
- Prevents mold, mildew stains and odors**
- Does NOT contain Bleach or Alcohol**
- Will not harm LCD screens**

ACTIVE INGREDIENTS:

Hydrogen Peroxide.....	24.0%
Peroxyacetic Acid.....	1.2%
Other Ingredients.....	74.8%
Total.....	100.0%

KEEP OUT OF REACH OF CHILDREN
DANGER PELIGRO

Si usted no entiende la etiqueta, busque a alguien para que se la explique a usted en detalle.
(If you do not understand the label, find someone to explain it to you in detail.)

See back panel for
additional precautionary
and first aid statements.

PN 810000015F_0

NET CONTENTS: 1 U.S. GALLON (3.8 Liters)



**PRECAUTIONARY STATEMENTS
HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

DANGER. CORROSIVE. Causes irreversible eye damage and skin burns. May be fatal if absorbed through the skin or inhaled. Harmful if swallowed. Do not get in eyes, on skin or on clothing. Do not breathe vapor or spray mist. Wear protective eyewear (goggles, face shield, or safety glasses), clothing, and chemical resistant gloves. Wash with soap and water thoroughly after handling and before eating, drinking or using tobacco. Remove contaminated clothing and wash clothing before reuse. Dilute the concentrated Peridox® to the appropriate use rate in a well-ventilated area or have available proper respiratory protection. **PHYSICAL AND CHEMICAL HAZARDS:** Strong Oxidizing Agent. Corrosive. Do not use in concentrated form. Mix only with water according to label instructions. Never bring concentrate in contact with other sanitizers, cleaners or organic substances. **ENVIRONMENTAL HAZARDS:** This pesticide is toxic to birds, fish and aquatic invertebrates. Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office at the EPA.

FIRST AID: IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice. **IF SWALLOWED:** Call a poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to by a poison control center or doctor. Do not give anything to an unconscious person. **IF ON SKIN OR CLOTHING:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or physician for treatment advice. **IF INHALED:** Move person to fresh air. If person is not breathing, call 911 or ambulance, and then give artificial respiration, preferably mouth-to-mouth, if possible. Call a poison control center or physician for further treatment advice. Have the product container with you when calling a poison control center or doctor or going for treatment. **FOR EMERGENCY MEDICAL INFORMATION CALL TOLL FREE: 1-800-222-1222.** NOTE TO PHYSICIAN: Probably mucosal damage may contraindicate the use of gastric lavage.

STORAGE AND DISPOSAL: Do not contaminate water, food or feed by storage or disposal. **PESTICIDE STORAGE:** Store in original container in a cool, well-ventilated area, away from direct sunlight. Do not allow product to become overheated in storage. This may cause increased degradation of the product, which will decrease product effectiveness. In case of spill, flood area with large quantities of water. **PESTICIDE DISPOSAL:** Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance. If material has spilled, an acceptable method of disposal is to dilute with at least 20 volumes of water followed by discharge into suitable treatment system in accordance with all local, state, and Federal environmental laws, rules, regulations, standards and other requirements. **CONTAINER DISPOSAL:** If empty, wrap container and put in trash or offer for recycling, if available. If partly filled, call your local solid waste agency for disposal instructions.

Peridox® is a registered trademark of BioMed Protect, St. Louis, MO, USA. Patents Pending

EPA Reg. No. 88089-2
EPA Est. No. 88089-MO-1^A, 87720-MO-1^B
67811-MO-1^C, 4313-TX-1^D

Superscript is first letter of lot number
See Lot No. on Bottle

For Emergency Medical Information: 1-800-222-1222
For Product Information: 1-800-691-7150

Product No. 85342



PN 810000015R_0

Manufactured for:



13475 LAKEFRONT DR, ST. LOUIS, MO 63045

file: 61639_book_art
folder: 61639_biomed_810000015R
colors: black, 381, 5405
barcode: upca, 80%, -1.54

created by: 02-11-13 ek
alt: 03-21-13 ek
size: 4.625" (w) x 5.5" (h)
 dp proofed

WARRANTY

This material is safe and effective under the conditions of application as specified in the Directions for Use. To the fullest extent permitted by law, the manufacturer shall not be liable for improper application or use of this product beyond those purposes referred to in the Directions for Use. NO OTHER EXPRESS OR IMPLIED WARRANTY OF FITNESS OR MERCHANTABILITY IS MADE.

Effective against a wide range of viruses, bacteria and fungi, including:

Bactericidal (2 minutes): *Enterobacter aerogenes* (ATCC 13048); *Listeria monocytogenes* (ATCC 19111); *Pseudomonas aeruginosa* (ATCC 15442); *Salmonella typhimurium* (ATCC 13311); *Staphylococcus aureus* (ATCC 6538); *Vibrio cholerae* (ATCC 14035); *Salmonella enterica* (ATCC 10708); *Acinetobacter baumannii* (ATCC 19606); *Campylobacter jejuni* (ATCC 29428); *E. coli* O157:H7 (ATCC 35150); *E. coli* ESBL (ATCC BAA-196); *Enterococcus faecalis*-Vancomycin Resistant (VRE) (ATCC 51299); *Enterococcus hirae* (ATCC 10541); *Haemophilus influenzae* (ATCC 19418); *Klebsiella pneumoniae* (ATCC 4352); *Legionella pneumophila* (ATCC 33153); *Proteus vulgaris* (ATCC 27973); *Serratia marcescens* (ATCC 13880); *Shigella sonnei* (ATCC 11060); *Staphylococcus aureus* – MRSA (ATCC 33592); Community Acquired - Methicillin Resistant *Staphylococcus aureus* (CA-MRSA) (Genotype 300), *Streptococcus pneumoniae* - Penicillin Resistant (PRSP) (ATCC 700671), *Streptococcus pyogenes* (ATCC 19615)

Tuberculocidal (5 minutes): *Mycobacterium bovis* (at 20°C)

Virucidal (2 minutes): *Hepatitis B Virus*; *Avian Influenza A (H3N2)* (ATCC VR 2072); *Reovirus* (ATCC VR 232); *Human coronavirus* (ATCC VR 740); *Adenovirus, Type 2* (ATCC 846); *Avian influenza H5N1* (NIBRG-14); *Herpes simplex, Type 1* (ATCC VR 260); *Herpes simplex, Type 2* (ATCC VR 734); *Human immunodeficiency virus Type 1* (HIV-1, Zeptomatrix Corp.); *Influenza A* (ATCC VR 544); *Influenza A (H1N1)* (ATCC VR 1469, Strain A/PR/8/34); *Influenza B* (ATCC VR 823); *Norovirus (F. Calicivirus)* (ATCC VR 782); *Poliovirus Type 1* (ATCC VR 1562); *Respiratory syncytial virus* (RSV) (ATCC VR 26); *Rhinovirus* (ATCC VR 338); *Rotavirus* (Strain WA)

Sporicidal Activity of Disinfectant (3 minutes): *Clostridium difficile* spores (ATCC 700792)

Fungicidal (1 minute hard surfaces): *Trichophyton mentagrophytes* (athlete's foot fungus) (ATCC 9533); *Candida albicans* (ATCC 10231); *Aspergillus niger* (ATCC 16404)

Disinfectant for use on hard, non-porous surfaces in Hospitals and Health Care Facilities such as Clinics, Dental Offices, Hospices, Laboratories, Nursing Homes, Acute Long-Term Care and Long-Term Care Facilities, Physical Therapy, Doctor's Offices, Radiology, Rehabilitation and Transport Vehicles; Critical Care Areas such as Critical Care Unit (CCU), Emergency Room, ICU, Neonatal Intensive Care Unit (NICU), Operating Rooms, Outpatient Surgery Centers, Pediatric Intensive Care Unit (PICU), Surgery, and Emergency Medical Services (EMS); Other Sites such as Schools, Colleges, Day-care Centers, Child-care Facilities, Correctional Facilities; Hospitality Establishments such as Hotels, Motels, and Cruise Ships; Veterinary Clinics, Animal Science Laboratories; Funeral Homes and Morgues; Research Laboratories; Industrial Facilities; Pharmaceutical Production Facilities; Cosmetic Processing Facilities; Consumer Home Use; Dining Areas (non-food contact surfaces); Office Buildings; Recreational Facilities; Retail and Wholesale Establishments; Prisons; Animal Care Facilities; Veterinary Facilities; Farms; Livestock Quarters; Poultry Premises; Poultry Houses and Hatcheries; and vehicles including ships, planes, automobiles, trucks, ambulances, trains, and farm vehicles.

Hospital, Healthcare, and Critical Use Sites: May be used on hard, non-porous surfaces of autoclaves, bed railings, blood glucose monitors, cabinets, carts, chairs, counters, exam tables, gurneys, isolettes, infant incubators and care cribs, IV poles, phlebotomy trays, PVC tubing, stethoscopes, stretchers, tables, bathrooms, sinks, faucets, toilets seats and rims, towel dispensers, hand railings, stall doors, bath tubs, showers, hampers, tiled walls, telephones, door knobs, vanities, floors, non-porous shelves, and display cases.

Hard, non-porous surfaces of ambulance equipment, diagnostic equipment, dialysis machines, mammography equipment, patient monitoring equipment, respiratory equipment, ultrasound transducers and probes, patient support and delivery equipment. Peridox® solution will not harm aluminum, LDPE, HDPE, vinyl, painted surfaces, polycarbonate, polypropylene, polyurethane varnish, PVC, silicon rubber, stainless steel, medical tubing, vinyl rubber, acrylic, brass, LCD screens, copper, Corian®, glass, laminate flooring and countertops, nickel, glazed porcelain, glass, and glazed tile.

DIRECTION FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Peridox® can be used for disinfection. This product must be used only as directed in this label.

Before using this product read the entire label.

Peridox® solution can be applied to hard, non-porous surfaces such as counter-tops, tables, furniture, cabinets, appliances, kitchen and bathroom fixtures, furnishings, hardware, floors, walls, equipment, tabletop, benchtops, drip pans, troughs, desks, locker room benches, and the use sites listed. For heavily soiled areas, pre-cleaning is required. To apply the Peridox® solution, the measured amount of Peridox® per the Directions for Use is added to a spray applicator or mixing container and diluted with the measured amount of water (up to 250 ppm hardness as CaCO₃). The "Use Dilutions" for Peridox® are based on the 24.0% hydrogen peroxide active ingredient. Apply the mixed Peridox® solution to the contaminated surface with a spray applicator device at an application rate of greater than 50 ml/m² (approximately 1.5 fluid oz. of the Peridox® solution per square yard), or by other means that ensure that the surface is wet (for example, by wiping with a clean cloth or sponge that is wet with Peridox® solution), or by soaking the contaminated surface and allowing it to remain on the surface the indicated contact time. Spray or soak the surface with a sufficient amount of Peridox® solution to remain on surface for the indicated contact time before evaporation. Air dry or remove with a dry clean cloth or wipe.

DISINFECTION

Use Peridox® solution as a disinfectant in Hospitals, Obstetrical Suites, Housekeeping Services, Physical Therapy Departments, Nursing Services, Autopsy Facilities, nursing homes, health-care facilities, schools, colleges, day-care centers, child-care facilities, veterinary clinics, animal life science laboratories, industrial facilities, dining areas (non-food contact surfaces), office buildings, recreational facilities, retail and wholesale establishments, and in transport vehicles including ships, planes, automobiles, trucks, trains, and EMS.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

Combination Disinfection and Cleaning: Peridox® disinfectant is effective against *Staphylococcus aureus*, *Salmonella enterica*, *Enterobacter aerogenes*, *Salmonella typhimurium*, *Listeria monocytogenes*, *Vibrio cholerae*, *Pseudomonas aeruginosa*, and other bacteria listed in the Chart at 4% (26 oz. Peridox®/gallon water) in hard water (250 ppm as CaCO₃) and a 5% organic load on hard, non-porous surfaces. For heavily soiled areas a pre-cleaning step is required. Apply Peridox® solution with spray applicator device or other means that ensures that the surface is thoroughly wet. Allow surfaces to remain wet for **2 minutes**; then, wipe with a dry clean cloth. Prepare a fresh solution daily or, if necessary, more frequently if the Peridox® solution becomes soiled (i.e., is not visibly clear) or diluted, replace and use fresh Peridox® solution. Always ensure the solution is fresh prior to application.

Virucidal: Peridox® disinfectant is effective against Influenza A (H3N2), Influenza A (H1N1), Norovirus, Poliovirus, Reovirus, Rhinovirus, Human Coronavirus and specific viruses listed on the Chart on hard, non-porous inanimate surfaces when used at 4% (26 oz. Peridox®/gallon water) in hard water (250 ppm CaCO₃) with a **2 minute** contact time in the presence of a 5% organic load. For heavily soiled areas a pre-cleaning step is required. Apply Peridox® solution with a spray applicator device or other means that ensures that the surface is thoroughly wet. Allow surfaces to remain wet for **2 minutes**; then, wipe with a dry, clean cloth. Prepare a fresh solution daily or, if necessary, more frequently if it becomes soiled (i.e., is not visibly clear) or diluted (i.e., the solution becomes further diluted below the 4% Use Dilution concentration). Always ensure the solution is fresh prior to application.

Fungicidal: Peridox® disinfectant can be used on hard, non-porous inanimate surfaces such as shower room floors, shower stalls, locker room benches, and hard, non-porous bath mats (plastic or rubber). At 4% (26 oz. Peridox®/gallon water), this product is effective against *Trichophyton mentagrophytes* (Athletes Foot Fungi), *Aspergillus niger* (mold), and *Candida albicans* (pathogenic yeast) in the presence of a 5% organic load in hard water (250 ppm CaCO₃) with a **1 minute** contact time. For heavily soiled areas a pre-cleaning step is required. Apply Peridox® solution with spray applicator device or other means that ensures that the surface is thoroughly wet. Allow surfaces to remain wet for **1 minute**; then, wipe with a dry, clean cloth. Prepare a fresh solution daily or, if necessary, more frequently if it becomes soiled (i.e., is not visibly clear) or diluted (i.e., the solution becomes further diluted below the 4% Use Dilution concentration). Always ensure the solution is fresh prior to application.

Sporicidal Disinfectant: Peridox® solution is effective against *Clostridium difficile* spores on hard, non-porous inanimate surfaces when used at full strength (4.4% concentration of hydrogen peroxide and 0.23% of peroxyacetic acid) with a **3 minute** contact time in the presence of a 5% organic load. Apply PERIDOX® solution with a spray applicator device, or other means that ensures that the surface is thoroughly wet. Allow surfaces to remain wet for **3 minutes**; then, wipe with a dry, clean cloth.

Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks and eye covering.

Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a microfiber cloth/pad, clean cloth, mop, and/or sponge saturated with the disinfectant product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.

Infectious Materials Disposal: Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

Specific Instructions for Blood-borne Pathogens: KILLS HIV AND HEPATITIS B VIRUS ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in healthcare settings and other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, on surfaces or materials soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) and *Hepatitis B Virus*.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING AGAINST HIV-1 (AIDS Virus) AND HEPATITIS B VIRUS OF SURFACE/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Disposable latex, nitrile, or vinyl gloves, gowns, face masks, or eye coverings as appropriate, must be worn during all cleaning of body fluids, blood, and decontamination.

Cleaning Procedure: Blood/body fluids must be thoroughly cleaned from surfaces/objects before application of this product.

Contact Time: Allow surface to remain wet for **2 minutes** to kill HIV-1 and *Hepatitis B Virus*. Use contact time listed for all other organisms claimed.

Disposal of Infectious Material: Blood/body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

Tuberculocidal: Peridox® microbiocide is also effective against *Mycobacterium bovis* (at 20°C) on hard, non-porous inanimate surfaces when used at 4% (26 oz. Peridox®/gallon water) with a **5 minute** contact time in the presence of hard water (250 ppm CaCO₃) and a 5% organic load. For heavily soiled areas a pre-cleaning step is required. Apply Peridox® solution with spray applicator device or other means that ensures that the surface is thoroughly wet. Allow surfaces to remain wet for **5 minutes**; then, wipe with a dry clean cloth. Prepare a fresh solution daily or, if necessary, more frequently if it becomes soiled (i.e., is not visibly clear) or diluted (i.e., the solution becomes further diluted below the 4% Use Dilution concentration).

Disinfection of Surfaces in Pharmaceutical and Cosmetic Processing Facilities: Peridox® microbiocide is recommended for use on hard, non-porous environmental surfaces such as floors, walls and processing equipment in pharmaceutical and cosmetic processing facilities. It is effective against *Staphylococcus aureus*, *Salmonella enterica*, *Enterobacter aerogenes*, *Salmonella typhimurium*, *Listeria monocytogenes*, *Vibrio cholerae*, and *Pseudomonas aeruginosa* at 4% (26 oz. Peridox®/gallon water) in hard water (250 ppm CaCO₃), and 5% blood serum. For heavily soiled areas a pre-cleaning step is required. Apply Peridox® solution with spray applicator device or other means that ensures that the surface is thoroughly wet with a contact time of **2 minutes**. Rinse with sterile water. Wipe with a dry, clean cloth. Air Dry. Prepare a fresh solution daily, or if necessary, more frequently if it becomes soiled (i.e., is not visibly clear) or diluted (i.e., the solution becomes further diluted below the 4% Use Dilution concentration). Always ensure the solution is fresh prior to application.

Disinfectant for Farm Premises: Peridox® microbiocide is recommended for use on hard, non-porous, environmental surfaces such as floors, walls, troughs, racks, appliances, and processing equipment on Farm Premises. It is effective against *Staphylococcus aureus*, *Salmonella enterica*, *Enterobacter aerogenes*, *Salmonella typhimurium*, *Listeria monocytogenes*, *Vibrio cholerae*, and *Pseudomonas aeruginosa* and other bacteria listed in the Chart at 4% (26 oz. Peridox®/gallon water) in hard water (250 ppm CaCO₃), and 5% blood serum. For heavily soiled areas, a pre-cleaning step is required. Always ensure the solution is fresh prior to application.

Remove all animals and feed from premises, vehicles, and enclosures. Remove all litter and manure from floors, walls and surfaces of barns, pens, stalls, chutes, and other facilities and fixtures occupied or traversed by animals. Empty all troughs, racks, and other feeding and watering appliances. Thoroughly clean all surfaces with soap or detergent and rinse with water. Apply Peridox® solution with a spray applicator device, or other means that ensures that the surface is thoroughly wet with a contact time of **2 minutes**. Rinse with clean water. Wipe with a dry, clean cloth. Air dry.

Immerse all halters, ropes, and other types of equipment used in handling and restraining animals, as well as forks, shovels, and scrapers used for removing litter and manure. If the Peridox® solution becomes soiled (i.e., is not visibly clear) or diluted, replace and use fresh Peridox® solution.

Ventilate buildings, cars, boats, and other closed spaces. Do not house livestock or employ equipment until treatment has been absorbed, set, or dried. Thoroughly scrub all treated feed racks, mangers, troughs, automatic feeders, fountains, and waterers with soap or detergent, and rinse with potable water before reuse.

Disinfectant for Poultry Houses and Poultry Hatcheries: Peridox® microbiocide is recommended for use on hard, non-porous, environmental surfaces such as floors, walls, troughs, racks, feeders, fountains, appliances, and processing equipment in Poultry Houses and Poultry Hatcheries. It is effective against *Staphylococcus aureus*, *Salmonella enterica*, *Enterobacter aerogenes*, *Salmonella typhimurium*, *Listeria monocytogenes*, *Vibrio cholerae*, and *Pseudomonas aeruginosa* and other bacteria listed in the Chart at 4% (26 oz. Peridox®/gallon water) in hard water (250 ppm CaCO₃), and 5% blood serum. For heavily soiled areas a pre-cleaning step is required. Always ensure the solution is fresh prior to application.

Remove all poultry and feeds from premises, trucks, coops, and crates. Remove all litter and droppings from floors, walls and surfaces of facilities occupied or traversed by poultry. Empty all troughs, racks, and other feeding and watering appliances. Thoroughly clean all surfaces with soap or detergent and rinse with water. Apply Peridox® solution with a spray applicator device, or other means, that ensures that the surface is thoroughly wet with a contact time of **2 minutes**. Rinse with clean water. Wipe with a dry, clean cloth. Air dry. Prepare a fresh solution daily or, if necessary, more frequently if it becomes soiled (i.e., is not visibly clear) or diluted (i.e., the solution becomes further diluted below the 4% Use Dilution concentration).

Ventilate buildings, coops, and other closed spaces. Do not house poultry or employ equipment until treatment has been absorbed, set, or dried. Thoroughly scrub treated feed racks, troughs, automatic feeders, fountains, and waterers with soap or detergent, and rinse with potable water before reuse.

NOTES

NOTES

NOTES



PRECAUTIONARY STATEMENTS HAZARDS TO HUMANS AND DOMESTIC ANIMALS

DANGER. CORROSIVE. Causes irreversible eye damage and skin burns. May be fatal if absorbed through the skin or inhaled. Harmful if swallowed. Do not get in eyes, on skin or on clothing. Do not breathe vapor or spray mist. Wear protective eyewear (goggles, face shield, or safety glasses), clothing, and chemical resistant gloves. Wash with soap and water thoroughly after handling and before eating, drinking or using tobacco. Remove contaminated clothing and wash clothing before reuse. Dilute the concentrated Peridox® to the appropriate use rate in a well-ventilated area or have available proper respiratory protection.

PHYSICAL AND CHEMICAL HAZARDS: Strong Oxidizing Agent. Corrosive. Do not use in concentrated form. Mix only with water according to label instructions. Never bring concentrate in contact with other sanitizers, cleaners or organic substances.

ENVIRONMENTAL HAZARDS: This pesticide is toxic to birds, fish and aquatic invertebrates. Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office at the EPA.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal

PESTICIDE STORAGE: Store in original container in a cool, well-ventilated area, away from direct sunlight. Do not allow product to become overheated in storage. This may cause increased degradation of the product, which will decrease product effectiveness. In case of spill, flood area with large quantities of water.

PESTICIDE DISPOSAL: Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance. If material has spilled, an acceptable method of disposal is to dilute with at least 20 volumes of water followed by discharge into suitable treatment system in accordance with all local, state, and Federal environmental laws, rules, regulations, standards and other requirements.

CONTAINER DISPOSAL: If empty, wrap container and put in trash or offer for recycling, if available. If partly filled, call your local solid waste agency for disposal instructions.

WARRANTY

This material is safe and effective under the conditions of application as specified in the Directions for Use. To the fullest extent permitted by law, the manufacturer shall not be liable for improper application or use of this product beyond those purposes referred to in the Directions for Use. NO OTHER EXPRESS OR IMPLIED WARRANTY OF FITNESS OR MERCHANTABILITY IS MADE.

Product No. 85342

Manufactured for: BioMed Protect, LLC
13475 Lakefront Dr.
St. Louis, MO 63045

EPA Reg. No. 88089-2
EPA Est. No. 88089-MO-1^A, 87720-MO-1^B
67811-MO-1^C, 4313-TX-1^D
Superscript is first letter of lot number
See Lot No. on Bottle