

Maximum Image Area Template
for Screen Print Decorating
on Priority Plastics'
76oz/120mm Canister
7/13/09

front panel image area

back panel image area



LICENSED

PERIOD 2015-2017 LIC. NO.

9021.19



SANICARE

Disinfecting Wipes

One Step Cleaner/Disinfectant/Deodorizer/Sanitizer

Active Ingredients:

Octyl decyl dimethyl ammonium chloride.....	0.069%
Dioctyl dimethyl ammonium chloride.....	0.028%
Didecyl dimethyl ammonium chloride.....	0.042%
Alkyl (C ₁₄ , 50%; C ₁₂ , 40%; C ₁₆ , 10%) dimethyl benzyl ammonium chloride.....	0.093%
Other Ingredients.....	99.768%
Total.....	100.000%

¹Includes weight of absorbent wipe

**KEEP OUT OF REACH OF CHILDREN
CAUTION**

See back panel for Precautionary Statements and First Aid
EPA Reg. No. 6836-313-559 EPA Est. No. 70930-PA-001

Net Contents: 120 disinfecting wipes
Net Contents: 2.0547 lbs (932 g)
7 in. x 8 in. (17.8 cm x 20.3 cm)

Product # 50960200

Manufactured for:

Buckeye International, Inc.

2700 Wagner Place • Maryland Heights, MO 63043 • USA
314.291.1900 • www.buckeyeinternational.com

Buckeye Sanicare Disinfecting Wipes help you accomplish things you never thought a wipe could do. Cleans, disinfects, and deodorizes in one labor saving step. Easy to use – no sprays, no sponges, no mess. Use **Buckeye Sanicare Disinfecting Wipes** on washable hard, nonporous surfaces of Athletic Equipment, Bathroom fixtures, Countertops, Desks, Doorknobs, Light switches, Tanning Beds, Toilet Seats, Tubs and Glazed Tiles. Cross-contamination is of major concern. **Buckeye Sanicare Disinfecting Wipes** have been formulated to aid in the reduction of cross-contamination on treated surfaces in schools, institutions and industry. Just pull, wipe and toss. Tough cleaning made easy.

Kills Pandemic 2009 H1N1 influenza A virus (formerly called swine flu).

Kills 99.9% of bacteria² in just 15 seconds

²Staphylococcus aureus (Staph), Klebsiella pneumoniae [Klebsiella], Campylobacter jejuni

Buckeye Sanicare Disinfecting Wipes are effective against:

Bacteria:

- Staphylococcus aureus (Staph)
- Salmonella enterica (Salmonella)
- Pseudomonas aeruginosa [Pseudomonas]
- Campylobacter jejuni [Campylobacter]
- Acinetobacter baumannii
- Burkholderia cepacia
- Corynebacterium ammoniagenes
- Enterobacter aerogenes
- Enterococcus faecalis [Enterococcus]
- Enterococcus faecalis – Vancomycin Resistant [VRE]
- Escherichia coli [E. coli]
- ESBL Escherichia coli – [Extended spectrum beta-lactamase producing E. coli]
- Klebsiella pneumoniae [Klebsiella]
- Legionella pneumophila
- Listeria monocytogenes
- Salmonella schottmuelleri [Salmonella]
- Salmonella typhi [Salmonella]
- Serratia marcescens [Serratia]
- Shigella dysenteriae [Shigella]
- Staphylococcus aureus – Community Associated Methillin-Resistant [CA-MRSA] [NRS384] [USA 300]
- Staphylococcus aureus – Community Associated Methillin-Resistant [CA-MRSA] [NRS123] [USA 400]

- Staphylococcus aureus – Methicillin-Resistant [MRSA]
- Staphylococcus aureus – Multi-Drug Resistant
- Staphylococcus aureus-Vancomycin Intermediate Resistant – [VISA]
- Streptococcus pyogenes [Strep] [a cause of scarlet fever]
- Vibrio cholerae

Viruses:

- *Hepatitis B Virus [HBV]
- *Hepatitis C Virus [HCV]
- *Herpes Simplex Virus Type 1
- *Herpes Simplex Virus Type 2
- *HIV-1 [AIDS]
- *Human Coronavirus
- *Influenza Type A Virus/Brazil [Influenza] [the Flu Virus]
- *Norwalk Virus [Norovirus]
- *Respiratory Syncytial Virus [RSV]
- *Rotavirus [Diarrhea-causing virus]
- *SARS Associated Coronavirus [SARS] [cause of Severe Acute Respiratory Syndrome]
- *Vaccinia Virus [Pox Virus]

Fungi:

- Aspergillus niger
- Trichophyton mentagrophytes [the Athlete's Foot Fungus] [A cause of Ringworm of the foot]

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

TO CLEAN AND DEODORIZE: Wipe surface with wipe. Let air dry. Discard used wipe in trash. Do not flush in toilet.

TO CLEAN AND DISINFECT / VIRUCIDAL* / MOLD AND MILDEW CONTROL DIRECTIONS: Wipe hard, nonporous surface with wipe until surface is visibly wet. Allow surface to remain wet for 10 minutes. Let air dry. If surfaces are extremely dirty, clean first with another wipe before disinfecting. Discard used wipe in trash. Do not flush in toilet. Rinse food contact surfaces with clean, potable water after disinfecting. Do not use to disinfect dishes, glassware, or utensils. Do not use as a diaper wipe or for personal cleansing. Note: Not recommended for unpainted wood, natural marble or brass.

For fungicidal activity, Buckeye Sanicare Disinfecting Wipes is an effective fungicide against Trichophyton mentagrophytes [the athlete's foot fungus] in areas such as locker rooms, dressing rooms, shower and bath areas and exercise facilities. **Buckeye Sanicare Disinfecting Wipes** is a one-step fungicide.

For mold and mildew, Buckeye Sanicare Disinfecting Wipes will effectively inhibit the growth of mold and mildew and the odors caused by them when applied to hard, nonporous surfaces. Follow disinfection instructions. Repeat treatment every seven days, or more often if new growth appears.

***KILLS HIV AND HBV AND HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) or Hepatitis B Virus (HBV) or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 OR HBV OR HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS. Personal protection: Clean-up must always be done wearing protective gloves, gowns, masks, and eye protection. **Cleaning procedures:** Blood and other body fluids containing HIV or HBV or HCV must be thoroughly cleaned from surfaces and objects before application of **Buckeye Sanicare Disinfecting Wipes**. Wipe can be used for this purpose. **Contact time:** Leave surface wet for 1 minute (60 seconds) for HIV-1 and 10 minutes for HBV and HCV. Use a 10-minute contact time for disinfection against all other viruses, bacteria and fungi claimed. **Disposal of infectious materials:** Blood, body fluids, cleaning materials and clothing must be autoclaved and disposed of according to local regulations for infectious waste disposal. **This product is for use on hard nonporous surfaces not conducive to treatment by immersion or excess liquid against Avian Influenza A.**

SANITIZING:

To sanitize hard, non-porous non-food contact surfaces; wipe surface, using enough wipes to allow treated surface to remain wet for 15 seconds for Staphylococcus aureus, Klebsiella pneumoniae, and Campylobacter jejuni. Let surface dry. For highly soiled surfaces, clean excess dirt first.

For plastic and painted surfaces, spot test on an inconspicuous area before use. This product will not harm most surfaces: including acrylic, sealed fiberglass and vinyl. Do not use on unpainted wood. Not for personal cleansing. This is not a baby wipe! Do not flush in toilet.

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

Caution. Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

FIRST AID

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling poison control center or doctor or going for treatment. For emergency information call 800-303-0441. Outside North America call 651-632-8956.

STORAGE AND DISPOSAL

Store in original container in areas inaccessible to children. Keep securely closed. Nonrefillable container. Do not reuse or refill this container. Discard empty container in trash or recycle if available.

To Open Package:




Remove cover. Pull up wipe corner from center of roll. Thread through hole in the cover. Replace cover. Open center cap on cover to remove wipe. When not in use, close center cap to prevent moisture loss.

Manufactured for:

Buckeye International, Inc.

2700 Wagner Place • Maryland Heights, MO 63043 • USA
314.291.1900 • www.buckeyeinternational.com



2404 (3/15)

PL

CL

PL

CL

PL

No print zone

No print zone

No print zone

5.312" high

IN/OUT.25"

7.28" deco area

.25"

7.28" deco area

.25" IN/OUT

15.315" circumference

.813" Minimum Baseline (lowest pass)