

**Healthcare Facilities, Academic Facilities, Dietary Areas, Office Buildings, Recreational Facilities, Retail and Wholesale Establishments, Institutional and Industrial Use**

*OxyCide Daily Disinfectant Cleaner* use solution cleans, disinfects and deodorizes hard, nonporous hospital/medical surfaces in one step with no rinsing required. This product evaporates completely and leaves no visible residue. This product is specifically designed to disinfect, deodorize and clean inanimate hard non-porous surfaces such as walls, floors, sink tops, furniture, patient beds, and operating tables. May be used to clean and disinfect finished floors. Cleans and disinfects without dulling gloss. May be used to clean and disinfect floor areas, sinks, faucets, bathrooms and tubs. In addition this product will deodorize those areas that are generally hard to keep fresh smelling such as garbage storage areas, empty garbage bins and cans, and any other areas that are prone to odors caused by microorganisms. Do not use on marble or un-sealed/un-coated terrazzo floors. This product use solution is a broad spectrum disinfectant cleaner and deodorizer designed for general cleaning and disinfecting of hard, nonporous surfaces. Quickly removes dirt, grime, blood and other organic matter commonly found in hospitals on medical surfaces. Use where housekeeping is of prime importance in controlling the hazard of cross contamination.

Using AOAC approved test methods (under Good Laboratory Practices), in the presence of 5% blood serum and hard water up to 400 ppm hardness (calculated as CaCO3), this product kills the following organisms on hard, nonporous inanimate surfaces.

<b>3 MINUTE CONTACT TIME:</b> <i>Pseudomonas aeruginosa</i> ; ATCC 15442 <i>Salmonella enterica</i> ; ATCC 10708 <i>Staphylococcus aureus</i> ; ATCC 6538 <i>Staphylococcus aureus</i> , (MRSA); ATCC 33592 <i>Staphylococcus aureus</i> , (Genotype USA 300) (CA-MRSA); ATCC BAA-1556 <i>Staphylococcus aureus</i> , (Genotype USA400) (CA-MRSA); ATCC BAA-1683 <i>Staphylococcus aureus</i> , (VISA); ATCC 700788 <i>Enterococcus faecalis</i> , (VRE); ATCC 51299 <i>Staphylococcus epidermidis</i> , (MRSE); ATCC 51625	<i>Streptococcus pneumoniae</i> ; ATCC 6303 <i>Streptococcus pyogenes</i> ; ATCC 19615 <i>Bordetella pertussis</i> , (Whooping Cough); ATCC 12743 <i>Escherichia coli</i> ; ATCC 11229 <i>Escherichia coli</i> , (ESBL); ATCC BAA-196 <i>Klebsiella pneumoniae</i> ; ATCC 4352 <i>Klebsiella pneumoniae</i> , (KPC); ATCC BAA-1705 <i>Acinetobacter baumannii</i> ; ATCC 19606 <i>Acinetobacter baumannii</i> , (Multi-drug Resistant - gentamicin, imipenem, ceftazidime); ATCC BAA-1605 <i>Proteus mirabilis</i> ; ATCC 7002 Human Immunodeficiency Virus Type 1 (AIDS virus)	Herpes Simplex Type I virus (F strain) ATCC VR-733 Herpes Simplex Type II virus (G strain); ATCC VR-734 Influenza A virus (Strain Hong Kong); ATCC VR-544 Respiratory Syncytial Virus, (RSV); ATCC VR-26 Vaccinia Virus, (Strain WR) (Pox Virus); ATCC VR-119 Norovirus (feline calicivirus tested surrogate) Rhinovirus (Type 37, Strain 151-1); ATCC VR-1147 Rotavirus (Strain WA) <i>Candida albicans</i> ; ATCC 10231
--	---	---

*OxyCide Daily Disinfectant Cleaner* has demonstrated effectiveness against influenza A virus and is expected to inactivate all influenza A viruses including Pandemic 2009 H1N1 influenza A virus.

**5 MINUTE CONTACT TIME:**  
Hepatitis B virus  
Hepatitis C virus

**DIRECTIONS FOR USE**  
It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.  
**COMBINATION DISINFECTANT AND CLEANING:**  
*OxyCide Daily Disinfectant Cleaner* is effective against the labeled organisms at 3.0 oz. per 1 gallon of water in hard water (400 ppm as CaCO3) and 5% blood serum on hard nonporous surfaces. For heavily soiled areas a precleaning step is required. Apply solution with mop, cloth, sponge, brush, scrubber, or coarse spray device or by soaking so as to wet all surfaces thoroughly. Allow to remain wet for required contact of 3 or 5 minutes and then either allow to air dry or if desired remove solution and entrapped soil with a clean wet mop, cloth, or wet vacuum pickup. Prepare a fresh solution daily or when it becomes soiled or diluted.

**DIRECTIONS FOR USE (CONTINUED)**  
**To Deodorize:** Apply *OxyCide Daily Disinfectant Cleaner* use solution to completely wet all surfaces. Let stand for 3 minutes to kill odor causing bacteria then wipe or allow to air dry. For heavily soiled areas, a precleaning is required.  
**To Clean and Deodorize Toilets:** To clean and deodorize toilet bowl, squirt liberally (1/2 cup) on toilet sides and upper toilet bowl rim. Swab or brush all surfaces let stand for 3 minutes and flush.  
**Dilution:** Disinfection 3.0 oz. per 1 gallon of water.  
*OxyCide Daily Disinfectant Cleaner* is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream, or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**SPORICIDAL**  
***Clostridium difficile***  
*OxyCide Daily Disinfectant Cleaner* kills and/or inactivates spores of *Clostridium difficile* on hard, non-porous surfaces. This product is effective against *Clostridium difficile* endospores after a 3 minute exposure time.  
**SPECIAL INSTRUCTIONS FOR CLEANING PRIOR TO DISINFECTION AGAINST *Clostridium difficile* SPORES**  
**Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks or eye covering.  
**Contact time:** Leave surface wet for 3 minutes with 3oz. per gallon use solution.  
**Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with the disinfectant product. This cleaning may be accomplished with any cleaning solution, including this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.  
**Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

**VIROCIDAL**  
*OxyCide Daily Disinfectant Cleaner* kills HIV and HBV and HCV on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces / objects with blood or body fluids, and in which the surfaces / objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) or Hepatitis B Virus (HBV) or Hepatitis C Virus (HCV).  
**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 OR HBV OR HCV ON SURFACES / OBJECTS SOILED WITH BLOOD / BODY FLUIDS.**  
**Personal protection:** Clean-up must always be done wearing protective gloves, gowns, masks and eye protection.  
**Cleaning procedure:** Blood and other body fluids containing HIV or HBV or HCV must be thoroughly cleaned from surfaces and objects before application of this product. This cleaning may be accomplished with any cleaning solution, including this product.  
**Contact time:** Leave surface wet for 3 minutes for HIV-1 and 5 minutes for HBV and HCV with 3oz. per gallon use-solution.  
**Disposal of infectious material:** Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.



**OxyCide™**  
**Daily Disinfectant Cleaner**

**LICENSED**  
PERIOD 2015-2017 LIC. NO. **9203.464**

**ONE-STEP DISINFECTANT / CLEANER / VIRUCIDE / SPORICIDE / DEODORIZER / KILLS GERMS / KILLS CLOSTRIDIUM DIFFICILE ENDOSPORES**

**LICENSED**  
PERIOD 2012-2014 LIC. NO.

- Kills *Clostridium difficile* spores in 3 minutes.
- Tough on C. diff spores but easy on surfaces and designed for daily use.
- Will not stain, discolor or bleach uniforms or fabrics.

**ACTIVE INGREDIENTS:**  
Hydrogen Peroxide . . . . . 27.5%  
Peroxyacetic Acid . . . . . 5.8%  
**INERT INGREDIENTS:** . . . . . 66.7%  
**TOTAL** . . . . . 100.00%

Contains no alkyl phenol ethoxylate detergents

**KEEP OUT OF REACH OF CHILDREN**  
**DANGER PELIGRO**  
See side panel for additional Precautionary and First Aid statements.

**CHLORINE BLEACH FREE**  
**FOR INSTUTIONAL USE**  
**STRONG OXIDIZING AGENT**

**XXXXXXXXX**  
Net Content:  
2.84 L (96 US fl oz)

**PRECAUTIONARY STATEMENTS**  
**HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

**DANGER: CORROSIVE:** Corrosive. Causes skin burns. Causes irreversible eye damage. May be fatal if inhaled. May be fatal if swallowed. Harmful if absorbed through the skin. Do not get in eyes, on skin, or on clothing. Wear coveralls worn over long-sleeved shirt and long pants, socks, chemical-resistant footwear, and gloves. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove and wash contaminated clothing before reuse. The following Personal Protective Equipment (PPE) should be used when handling the concentrated product: coveralls over long-sleeved shirt and long pants, socks and chemical-resistant footwear, goggles or face shield, chemical-resistant gloves (such as rubber or made out of any waterproof material), chemical-resistant apron. Wear a mask or pesticide respirator jointly approved by Mine Safety and Health Administration and the National Institute for Occupational Safety and Health. After product has been diluted according to label directions PPE is not required.  
**PHYSICAL AND CHEMICAL HAZARDS:** Strong oxidizing agent. Corrosive. Do not use in concentrated form. Mix only with water according to label instructions. Never bring concentrate in contact with other sanitizers, cleaners or organic substances.

**ENVIRONMENTAL HAZARDS:** This pesticide is toxic to birds, fish, and aquatic invertebrates. Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA.

**FIRST AID**  
**IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.  
**IF ON SKIN OR CLOTHING:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.  
**IF SWALLOWED:** Call a poison control center or doctor for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by a poison control center or doctor. Do not give anything by mouth to an unconscious person. Have the product container or label with you when calling a poison control center or going for treatment.  
**IF INHALED:** Move person to fresh air. If the person is not breathing, call 911 or ambulance, then give artificial respiration, preferably mouth-to-mouth, if possible. Call a Poison Control Center or doctor for further treatment advice.  
**FOR EMERGENCY MEDICAL INFORMATION, CALL TOLL-FREE 1.800.328.0026.**  
Have the product container or label with you when calling a poison control center or doctor, or going for treatment.  
**NOTE TO PHYSICIAN:** Probable mucosal damage may contraindicate the use of gastric lavage.

**STORAGE AND DISPOSAL**  
Do not contaminate water, food or feed by storage or disposal.  
**PESTICIDE STORAGE:** Product should be kept cool and in a vented container to avoid any explosion hazard.  
**PESTICIDE DISPOSAL:** Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance.  
**CONTAINER DISPOSAL:** Nonrefillable sealed container. Do not reuse or refill this container. Wrap empty container and put in trash.

For product information, call 1.866.781.8787.  
Ecolab Inc.  
370 Wabasha Street N • St Paul MN 55102-1390 USA  
© 2012 Ecolab USA Inc • All rights reserved  
Made in United States • OxyCide-B/EPA/0112

EPA Reg No. 1677-237  
EPA Est. 60156-IL-1(SI), 72806-OK-1 (AD), 1677-IL-2(J), 1677-TX-1(D), 1677-GA-1(M), 1677-CA-1(S), 1677-CA-2(R), 1677-MN-1(P), 1677-WV-1(V), 70271-CA-2(A), 303-IN-1(L).  
The superscript refers to first letter of date code.