



LICENSED

PERIOD 2014-2016 LIC. NO. 9308.9



LICENSED

PERIOD 2011-2013 LIC. NO.



LICENSED

PERIOD 2008-2010 LIC. NO.

BACTERICIDAL
VIRUCIDAL
TUBERCULOCIDAL
FUNGICIDAL
PSEUDOMONACIDAL

Manufactured by:
Medical Chemical Corporation
19430 Van Ness Avenue
Torrance, CA 90501
800-424-9394

U.S. Patents: 5,004,757;
5,219,890 and pending
EPA Reg. No 15136-9
EPA Est. No. 15136-CA-01

Wavicide-06 PlusTM
HARD SURFACE DISINFECTANT SPRAY
REVISED DEC. 9, 2008

ACTIVE INGREDIENTS:

Glutaraldehyde (1.5 Pentanedial).....0.25%
Ethanol 50.00%
Inert Ingredients: 49.75%
Contains 3 corrosion inhibitor

**-DANGER-
KEEP OUT OF REACH OF
CHILDREN**

Net Contents: 1 US gallon (3.79 liters)
FOR PROFESSIONAL AND INDUSTRIAL USE ONLY

DIRECTIONS FOR USE

It is a violation of federal law to use this product in a manner inconsistent with its labeling.

Recommended for the disinfection of hard non-porous surfaces such as countertops, tables, chairs, cages and other environmental surfaces in hospitals, laboratories and other institutions. Wavicide-06 Plus is compatible with rubber, latex, aluminum and stainless steel when used according to instructions.

Preclean all surfaces to be disinfected. This will insure maximum contact with the disinfection solution. DO NOT DILUTE. For effective disinfection, spray non-critical devices with full strength Wavicide-06 Plus at room temperature until saturated. Surfaces should be kept moist for 10 minutes. Dry surfaces with a clean cloth. This procedure will destroy vegetative pathogens including Staphylococcus aureus, Salmonella choleraesuis, Pseudomonas aeruginosa and pathogenic fungi (such as Trichophyton mentagrophytes) on hard, non-porous, inanimate surfaces. A ten minute contact time with this procedure will also destroy the following viruses on hard non-porous inanimate surfaces: Polio virus type 2, HIV-1, Adenovirus type 2 and Herpes simplex virus type 1.

Tuberculocidal at room temperature (22°C) and a 10 minute contact time (tested against M. bovis using the quantitative test method).

Special instructions for cleaning and decontaminating against HIV on surfaces/objects soiled with blood/body fluids. Observe universal precautions as defined by the CDC for surfaces soiled with blood/body fluids. Observe universal precautions as defined by the CDC for Medical Instruments or equipment soiled with blood/body fluids. Personal Protection: Disposable nitrile gloves, eye protection, and face masks should be worn when cleaning and disinfecting soiled equipment. Cleaning Procedure: Blood and other bodily fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant. Disinfection: Treat with Wavicide-06 Plus for ten minutes to destroy HIV-1 virus. Disposal of Infectious Material: Blood and other bodily fluids should be disposed of according to approved hospital procedures.

STORAGE AND DISPOSAL: Store at room temperature. **CONTAINER DISPOSAL:** Thoroughly rinse container and discard in trash. **PESTICIDE DISPOSAL:** Wastes resulting from the use of this product may be disposed of on site or at an approved waste disposal facility.

PRECAUTIONARY STATEMENTS: Hazardous to humans and domestic animals. This product is a sensitizer. **DANGER.** Corrosive. Causes irreversible eye damage. Do not get into eyes or on clothing. Wear goggles or face shield. Prolonged or repeated skin contact may cause allergic reaction in sensitive individuals. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash before reuse. **FLAMMABLE FLUID.** Contains alcohol. Keep away from heat, flame or other sources of ignition.

FIRST AID

IF SWALLOWED: Call poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not give anything by mouth to an unconscious person.

IF ON SKIN: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control for treatment advice.

IF INHALED: Move person to fresh air. If person is not breathing call 911 or an ambulance, then give artificial respiration, preferably mouth-to-mouth if possible. Call a poison control center or doctor for further treatment advice.

IF IN EYES: Hold eye open and rinse gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

This product is not to be used as a terminal sterilant/high level disinfectant on any surfaces or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.