

STERIS®**Coverage® Spray TB Plus**

Ready to Use

Order Number:
1629-B4Net Contents:
24 fl oz (710 mL)

Hospital / Medical / Surgical Hard and Soft Surface Sanitizer** & Disinfectant / Cleaner

Bactericidal • Fungicidal • Virucidal* • Tuberculocidal

EPA Reg. No. 70144-1-1043 EPA Est. 1130-IL-1

ACTIVE INGREDIENTS:

n-Alkyl (60% C ₁₄ , 30% C ₁₆ , 5% C ₁₂ , 5% C ₁₈) dimethyl benzyl ammonium chloride		0.154%
n-Alkyl (68% C ₁₂ , 32% C ₁₄) dimethyl ethylbenzyl ammonium chloride		0.154%
Isopropanol	Department of Agriculture STATE OF HAWAII	21.000%
OTHER INGREDIENTS		78.692%
TOTAL	LICENSED	100.000%

PERIOD **2013-2015** LIC. NO.**KEEP OUT OF REACH OF CHILDREN
CAUTION**

See back panel for First Aid and additional Precautionary Statements.

9445.118 Disc.

1629-61J(F)(0713)

62377

Made in the U.S.A.

NEST24S-82368B-0513

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

Apply this undiluted product by sponge, rag, cloth, mop, spray or immersion to wet all surfaces thoroughly.

When used as directed this product effectively kills the following microorganisms at room temperature with a: 10-second sanitization** contact time against *Enterobacter aerogenes*** , *Staphylococcus aureus*** , and a 2-minute disinfectant contact time against *Acinetobacter baumannii* , Avian Influenza A (H3N2) Reassortant Virus* , *Escherichia coli* (ESBL) (*E. coli*) , Hepatitis B Virus* (HBV) , Hepatitis C Virus* (HCV) , Herpes Simplex Virus* (HSV-2) , Human Immunodeficiency Virus* (HIV-1) (associated with AIDS) , Influenza A Virus* , *Klebsiella pneumoniae*-Carbapenem-resistant (KPC) , *Listeria monocytogenes* , Methicillin Resistant *Staphylococcus aureus* (MRSA) , *Mycobacterium bovis* BCG (tuberculosis) (TB) , *Pseudomonas aeruginosa* , Rotavirus* , *Salmonella enterica* , *Serratia marcescens* , *Staphylococcus aureus* , Swine Influenza A Virus* (H1N1) (formerly known as 'swine flu') , *Trichophyton mentagrophytes* (Ringworm) and Vancomycin-Resistant *Enterococcus faecalis* (VRE) . For disinfection of Avian Influenza A, use this product as a spray only for surfaces that are not conducive to treatment by immersion or excess liquid.

Use this product to clean, disinfect and sanitize** soft and hard non-food contact surfaces, food contact surfaces***, non-porous inanimate surfaces such as those made of acrylic, metals (aluminum, brass, stainless steel, chrome), glass, plastics (polycarbonate, polypropylene, polystyrene and polyvinylchloride), clear plastics, Plexiglas®, vinyl and leather surfaces as well as laminated and painted surfaces. (** Rinse with potable water before use.)

Hospital, Surgical & Home Healthcare Use

This product is effective in clinical areas, child care facilities, dialysis centers, drug preparation facilities***, emergency medical service (EMS) units, exam rooms, isolation rooms, infant care facilities, laboratories, neonatal intensive care units (NICUs), operating rooms, ophthalmology facilities, patient care rooms, pediatric units, personal care service facilities, pharmacies, physical therapy facilities, plasma pheresis facilities, primary care facilities, GI labs, sick rooms and in other similar establishments.

May be used on items such as anesthesia machines, monitors, beds, bedrails, bedside devices, biohazard team equipment, birthing room fixtures, blood bank equipment surfaces, blood pressure equipment, chest pumps, coated mattresses, computer keyboards and screens, colon equipment surfaces, drug preparation equipment, empty whirlpool tanks***, endoscopes, emergency room stretchers, exam room surfaces, floors, glucometers, gonial lenses, empty hydrotherapy units, infant care equipment, Isolettes®, infusion (IV) pump surfaces, lab benches, mats, mattress covers, operating room surfaces, oxygen hoods, ▶

Manufactured in the USA for: STERIS Corporation
7501 Page Avenue ■ St. Louis, MO 63133 ■ USA ■ 800-548-4873
1629-61JB(0713) 62378

NEST24S-82369B-0513



(01)00724995 083793

patient beds, patient curtains, physician office surfaces, respiratory therapy equipment, skin care equipment surfaces, telemetry equipment, tonometer equipment surfaces and other similar items.

SURFACES:

For surface cleaning: Prior to surface disinfection or sanitization**, pre-clean surfaces, making sure to remove all gross filth and heavy soil, by saturating surfaces with this product then immediately wipe surface clean using a fresh paper or cloth towel. Repeat, if necessary, until surfaces are visibly clean. Discard paper towels or laundry cloth towels before reusing.

For surface sanitization: Apply this product to clean surfaces by sponge, cloth, mop, spray or immersion to wet all surfaces thoroughly. Allow to remain visibly wet for 10 seconds at room temperature, then wipe surfaces using a clean mop, paper or cloth towel, or rinse surface using potable water or wipe surface dry or allow to air dry. Discard or launder spent materials before reuse.

For surface disinfection: Apply this product to pre-cleaned surfaces by sponge, cloth, mop, spray or immersion to wet all surfaces thoroughly. Allow to remain wet for 2 minutes, then wipe dry using a clean mop, paper or cloth towel, or rinse surface using potable water and wipe surface dry or allow to air dry. Discard or launder spent materials before reuse.

INSTRUMENTS / DEVICES:

For use as a pre-cleaning spray: Place item into a suitable container. Thoroughly spray item with this product until drenched, thoroughly rinse. Follow with an appropriate terminal cleaning and sterilization/disinfection process.

For use as a protective eye wear/loupes cleaner/disinfectant/sanitizer:** Thoroughly wet entire surface and wipe clean using a soft cloth. Reapply and allow to remain wet for 10 seconds for sanitization** or 2 minutes for disinfection and re-wipe using a fresh cloth. This product should not be used on contact lenses.

HBV, HCV and HIV-1 Precautions: This product kills HBV, HCV and HIV-1 on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects can be associated with the potential for transmission of HBV, HCV and/or the human immunodeficiency virus Type-1 (HIV-1) associated with AIDS. **Special instructions for cleaning and decontamination against HBV, HCV and HIV-1 of surfaces/objects soiled with blood/body fluids:** **Personal Protection:** Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks and eye coverings. **Cleaning Procedure:** Blood and body fluids must be thoroughly cleaned from surfaces/objects before disinfection. **Contact Time:** Allow the surface to remain wet for 2 minutes. **Infectious Materials Disposal:** Cleaning materials used that may contain blood/body fluids must be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

This product is not to be used as a terminal sterilant or high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter sterile areas of the body. This product may be used to pre-clean or decontaminate critical and semi-critical medical devices prior to sterilization or high-level disinfection.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

Pesticide Storage: Store in original container in a cool, well-ventilated area inaccessible to children and pets. **Pesticide Disposal:** Wastes resulting from the use of this product must be disposed of on site or at an approved waste disposal facility. **Container Disposal:** Refillable container. Refill this container with only the product originally supplied in this container. Do not reuse this container for any other purpose. When empty, fill the container about 10 percent full of water. Agitate vigorously and dispose of contents on site or at an approved waste disposal facility. Offer container for recycling, if available; otherwise dispose of container in trash.

FIRST AID

If Inhaled: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible. **If on Skin or Clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. **If in Eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. You may also contact 1-800-548-4873 for emergency medical treatment information.

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

CAUTION: May be harmful if inhaled or absorbed through skin. May cause moderate eye irritation. Avoid breathing vapors and contact with eyes, skin or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove contaminated clothing and wash clothing before reuse.

Physical or Chemical Hazards

Do not use or store near heat or open flame.