

Reorder Number: 1601-D2    Contents: 59.2 mL (2 fl. oz.)

**FOR SURGICAL INSTRUMENTS AND APPARATUS**

**Active Ingredient:**  
Bacterial Protease ..... 0.20%  
**Inert Ingredients** ..... 99.80%  
**Total** ..... 100.00%

CA Reg. No. 1043-50001

**SEE BACK PANEL FOR USE INSTRUCTIONS AND ADDITIONAL PRECAUTIONS**

**KEEP OUT OF REACH OF CHILDREN DANGER**



STATE OF HAWAII  
Department of Agriculture  
**ACCEPTED**  
LICENSE NO. **9445.80**

Reorder Number: 1601-D2    Contents: 59.2 mL (2 fl. oz.)

**EFFECTIVE IN THE PRESENCE OF SERUM AND HARD WATER**  
For Interim Instrument Decontamination prior to Terminal Cleaning and Sterilization of Surgical Instruments and Apparatus. Use only with **Vesta-Syde Enzyme Presoak**. Do not add other chemicals. Do not store instruments in the decontamination solution.

**FOR INSTITUTIONAL USE**

**Active Ingredients:**  
o-phenylphenol ..... 9.09%  
p-tertiary amylphenol ..... 7.66%  
**Inert Ingredients** ..... 83.25%  
**Total** ..... 100.00%

EPA Reg. No. 1043-114    EPA Est. No. 11563-MI-1

**KEEP OUT OF REACH OF CHILDREN DANGER**

FIRST AID	
If in eyes	<ul style="list-style-type: none"> <li>Hold eye open and rinse slowly and gently with water for 15-20 minutes.</li> <li>Remove contact lenses, if present, after the first 5 minutes, then continue rinsing the eye.</li> <li>Call a poison control center or doctor for treatment advice.</li> </ul>
If on skin or clothing	<ul style="list-style-type: none"> <li>Take off contaminated clothing.</li> <li>Rinse skin immediately with plenty of water for 15-20 minutes.</li> <li>Call a poison control center or doctor for treatment advice.</li> </ul>
If swallowed	<ul style="list-style-type: none"> <li>Call poison control center or doctor immediately for treatment advice.</li> <li>Have person sip a glass of water if able to swallow.</li> <li>Do not induce vomiting unless told to do so by the poison control center or doctor.</li> <li>Do not give anything by mouth to an unconscious person.</li> </ul>
If inhaled	<ul style="list-style-type: none"> <li>Move person to fresh air.</li> <li>If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible.</li> <li>Call a poison control center or doctor for further treatment advice.</li> </ul>
HOT LINE NUMBER	
Have the product container or label with you when calling a poison control center or doctor, or going for treatment.	
NOTE TO PHYSICIAN	
Probable mucosal damage may contraindicate the use of gastric lavage.	
1C19-67P(F)(313)	62214

FIRST AID	
If in eyes	<ul style="list-style-type: none"> <li>Hold eye open and rinse slowly and gently with water for 15-20 minutes.</li> <li>Remove contact lenses, if present, after the first 5 minutes, then continue rinsing the eye.</li> <li>Call a poison control center or doctor for treatment advice.</li> </ul>
If on skin or clothing	<ul style="list-style-type: none"> <li>Take off contaminated clothing.</li> <li>Rinse skin immediately with plenty of water for 15-20 minutes.</li> <li>Call a poison control center or doctor for treatment advice.</li> </ul>
If swallowed	<ul style="list-style-type: none"> <li>Call poison control center or doctor immediately for treatment advice.</li> <li>Have person sip a glass of water if able to swallow.</li> <li>Do not induce vomiting unless told to do so by the poison control center or doctor.</li> <li>Do not give anything by mouth to an unconscious person.</li> </ul>
If inhaled	<ul style="list-style-type: none"> <li>Move person to fresh air.</li> <li>If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible.</li> <li>Call a poison control center or doctor for further treatment advice.</li> </ul>
HOT LINE NUMBER	
Have the product container or label with you when calling a poison control center or doctor, or going for treatment.	
NOTE TO PHYSICIAN	
Probable mucosal damage may contraindicate the use of gastric lavage.	
1C18-67P(F)(313)	

5-3/8"

4-3/4"

2-1/4"

2-1/4"

**Vesta-Syde® Enzyme Presoak**

**FOR USE WITH VESTA-SYDE INTERIM INSTRUMENT DECONTAMINATION SOLUTION**

**Vesta-Syde Enzyme Presoak** is formulated for use on all surgical instruments and apparatus where blood, protein, and mucous form hard to remove soils or where scrubbing is impractical. This material is formulated for use as the first of two steps, for interim instrument decontamination, once surgery is completed. **Enzyme Presoak** has no biocidal efficacy, and is formulated for use as a digesting agent in proteinaceous soils. Equipment should be sterilized or disinfected by appropriate methods in separate subsequent steps.

**DIRECTIONS:** Dilute 2 ounces of **Vesta-Syde Enzyme Presoak** per gallon of water in the dirty instrument catch basin or instrument tray. Allow a minimum of 20 minutes contact time as a presoak for optimum results. After soaking for 20 minutes, add 2 ounces of **Vesta-Syde Interim Decontamination Solution** and allow an additional 20 minutes for the decontamination step to take place. If ultrasonic cleaning is desired, **Vesta-Syde Enzyme Presoak** may be used in an ultrasonic cleaning system. As with all ultrasonic cleansing procedures, the ultrasonic system should be covered. Allow a minimum of 20 minutes agitation time when using. Use only in conjunction with **Vesta-Syde Interim Instrument Decontamination Solution**. Do not add other chemicals. Added chemicals may destroy the effectiveness of this product.

**PRECAUTIONARY STATEMENT**

**DANGER:** Corrosive. Causes eye and skin damage. Wear goggles or face shield and rubber gloves when handling. Do not get in eyes, on skin, or clothing. Remove contaminated clothing; rinse and launder immediately. Clean up spills immediately. Wash hands thoroughly after use. Keep containers closed when not in use. Keep out of reach of children.

**STORAGE AND DISPOSAL**

Store in a cool dry place. Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available. Best if used before date on the pouch.

**FOR INSTITUTIONAL USE ONLY. NOT FOR HOME USE.**



(01)00724995 013462

Manufactured For  
STERIS Corporation

7501 Page Avenue ■ St. Louis, MO 63133 ■ USA  
800-548-4873 ■ www.steris.com

1C19-67P(B)(313)

Product Made in U.S.A.

**PRECAUTIONARY STATEMENTS HAZARD TO HUMANS AND DOMESTIC ANIMALS: DANGER**

Corrosive. Causes irreversible eye damage and skin burns. Harmful if swallowed, inhaled or absorbed through the skin. Do not get in eyes, on skin, or on clothing. Wear protective eyewear (goggles or face shield), protective clothing and rubber gloves. Wash thoroughly with soap and water after handling and before eating, drinking, using tobacco or using the toilet.

**This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.**

**DIRECTION FOR USE:** It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

**DECONTAMINATING SOILED SURGICAL INSTRUMENTS AND APPARATUS:** This product is formulated for use in conjunction with **Vesta-Syde Enzyme Presoak**. **Vesta-Syde Enzyme Presoak** is to be diluted 2 ounces per gallon of water (1:64) and a minimum 20 minute exposure time is to be allowed for **Vesta-Syde Enzyme Presoak** to work sufficiently to loosen proteinaceous soils. After a minimum 20 minute exposure to **Vesta-Syde Enzyme Presoak**, 2 ounces per gallon (1:64) of **Vesta-Syde Instrument Decontamination Solution** is to be deposited in the same vessel holding the soiled instruments and enzyme presoak solution. Gently stir the solution to provide a uniform mixture and thorough contact with the treated surfaces. Allow a minimum 20 minutes contact time to achieve interim instrument decontamination on all soiled instruments and apparatus deposited in the vessel. Avoid splashing and cover, where possible, when transporting soiled instruments.

**STORAGE AND DISPOSAL**

Do not contaminate water, food or feed by storage or disposal.

**PESTICIDE STORAGE:** Store in original container in areas inaccessible to persons unfamiliar with its use. **PESTICIDE DISPOSAL:** Pesticide waste may be acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal Law. If these waste cannot be disposed of by use according to label instructions, contact your State Pesticide or ENVIRONMENTAL Control Agency, or the Hazardous Waste Representative at the nearest EPA Regional Office for guidance. **CONTAINER DISPOSAL:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available. If not, puncture and dispose in sanitary landfill or by other procedures approved by state and local authorities.

Manufactured For  
STERIS Corporation

7501 Page Avenue ■ St. Louis, MO 63133 ■ USA  
800-548-4873 ■ www.steris.com

1C18-67P(B)(313)

62215

Product Made in U.S.A.

6-1/2"

5-3/16"



**Vesta-Syde®  
Enzyme Presoak**

For Use With VESTA-SYDE INTERIM INSTRUMENT  
DECONTAMINATION SOLUTION

Vesta-Syde Enzyme Presoak is formulated for use on all surgical instruments and apparatus where blood, protein, and mucous form hard to remove soils or where scrubbing is impractical. This material is formulated for use as the first of two steps, for interim instrument decontamination, once surgery is completed. Enzyme Presoak has no biocidal efficacy, and is formulated for use as a digesting agent in proteinaceous soils. Equipment should be sterilized or disinfected by appropriate methods in separate subsequent steps.

**DIRECTIONS:** Dilute 2 ounces of Vesta-Syde Enzyme Presoak per gallon of water in the dirty instrument catch basin or instrument tray. Allow a minimum of 20 minutes contact time as a presoak for optimum results. After soaking for 20 minutes, add 2 ounces of Vesta-Syde Interim Decontamination Solution and allow an additional 20 minutes for the decontamination step to take place. If ultrasonic cleaning is desired, Vesta-Syde Enzyme Presoak may be used in an ultrasonic cleaning system. As with all ultrasonic cleaning procedures, the ultrasonic system should be covered. Allow a minimum of 20 minutes agitation time when using. Use only in conjunction with Vesta-Syde Interim Instrument Decontamination Solution. Do not add other chemicals. Added chemicals may destroy the effectiveness of this product.

**PRECAUTIONARY STATEMENT**

**DANGER:** Corrosive. Causes eye and skin damage. Wear goggles or face shield and rubber gloves when handling. Do not get in eyes, on skin, or clothing. Remove contaminated clothing; rinse and launder immediately. Clean up spills immediately. Wash hands thoroughly after use. Keep containers closed when not in use. Keep out of reach of children.

**STORAGE AND DISPOSAL:** Store in a cool dry place. Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available.

**FOR INSTITUTIONAL USE.**

STERIS Corporation  
7501 Page Avenue  
St. Louis, MO 63133 ■ USA  
800-JIT-4-USE (800-548-4873)  
www.steris.com

160144Q(1113)

Product Made in U.S.A.

(Continued from Front Panel)

VESTA-SYDE INTERIM INSTRUMENT DECONTAMINATION SOLUTION represents a research-developed, state of the art compound for use with Vesta-Syde Enzyme Presoak as part of an interim instrument decontaminating step for soiled surgical instruments and apparatus. It substantially reduces the risk of exposure once surgery is completed. This product is effective in both hard water and in the presence of 5% blood serum. It is formulated as a concentrate for use in conjunction with the Enzyme Presoak. Its unique formulation allows for the addition of this product into the soiled instrument container in surgical suites or other procedure areas.

Although efficacy was shown at various dilutions and contact times, use this product at a 1:64 dilution for a 20 minute contact time for use against all of the organisms claimed.

**GERMICIDAL:** Passes A.O.A.C. Germicidal Use-Dilution Method (S. aureus, S. enterica, Ps. aeruginosa) when diluted with 400 ppm A.O.A.C. hard water to make a 1:128 (1 ounce per 1 gallon) solution, in the presence of 5% organic soil (serum), 10 minutes at 20°C.

**INTERIM INSTRUMENT DECONTAMINATION/DISINFECTION PRIOR TO TERMINAL CLEANING AND STERILIZATION TO REDUCE THE RISK OF CROSS CONTAMINATION:** When tested according to a protocol which comprised simulated contamination of instruments with a slurry containing proteinaceous matter derived from beef tissue and sheep's blood (50% by volume) containing A.O.A.C. specified concentrations of S. aureus, S. choleraesuis, and Ps. aeruginosa dried then exposed to Vesta-Syde Enzyme Presoak diluted 1:64 dilution in 400 ppm A.O.A.C. hard water for 20 minutes followed by the addition of Vesta-Syde Interim Instrument Decontamination Solution at a 1:64 dilution in the Enzyme Presoak solution completely inactivated S. aureus, S. choleraesuis, and Ps. aeruginosa in replicate tests in 20 minutes at 25°C (room temperature).

**BROAD SPECTRUM DATA:** In addition, the following organisms pass the A.O.A.C. Use-Dilution Test in 400 ppm A.O.A.C. hard water at a dilution of 1:128 in the presence of 5% organic soil (serum), 10 minutes at 20°C.

Acinetobacter calcoaceticus, ATCC 19606	Proteus vulgaris, ATCC 13315
Candida albicans, Clinical Isolate	Pseudomonas aeruginosa, ATCC 27853
Candida parapsilosis, Clinical Isolate	Pseudomonas cepacia, ATCC 25416
Citrobacter freundii, ATCC 8090	Salmonella typhi, ATCC 6539
Enterobacter aerogenes, ATCC 13048	Salmonella typhimurium, ATCC 14028
Enterobacter cloacae, ATCC 23355	Serratia marcescens, ATCC 8100
Escherichia coli, ATCC 25922	Shigella flexneri, ATCC 12022
Klebsiella pneumoniae, ATCC 13883	Shigella sonnei, ATCC 25931
Staphylococcus aureus, (MRSA),	Staphylococcus aureus, ATCC 25923
Multiply (Methicillin)-Resistant Clinical Isolate	Staphylococcus epidermidis, ATCC 12228
Proteus mirabilis, Clinical Isolate	Streptococcus faecalis, ATCC 19433
	Streptococcus pyogenes, ATCC 19615

Vesta-Syde Interim Instrument Decontamination Solution is effective in 3 minutes at 20°C against Pseudomonas aeruginosa ATCC 13388, according to the A.O.A.C. Use Dilution Test when diluted with 400 ppm hard water to make a 1:64 (2 ounces per gallon) solution, in the presence of 5% added organic soil (serum).

**FUNGICIDAL:** Passes A.O.A.C. Fungicidal Test (T. Mentagrophytes) when diluted with 400 ppm A.O.A.C. hard water to make a 1:128 (1 ounce per 1 gallon) solution, in the presence of 5% organic soil (serum), 10 minutes at 20°C.

**TUBERCULOCIDAL:** Passes A.O.A.C. Tuberculocidal Test (Mycobacterium tuberculosis, var. bovis) when diluted with 400 ppm A.O.A.C. hard water to make a 1:128 (1 ounce per 1 gallon) solution, in the presence of 5% organic soil (serum), 10 minutes at 20°C.

**VIROCIDAL:** Passes Virucidal Assay (AOAC Use Dilution Method) [Influenza A2 (Japan) virus, Herpes Simplex Type 2 virus, Vaccinia virus, Adenovirus Type 2 virus, Avian Infectious Bronchitis virus, Avian Laryngotracheitis virus, Avian Nucleocystic Disease virus, and Porcine Pseudorabies virus], when diluted with 400 ppm A.O.A.C. hard water to make a 1:128 (1 ounce per 1 gallon) solution in the presence of 5% organic soil (serum) in 10 minutes at 20°C.

When tested by an AOAC Use Dilution Method, the HIV-1(AIDS) virus, with added 10% organic soil (serum) was completely inactivated by a 1:128 (1 ounce per 1 gallon) solution in 400 ppm A.O.A.C. hard water in 60 seconds at 20-25°C. Although efficacy at 1 minute contact time has been shown to be adequate against HIV-1, this time would not be sufficient for other organisms. Use a 20-minute contact time for disinfection against all of the organisms claimed.

Kills Pandemic 2009 H1N1 influenza A virus formerly called swine flu.

**ODOR CONTROL:** It eliminates many odors by killing odor-causing bacteria while simultaneously chemically neutralizing their odors and leaving a pleasant fragrance. (Continued on Back Panel)



(01)10724995 013469

TEAR OFF TOP ON DOTTED LINE

52529



Reorder Number: 1601-D2 Net Contents: 24-2 fl. oz. (59.2 mL) packets EPA Est. No. 11563-MH-1 EPA Reg. No. 1043-114

**Interim Presoak**

**FOR SURGICAL INSTRUMENTS AND APPARATUS**

<b>Active Ingredient:</b>	
Bacterial Protease .....	0.20%
Inert Ingredients .....	99.80%
<b>Total .....</b>	<b>100.00%</b>

CA Reg. No. 1043-50001

**KEEP OUT OF REACH OF CHILDREN  
DANGER**

**SEE LEFT PANEL FOR USE INSTRUCTIONS AND ADDITIONAL PRECAUTIONS**

FIRST AID	
If in eyes	<ul style="list-style-type: none"> <li>Hold eye open and rinse slowly and gently with water for 15-20 minutes.</li> <li>Remove contact lenses, if present, after the first 5 minutes, then continue rinsing the eye.</li> <li>Call a poison control center or doctor for treatment advice.</li> </ul>
If on skin or clothing	<ul style="list-style-type: none"> <li>Take off contaminated clothing.</li> <li>Rinse skin immediately with plenty of water for 15-20 minutes.</li> <li>Call a poison control center or doctor for treatment advice.</li> </ul>
If swallowed	<ul style="list-style-type: none"> <li>Call poison control center or doctor immediately for treatment advice.</li> <li>Have person sip a glass of water if able to swallow.</li> <li>Do not induce vomiting unless told to do so by the poison control center or doctor.</li> <li>Do not give anything by mouth to an unconscious person.</li> </ul>
If inhaled	<ul style="list-style-type: none"> <li>Move person to fresh air.</li> <li>If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible.</li> <li>Call a poison control center or doctor for further treatment advice.</li> </ul>

**HOT LINE NUMBER**  
Have the product container or label with you when calling a poison control center or doctor, or going for treatment.  
**NOTE TO PHYSICIAN**  
Probable mucosal damage may contraindicate the use of gastric lavage.

**Interim Instrument Decontamination Solution**

**EFFECTIVE IN THE PRESENCE OF SERUM AND HARD WATER**  
For Interim Instrument Decontamination prior to Terminal Cleaning and Sterilization of Surgical Instruments and Apparatus. Use only with Vesta-Syde Enzyme Presoak. Do not add other chemicals. Do not store instruments in decontamination solution.

<b>Active Ingredients:</b>	
o-phenylphenol .....	9.09%
p-tertiary amylphenol .....	7.66%
Inert Ingredients .....	83.25%
<b>Total .....</b>	<b>100.00%</b>

**KEEP OUT OF REACH OF CHILDREN  
DANGER**

**SEE BACK PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS**

FIRST AID	
If in eyes	<ul style="list-style-type: none"> <li>Hold eye open and rinse slowly and gently with water for 15-20 minutes.</li> <li>Remove contact lenses, if present, after the first 5 minutes, then continue rinsing the eye.</li> <li>Call a poison control center or doctor for treatment advice.</li> </ul>
If on skin or clothing	<ul style="list-style-type: none"> <li>Take off contaminated clothing.</li> <li>Rinse skin immediately with plenty of water for 15-20 minutes.</li> <li>Call a poison control center or doctor for treatment advice.</li> </ul>
If swallowed	<ul style="list-style-type: none"> <li>Call poison control center or doctor immediately for treatment advice.</li> <li>Have person sip a glass of water if able to swallow.</li> <li>Do not induce vomiting unless told to do so by the poison control center or doctor.</li> <li>Do not give anything by mouth to an unconscious person.</li> </ul>
If inhaled	<ul style="list-style-type: none"> <li>Move person to fresh air.</li> <li>If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible.</li> <li>Call a poison control center or doctor for further treatment advice.</li> </ul>

**HOT LINE NUMBER**  
Have the product container or label with you when calling a poison control center or doctor, or going for treatment.  
**NOTE TO PHYSICIAN**  
Probable mucosal damage may contraindicate the use of gastric lavage.

(Continued from Front Panel)



**PRECAUTIONARY STATEMENTS  
HAZARD TO HUMANS AND DOMESTIC ANIMALS:  
DANGER**

Corrosive. Causes irreversible eye damage and skin burns. Harmful if swallowed, inhaled or absorbed through the skin. Do not get in eyes, on skin, or on clothing. Wear protective eyewear (goggles or face shield), protective clothing and rubber gloves. Wash thoroughly with soap and water after handling and before eating, drinking, using tobacco or using the toilet.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

**DIRECTIONS FOR USE:** It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

**DECONTAMINATING SOILED SURGICAL INSTRUMENTS AND APPARATUS:** This product is formulated for use in conjunction with Vesta-Syde Enzyme Presoak. The Enzyme Presoak is to be diluted 2 ounces per gallon of water (1:64) and a minimum 20 minute exposure time should be allowed for the Enzyme Presoak to work sufficiently to loosen proteinaceous soils. After a minimum 20 minute exposure to Enzyme Presoak, 2 ounces per gallon (1:64) of Vesta-Syde Interim Decontamination Solution is to be deposited in the same vessel holding the soiled instruments and enzyme presoak solution. Gently stir the solution to provide a uniform mixture and thorough contact with the treated surfaces. Allow a minimum of 20 minutes contact time to achieve interim instrument decontamination on all soiled instruments and apparatus deposited in the vessel. Avoid splashing and cover, where possible, when transporting soiled instruments.

**KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 (HUMAN IMMUNODEFICIENCY VIRUS OR AIDS VIRUS) OF SURFACES/OBJECTS SOILED WITH BLOOD/ BODY FLUIDS:**

**Personal Protection:** Wear appropriate barrier protection such as latex gloves, gowns, masks or eye coverings.

**Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of a 1:128 solution (1 fl. oz. per gallon). Prepare and apply solution as directed in paragraph above.

**Contact Time:** While the HIV-1 virus is inactivated in 1 minute, use a 20-minute contact time for disinfection of all organisms on this label. **Infectious Materials Disposal:** Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious disease.

**STORAGE AND DISPOSAL:** Do not contaminate water, food or feed by storage or disposal. **PESTICIDE STORAGE:** Store in original container in areas inaccessible to persons unfamiliar with its use. **PESTICIDE DISPOSAL:** Pesticide waste may acutely be hazardous. Improper disposal of excess pesticide, spray mixture or rinsate is a violation of Federal Law. If these waste cannot be disposed of by use according to label instructions, contact your State Pesticide or ENVIRONMENTAL Control Agency, or the Hazardous Waste Representative at the nearest EPA Regional Office for guidance. **CONTAINER DISPOSAL:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available. If not, puncture and dispose in sanitary landfill or by other procedures approved by state and local authorities.