

DARBY SURFACE DISINFECTANT CLEANER EFFECTIVELY KILLS HUMAN IMMUNODEFICIENCY VIRUS TYPE 1 (HIV-1)\* AND HEPATITIS B VIRUS (HBV)\* ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS IN HEALTHCARE SETTINGS OR OTHER SETTINGS IN WHICH THERE IS AN EXPECTED LIKELIHOOD OF SOILING OF INANIMATE SURFACES/OBJECTS WITH BLOOD/ BODY FLUIDS, AND IN WHICH THE SURFACES/ OBJECTS CAN BE ASSOCIATED WITH THE POTENTIAL FOR TRANSMISSION OF HIV-1\* AND HBV\*.

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1\* AND HBV\* AND OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUID:**

**Personal Protection:** Wear appropriate barrier protection such as protective gloves, gowns, masks and eye coverings.

**Cleaning Procedure:** Blood and other bodily fluids must be thoroughly cleaned from surfaces and objects before disinfection with Darby Surface Disinfectant Cleaner.

**Contact Time:** A contact time of 2 minutes is required for HIV-1\* and HBV\*. This contact time will not be sufficient for some of the organisms listed on the label. Refer to label for contact times against other organisms.

**Infectious Materials Disposal:** Cleaning materials used that may contain blood or other bodily fluids must be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

This product is not to be used as a terminal sterilant/ high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or

decontaminate critical or semi-critical medical devices prior to sterilization/high-level disinfection.

**FIRST AID**

**If in eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

**If on skin:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

For Chemical Emergency ONLY (spill, leak, fire, exposure or accident), call CHEMTREC at 1-800-424-9300. For ALL other inquiries about this product, call Customer Service at 1-800-645-2310, Monday through Friday, between 9 a.m. and 7 p.m. EST.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

**STORAGE AND DISPOSAL**

Do not contaminate water, food, or feed by storage and disposal.

**Pesticide Storage:** Store in a cool place.

**Pesticide Disposal:** Wastes resulting from the use of this product may be disposed of according to applicable Federal, State or local procedures.

**Container Disposal:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available.

To reorder: 800-645-2310 or darby.com

70-3510-1



9430504



DARBY DENTAL SUPPLY, LLC

# Darby Surface Disinfectant Cleaner

Surface Disinfectant/Decontaminant Cleaner

Bactericidal | Virucidal\* | Fungicidal | Tuberculocidal\*\*

**ACTIVE INGREDIENTS:**

Diisobutylphenoxethoxyethyl dimethyl benzyl ammonium chloride.....0.28%

Isopropanol.....17.20%

**INERT INGREDIENTS:**.....82.52%

TOTAL.....100.00%

KEEP OUT OF REACH OF CHILDREN

## CAUTION

See other panel for additional precautionary statements.

**Precautionary Statements**

**Hazards to Humans and Domestic Animals.**

**CAUTION**

Harmful if absorbed through skin. Causes moderate eye irritation. Avoid contact with eyes, skin or clothing. Wash hands before eating, drinking, chewing gum, using tobacco or using the toilet. Remove contaminated clothing and wash clothing before reuse.

Physical or Chemical Hazards: Combustible. Do not use or store near heat or open flame.

**Net Contents:**

24 fl oz/710 ml

(Net Weight 1.5 lbs./680 g.)

Distributed by:

Darby Dental Supply, LLC

300 Jericho Quadrangle

Jericho, NY 11753

Made in the USA

EPA Reg. No. 46781-6-90812

EPA Est. No. 46781-MI-1

Department of Agriculture  
STATE OF HAWAII

LICENSED

9430504

**EFFECTIVE AGAINST**

- *Mycobacterium tuberculosis* var: *bovis* (BCG)\*\* • *Staphylococcus aureus* • *Pseudomonas aeruginosa* • *Salmonella enterica* • *Trichophyton mentagrophytes* • Methicillin Resistant *Staphylococcus aureus* (MRSA) • Vancomycin Resistant *Enterococcus faecalis* (VRE) • Hepatitis B Virus (HBV)\* • Herpes Simplex Virus Types 1 and 2\* • Human Immunodeficiency Virus (HIV-1)\*

**DESCRIPTION**

Darby Surface Disinfectant Cleaner is a multi-purpose disinfectant/decontaminant cleaner for use on hard, non-porous inanimate surfaces.

It is useful in healthcare settings such as dental operatories, surgical suites and other critical care areas where control of cross contamination between treated surfaces is important.

Darby Surface Disinfectant Cleaner will effectively clean and disinfect, when used as directed, such items as dental unit instrument trays, operatory light switches, exterior surfaces of amalgamators and dental curing lights, exterior surfaces of anesthesia machines and respiratory therapy equipment, exterior surfaces of dental operatories, exterior surfaces of endodontic equipment, such as apex locators, pulp testers and motors, stools, floors, walls, handrails, door knobs, bathroom fixtures, exterior toilet surfaces, cabinets, keyboards and telephones that are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.

**DIRECTIONS FOR USE**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**FOR PRECLEANING NON-INSTRUMENT SURFACES PRIOR TO DISINFECTION:**

Remove gross filth and heavy soil loads. Spray Darby Surface Disinfectant Cleaner directly onto surface. Allow to remain visibly wet for 30 seconds. Wipe surface using a towel.

**FOR PRECLEANING INSTRUMENTS PRIOR TO DISINFECTION:**

Remove gross filth and heavy soil loads.

**For use as a precleaning decontaminant spray:** Place instruments in a container. Spray Darby Surface Disinfectant Cleaner so as to thoroughly drench all instrument surfaces and allow to remain visibly wet for 30 seconds. Remove and rinse instruments. Follow with

appropriate cleaning and disinfection process.

**FOR USE AS A VIRUCIDE\* ON NON-INSTRUMENT SURFACES, INCLUDING HIV-1\* AND HBV\* APPLICATIONS:**

Spray Darby Surface Disinfectant Cleaner directly onto precleaned surface, thoroughly wetting area to be treated. Allow surface to remain visibly wet for 2 minutes at room temperature (68°F/20°C).

**FOR USE AS A DISINFECTANT ON NON-INSTRUMENT SURFACES:**

Spray Darby Surface Disinfectant Cleaner directly onto precleaned surface, thoroughly wetting area to be disinfected. Allow surface to remain visibly wet for 3 minutes at room temperature (68°F/20°C) for Tuberculocidal\*\* activity and effectiveness against *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes* and *Salmonella enterica*.

Wipe surface using a towel or allow to air dry.

**FOR USE AS A VIRUCIDE\* ON PRECLEANED NON-CRITICAL MEDICAL DEVICES\*, INSTRUMENTS AND IMPLEMENTS INCLUDING HIV-1\* AND HBV\* APPLICATIONS:**

Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Using either a tray or an ultrasonic unit, immerse instruments in Darby Surface Disinfectant Cleaner for 2 minutes at room temperature (68°F/20°C).

**FOR USE AS A DISINFECTANT ON PRECLEANED NON-CRITICAL MEDICAL DEVICES\*, INSTRUMENTS AND IMPLEMENTS:**

Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Using either a tray or an ultrasonic unit, immerse instruments in Darby Surface Disinfectant Cleaner for 3 minutes at room temperature (68°F/20°C) for Tuberculocidal\*\* activity and effectiveness against *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes* and *Salmonella enterica*. Remove and rinse instruments. Wipe dry prior to use. Discard solution after each use.

\*Non-critical medical devices are items that come in contact only with intact skin.

Note: Critical and semi-critical devices must be followed by appropriate terminal sterilization/high-level disinfection process.



# Darby Surface Disinfectant Cleaner

Surface Disinfectant/Decontaminant Cleaner

**Bactericidal | Virucidal\* | Fungicidal | Tuberculocidal\*\***

## ACTIVE INGREDIENTS:

Diisobutylphenoxethoxyethyl dimethyl benzyl ammonium chloride .....	0.28%
Isopropanol.....	17.20%
<b>INERT INGREDIENTS</b> .....	<b>82.52%</b>
<b>TOTAL</b> .....	<b>100.00%</b>

KEEP OUT OF REACH OF CHILDREN

## CAUTION

### Precautionary Statements

**Hazards to Humans and Domestic Animals.**

### CAUTION

Harmful if absorbed through skin. Causes moderate eye irritation. Avoid contact with eyes, skin or clothing. Wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove contaminated clothing and wash clothing before reuse. See side panel for additional precautionary statements.

Physical or Chemical Hazards: Combustible. Do not use or store near heat or open flame.

EPA REG. NO.: 46781-6-90812

EPA EST. NO.: 46781-MI-1

**Reorder Number: 9430504**

**Net Contents: 12 x 24 oz Bottles (spray)**

To Reorder: **800-645-2310** or **Darby.com**

Distributed by:  
Darby Dental Supply, LLC  
300 Jericho Quadrangle  
Jericho, NY 11753

Made in the USA

## EFFECTIVE AGAINST

- *Mycobacterium tuberculosis* var. *bovis* (BCG)\*\* • *Staphylococcus aureus* • *Pseudomonas aeruginosa* • *Salmonella enterica* • *Trichophyton mentagrophytes* • Methicillin Resistant *Staphylococcus aureus* (MRSA) • Vancomycin Resistant *Enterococcus faecalis* (VRE) • Hepatitis B Virus (HBV)\* • Herpes Simplex Virus Types 1 and 2\* • Human Immunodeficiency Virus (HIV-1)\*

## DESCRIPTION

**Darby Surface Disinfectant Cleaner** is a multi-purpose disinfectant/decontaminant cleaner for use on hard, non-porous inanimate surfaces.

It is useful in healthcare settings such as dental operatories, surgical suites and other critical care areas where control of cross contamination between treated surfaces is important.

Darby Surface Disinfectant Cleaner will effectively clean and disinfect, when used as directed, such items as dental unit instrument trays, operatory light switches, exterior surfaces of amalgamators and dental curing lights, exterior surfaces of anesthesia machines and respiratory therapy equipment, exterior surfaces of dental operatories, exterior surfaces of endodontic equipment, such as apex locators, pulp testers and motors, stools, floors, walls, handrails, door knobs, bathroom fixtures, exterior toilet surfaces, cabinets, keyboards and telephones that are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.

## DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**FOR PRECLEANING NON-INSTRUMENT SURFACES PRIOR TO DISINFECTION:** Remove gross filth and heavy soil loads. Spray Darby Surface Disinfectant Cleaner directly onto surface. Allow to remain visibly wet for 30 seconds. Wipe surface using a towel.

**FOR PRECLEANING INSTRUMENTS PRIOR TO DISINFECTION:** Remove gross filth and heavy soil loads.

**For use as a precleaning decontaminant spray:** Place instruments in a container. Spray Darby Surface Disinfectant Cleaner so as to thoroughly drench all instrument surfaces and allow to remain visibly wet for 30 seconds. Remove and rinse instruments. Follow with appropriate cleaning and disinfection process.

**FOR USE AS A VIRUCIDE\* ON NON-INSTRUMENT SURFACES, INCLUDING HIV-1\* AND HBV\* APPLICATIONS:** Spray Darby Surface Disinfectant Cleaner directly onto precleaned surface, thoroughly wetting area to be treated. Allow surface to remain visibly wet for 2 minutes at room temperature (68°F/20°C).

**FOR USE AS A DISINFECTANT ON NON-INSTRUMENT SURFACES:** Spray Darby Surface Disinfectant Cleaner directly onto precleaned surface, thoroughly wetting area to be disinfected. Allow surface to remain visibly wet for 3 minutes at room temperature (68°F/20°C) for Tuberculocidal\*\* activity and effectiveness against *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes* and

*Salmonella enterica*.

Wipe surface using a towel or allow to air dry.

**FOR USE AS A VIRUCIDE\* ON PRECLEANED NON-CRITICAL MEDICAL DEVICES\*, INSTRUMENTS AND IMPLEMENTS INCLUDING HIV-1\* AND HBV\* APPLICATIONS:** Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Using either a tray or an ultrasonic unit, immerse instruments in Darby Surface Disinfectant Cleaner for 2 minutes at room temperature (68°F/20°C).

**FOR USE AS A DISINFECTANT ON PRECLEANED NON-CRITICAL MEDICAL DEVICES\*, INSTRUMENTS AND IMPLEMENTS:** Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Using either a tray or an ultrasonic unit, immerse instruments in Darby Surface Disinfectant Cleaner for 3 minutes at room temperature (68°F/20°C) for Tuberculocidal\*\* activity and effectiveness against *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes* and *Salmonella enterica*. Remove and rinse instruments. Wipe dry prior to use. Discard solution after each use.

\*Non-critical medical devices are items that come in contact only with intact skin.

Note: Critical and semi-critical devices must be followed by appropriate terminal sterilization/high-level disinfection process.

**DARBY SURFACE DISINFECTANT CLEANER EFFECTIVELY KILLS HUMAN IMMUNO DEFICIENCY VIRUS TYPE 1 (HIV-1)\* AND HEPATITIS B VIRUS (HBV)\* ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS IN HEALTHCARE SETTINGS OR OTHER SETTINGS IN WHICH THERE IS AN EXPECTED LIKELIHOOD OF SOILING OF INANIMATE SURFACES/OBJECTS WITH BLOOD/BODY FLUIDS, AND IN WHICH THE SURFACES/OBJECTS CAN BE ASSOCIATED WITH THE POTENTIAL FOR TRANSMISSION OF HIV-1\* AND HBV\*.**

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1\* AND HBV\* OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUID: Personal Protection:** Wear appropriate barrier protection such as protective gloves, gowns, masks and eye coverings.

**Cleaning Procedure:** Blood and other bodily fluids must be thoroughly cleaned from surfaces and objects before disinfection with Darby Surface Disinfectant Cleaner.

**Contact Time:** A contact time of 2 minutes is required for HIV-1\* and HBV\*. This contact time will not be sufficient for some of the organisms listed on the label. Refer to label for contact times against other organisms.

**Infectious Materials Disposal:** Cleaning materials used that may contain blood or other bodily fluids must be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal. This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body,

either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization/high-level disinfection.

## FIRST AID

**If in eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

**If on skin:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

For Chemical Emergency ONLY (spill, leak, fire, exposure or accident), call CHEMTREC at 1-800-424-9300. For ALL other inquiries about this product, call Customer Service at 1-800-645-2310, Monday through Friday, between 9 a.m. and 7 p.m. EST.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

## STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

**Pesticide Storage:** Store in a cool place.

**Pesticide Disposal:** Wastes resulting from the use of this product may be disposed of according to applicable Federal, State or local procedures.

**Container Disposal:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling, if available.



## Gallon Die Line

### 4.75" x 17.8125" with 0.125" radius corners

DARBY SURFACE DISINFECTANT CLEANER EFFECTIVELY KILLS HUMAN IMMUNODEFICIENCY VIRUS TYPE 1 (HIV-1)\* AND HEPATITIS B VIRUS (HBV)\* ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS IN HEALTHCARE SETTINGS OR OTHER SETTINGS IN WHICH THERE IS AN EXPECTED LIKELIHOOD OF SOILING OF INANIMATE SURFACES/OBJECTS WITH BLOOD/BODY FLUIDS, AND IN WHICH THE SURFACES/OBJECTS CAN BE ASSOCIATED WITH THE POTENTIAL FOR TRANSMISSION OF HIV-1\* AND HBV\*.

#### SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1\* AND HBV\* OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUID:

**Personal Protection:** Wear appropriate barrier protection such as protective gloves, gowns, masks and eye coverings.

**Cleaning Procedure:** Blood and other bodily fluids must be thoroughly cleaned from surfaces and objects before disinfection with Darby Surface Disinfectant Cleaner.

**Contact Time:** A contact time of 2 minutes is required for HIV-1\* and HBV\*. This contact time will not be sufficient for some of the organisms listed on the label. Refer to label for contact times against other organisms.

**Infectious Materials Disposal:** Cleaning materials used that may contain blood or other bodily fluids should be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

This product is not to be used as a terminal sterilant/ high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization/high-level disinfection.

#### FIRST AID

**If in eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

**If on skin:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

For Chemical Emergency ONLY (spill, leak, fire, exposure or accident), call CHEMTREC at 1-800-424-9300. For ALL other inquiries about this product, call Customer Service at 1-800-645-2310, Monday through Friday, between 9 a.m. and 7 p.m. EST.

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

#### STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

**Pesticide Storage:** Store in a cool place.

**Pesticide Disposal:** Wastes resulting from the use of this product may be disposed of according to applicable federal, state, or local procedures.

**Container Disposal:** Triple rinse. Then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration or, if allowed by state and local authorities, by burning. If burned, keep out of smoke. Nonrefillable container. Do not reuse or refill this container.

#### EFFECTIVE AGAINST

• *Mycobacterium tuberculosis* var. *bovis* (BCG)\*\* • *Staphylococcus aureus* • *Pseudomonas aeruginosa* • *Salmonella enterica* • *Trichophyton mentagrophytes* • Methicillin Resistant *Staphylococcus aureus* (MRSA) • Vancomycin Resistant *Enterococcus faecalis* (VRE) • Hepatitis B Virus (HBV)\* • Herpes Simplex Virus Types 1 and 2\* • Human Immunodeficiency Virus (HIV-1)\*

#### DESCRIPTION

**Darby Surface Disinfectant Cleaner** is a multi-purpose disinfectant/decontaminant cleaner for use on hard, non-porous inanimate surfaces.

It is useful in healthcare settings such as dental operatories, surgical suites and other critical care areas where control of cross contamination between treated surfaces is important.

**Darby Surface Disinfectant Cleaner** will effectively clean and disinfect, when used as directed, such items as dental unit instrument trays, operatory light switches, exterior surfaces of amalgamators and dental curing lights, exterior surfaces of anesthesia machines and respiratory therapy equipment, exterior surfaces of dental operatories, exterior surfaces of endodontic equipment, such as apex locators, pulp testers and motors, stools, floors, walls, handrails, door knobs, bathroom fixtures, exterior toilet surfaces, cabinets, keyboards and telephones which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.



9430506

703512-1



DARBY DENTAL SUPPLY, LLC

# Darby Surface Disinfectant Cleaner

## Surface Disinfectant/Decontaminant Cleaner

**Bactericidal | Virucidal\* | Fungicidal | Tuberculocidal\*\***

#### ACTIVE INGREDIENTS:

Disobutylphenoxethoxyethyl dimethyl benzyl ammonium chloride ..... 0.28%

Isopropanol ..... 17.20%

INERT INGREDIENTS: ..... 82.52%

TOTAL ..... 100.00%

#### Precautionary Statements

**Hazards to Humans and Domestic Animals.**

#### CAUTION

Harmful if absorbed through skin. Causes moderate eye irritation. Avoid contact with eyes, skin or clothing. Wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove contaminated clothing and wash clothing before reuse.

Physical or Chemical Hazards: Combustible. Do not use or store near heat or open flame.

Net Contents: ONE GALLON  
(128 fluid oz.) (3.785 liters)



LICENSED

PERIOD 2015-2017 LIC. NO.

9591.74

9430506

KEEP OUT OF REACH  
OF CHILDREN  
**CAUTION**

See other panel for additional precautionary statements.

Distributed by:  
Darby Dental Supply, LLC  
300 Jericho Quadrangle  
Jericho, NY 11753

Made in the USA

#### DIRECTIONS FOR USE

If is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**FOR PRECLEANING NON-INSTRUMENT SURFACES PRIOR TO DISINFECTION:**

Remove gross filth and heavy soil loads. Spray Darby Surface Disinfectant Cleaner directly onto surface. Allow to remain visibly wet for 30 seconds. Wipe surface using a towel.

**FOR PRECLEANING INSTRUMENTS PRIOR TO DISINFECTION:**

Remove gross filth and heavy soil loads.

**For use as a precleaning decontaminant immersion solution:** Place

instruments in a container with a sufficient amount of Darby Surface Disinfectant Cleaner so as to allow for complete submersion, cover and immerse for 10 minutes at room temperature (68°F/20°C). Remove and rinse instruments. Follow with appropriate cleaning and disinfection process. Discard solution after each use.

**For use as a precleaning decontaminant spray:** Place instruments in a container.

Spray Darby Surface Disinfectant Cleaner so as to thoroughly drench all instrument surfaces and allow to remain visibly wet for 30 seconds. Remove and rinse instruments. Follow with appropriate cleaning and disinfection process.

**For use as an ultrasonic precleaning solution:** Thoroughly rinse instruments under

running water to remove visible gross debris. Add 1 ounce Darby Surface Disinfectant Cleaner per liter of water to ultrasonic unit. Immerse instruments in solution and activate ultrasonic unit for 5 minutes or longer, if necessary. Remove and rinse instruments. Follow with appropriate cleaning and disinfection process. Discard solution after each use.

**For use as a manual precleaner:** Thoroughly rinse instruments under running water

to remove visible gross debris. Place instruments in a container with 1 ounce Darby Surface Disinfectant Cleaner per liter of water. Scrub instruments while submerged for 30 seconds using a stiff bristle brush until visibly clean. Remove and rinse instruments. Follow with appropriate cleaning and disinfection process. Discard solution after each use.

**FOR USE AS A VIRUCIDE\* ON NON-INSTRUMENT SURFACES,**

**INCLUDING HIV-1\* AND HBV\* APPLICATIONS:**

Spray Darby Surface Disinfectant Cleaner directly onto precleaned surface, thoroughly wetting area to be treated. Allow surface to remain visibly wet for 2 minutes at room temperature (68°F/20°C).

**FOR USE AS A DISINFECTANT ON NON-INSTRUMENT SURFACES:**

Spray Darby Surface Disinfectant Cleaner directly onto precleaned surface, thoroughly wetting area to be disinfected. Allow surface to remain visibly wet for 3 minutes at room temperature (68°F/20°C) for Tuberculocidal\*\* activity and effectiveness against

*Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes* and *Salmonella enterica*.

Wipe surface using a towel or allow to air dry.

**FOR USE AS A VIRUCIDE\* ON PRECLEANED NON-CRITICAL MEDICAL DEVICES', INSTRUMENTS AND IMPLEMENTS INCLUDING HIV-1\* AND HBV\* APPLICATIONS:**

Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Using either a tray or an ultrasonic unit, immerse instruments in Darby Surface Disinfectant Cleaner for 2 minutes at room temperature (68°F/20°C).

**FOR USE AS A DISINFECTANT ON PRECLEANED NON-CRITICAL MEDICAL DEVICES', INSTRUMENTS AND IMPLEMENTS:**

Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Using either a tray or an ultrasonic unit, immerse instruments in Darby Surface Disinfectant Cleaner for 3 minutes at room temperature (68°F/20°C) for Tuberculocidal\*\* activity and effectiveness against *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes* and *Salmonella enterica*.

Remove and rinse instruments. Wipe dry prior to use. Discard solution after each use.

\*Non-critical medical devices are items that come in contact only with intact skin.

Note: Critical and semi-critical devices must be followed by appropriate terminal sterilization/high-level disinfection process.

To reorder: **800-645-2310** or **darby.com**

EPA Reg. No. 46781-6-90812

EPA Est. No. 46781-WF-1

# Die #: 15761



## Darby Surface Disinfectant Cleaner

Surface Disinfectant/Decontaminant Cleaner

**Bactericidal | Virucidal\* | Fungicidal | Tuberculocidal\*\***  
**ACTIVE INGREDIENTS:**  
Diisobutylphenoxyethoxyethyl dimethyl benzyl ammonium chloride.....0.28%  
Isopropanol.....17.20%  
**INERT INGREDIENTS.....82.52%**  
**TOTAL.....100.00%**

**KEEP OUT OF REACH OF CHILDREN**

### CAUTION

#### Precautionary Statements

**Hazards to Humans and Domestic Animals.**

#### CAUTION

Harmful if absorbed through skin. Causes moderate eye irritation. Avoid contact with eyes, skin or clothing. Wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove contaminated clothing and wash clothing before reuse. See side panel for additional precautionary statements. Physical or Chemical Hazards: Combustible. Do not use or store near heat or open flame.

EPA REG. NO.: 46781-6-90812  
EPA EST. NO.: 46781-MI-1

**Reorder Number: 9430506**  
**Net Contents: 4 x 1 Gallon**

To reorder: **800-645-2310** or **darby.com**

Distributed by:  
Darby Dental Supply, LLC  
300 Jericho Quadrangle  
Jericho, NY 11753

Made in the USA

#### EFFECTIVE AGAINST

- *Mycobacterium tuberculosis* var: bovis (BCG)\*\* • *Staphylococcus aureus*
- *Pseudomonas aeruginosa* • *Salmonella enterica* • *Trichophyton mentagrophytes* • Methicillin Resistant *Staphylococcus aureus* (MRSA)
- Vancomycin Resistant *Enterococcus faecalis* (VRE) • Hepatitis B Virus (HBV)\*
- Herpes Simplex Virus Types 1 and 2\*\* • Human Immunodeficiency Virus (HIV-1)\*

#### DESCRIPTION

**Darby Surface Disinfectant Cleaner** is a multi-purpose disinfectant/decontaminant cleaner for use on hard, non-porous inanimate surfaces. It is useful in healthcare settings such as dental operatories, surgical suites and other critical care areas where control of cross contamination between treated surfaces is important.

**Darby Surface Disinfectant Cleaner** will effectively clean and disinfect, when used as directed, such items as dental unit instrument trays, operator light switches, exterior surfaces of amalgamators and dental curing lights, exterior surfaces of anesthesia machines and respiratory therapy equipment, exterior surfaces of dental operatories, exterior surfaces of endodontic equipment, such as apex locators, pulp testers and motors, stools, floors, walls, handrails, door knobs, bathroom fixtures, exterior toilet surfaces, cabinets, keyboards and telephones which are made of plastic (i.e., polycarbonate, polypropylene, polyvinylchloride, polystyrene), and use sites listed on this label consisting of non-porous vinyl, stainless steel, painted surfaces, Plexiglas®, Formica® and glass.

#### DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

#### FOR PRECLEANING NON-INSTRUMENT SURFACES PRIOR TO DISINFECTION:

Remove gross filth and heavy soil loads. Spray Darby Surface Disinfectant Cleaner directly onto surface. Allow to remain visibly wet for 30 seconds. Wipe surface using a towel.

#### FOR PRECLEANING INSTRUMENTS PRIOR TO DISINFECTION:

Remove gross filth and heavy soil loads.

**For use as a precleaning decontaminant immersion solution:** Place instruments in a container with a sufficient amount of Darby Surface Disinfectant Cleaner so as to allow for complete submersion, cover and immerse for 10 minutes at room temperature (68°F/20°C). Remove and rinse instruments. Follow with appropriate cleaning and disinfection process. Discard solution after each use.

**For use as a precleaning decontaminant spray:** Place instruments in a container. Spray Darby Surface Disinfectant Cleaner so as to thoroughly drench all instrument surfaces and allow to remain visibly wet for 30 seconds. Remove and rinse instruments. Follow with appropriate cleaning and disinfection process.

**For use as an ultrasonic precleaning solution:** Thoroughly rinse instruments under running water to remove visible gross debris. Add 1 ounce Darby Surface Disinfectant Cleaner per liter of water to ultrasonic unit. Immerse instruments in solution and activate ultrasonic unit for 5 minutes or longer, if necessary. Remove and rinse instruments. Follow with appropriate cleaning and disinfection process. Discard solution after each use.

**For use as a manual pre-cleaner:** Thoroughly rinse instruments under running water to remove visible gross debris. Place instruments in a container with 1 ounce Darby Surface Disinfectant Cleaner per liter of water. Scrub instruments while submerged for 30 seconds using a stiff bristle brush until

visibly clean. Remove and rinse instruments. Follow with appropriate cleaning and disinfection process. Discard solution after each use.

#### FOR USE AS A VIRUCIDE\* ON NON-INSTRUMENT SURFACES, INCLUDING HIV-1\* AND HBV\* APPLICATIONS:

Spray Darby Surface Disinfectant Cleaner directly onto precleaned surface, thoroughly wetting area to be treated. Allow surface to remain visibly wet for 2 minutes at room temperature (68°F/20°C).

#### FOR USE AS A DISINFECTANT ON NON-INSTRUMENT SURFACES:

Spray Darby Surface Disinfectant Cleaner directly onto precleaned surface, thoroughly wetting area to be disinfected. Allow surface to remain visibly wet for 3 minutes at room temperature (68°F/20°C) for Tuberculocidal\*\* activity and effectiveness against *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes* and *Salmonella enterica*. Wipe surface using a towel or allow to air dry.

#### FOR USE AS A VIRUCIDE\* ON PRECLEANED NON-CRITICAL MEDICAL DEVICES<sup>1</sup>, INSTRUMENTS AND IMPLEMENTS INCLUDING HIV-1\* AND HBV\* APPLICATIONS:

Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Using either a tray or an ultrasonic unit, immerse instruments in Darby Surface Disinfectant Cleaner for 2 minutes at room temperature (68°F/20°C).

#### FOR USE AS A DISINFECTANT ON PRECLEANED NON-CRITICAL MEDICAL DEVICES<sup>1</sup>, INSTRUMENTS AND IMPLEMENTS:

Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Thoroughly clean and rinse lumens of hollow instruments. Using either a tray or an ultrasonic unit, immerse instruments in Darby Surface Disinfectant Cleaner for 3 minutes at room temperature (68°F/20°C) for Tuberculocidal\*\* activity and effectiveness against *Staphylococcus aureus*, *Pseudomonas aeruginosa*, *Trichophyton mentagrophytes* and *Salmonella enterica*. Remove and rinse instruments. Wipe dry prior to use. Discard solution after each use.

<sup>1</sup>Non-critical medical devices are items that come in contact only with intact skin. Note: Critical and semi-critical devices must be followed by appropriate terminal sterilization/high-level disinfection process.

DARBY SURFACE DISINFECTANT CLEANER EFFECTIVELY KILLS **HUMAN IMMUNODEFICIENCY VIRUS TYPE 1 (HIV-1)\* AND HEPATITIS B VIRUS (HBV)\*** ON PRECLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS IN HEALTHCARE SETTINGS OR OTHER SETTINGS IN WHICH THERE IS AN EXPECTED LIKELIHOOD OF SOILING OF INANIMATE SURFACES/OBJECTS WITH BLOOD/BODY FLUIDS, AND IN WHICH THE SURFACES/OBJECTS CAN BE ASSOCIATED WITH THE POTENTIAL FOR TRANSMISSION OF HIV-1\* AND HBV\*.

#### SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION

#### TION AGAINST HIV-1\* AND HBV\* OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUID:

**Personal Protection:** Wear appropriate barrier protection such as protective gloves, gowns, masks and eye coverings.

**Cleaning Procedure:** Blood and other bodily fluids must be thoroughly cleaned from surfaces and objects before disinfection with Darby Surface Disinfectant Cleaner.

**Contact Time:** A contact time of 2 minutes is required for HIV-1\* and HBV\*. This contact time will not be sufficient for some of the organisms listed on the label. Refer to label for contact times against other organisms.

**Infectious Materials Disposal:** Cleaning materials used that may contain blood or other bodily fluids should be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

This product is not to be used as a terminal sterilant/ high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization/high-level disinfection.

#### FIRST AID

**If in eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice.

**If on skin:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

For Chemical Emergency ONLY (spill, leak, fire, exposure or accident), call CHEMTREC at 1-800-424-9300. For ALL other inquiries about this product, call Customer Service at 1-800-645-2310, Monday through Friday, between 9 a.m. and 7 p.m. EST. Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

#### STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

**Pesticide Storage:** Store in a cool place.

**Pesticide Disposal:** Wastes resulting from the use of this product may be disposed of according to applicable federal, state, or local procedures.

**Container Disposal:** Triple rinse. Then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, or incineration or, if allowed by state and local authorities, by burning. If burned, keep out of smoke. Nonrefillable container. Do not reuse or refill this container.

5 in

8 in

Press Web Direction →