

READ LABEL BEFORE USING
LOOSEN CAP TO FILL RESERVOIR - DO NOT SQUEEZE BOTTLE WHILE POURING



STATE OF HAWAII
Department of Agriculture

ACCEPTED
LICENSE NO. 9791.3

ISABEL CRISTINA

LET'S TOUCH®

HOSPITAL AND TUBERCULOCIDAL METAL SALON DISINFECTANT

DOES NOT RUST • DOES NOT DULL

**For Disinfection and Storage of All
Quality Metal Instruments Used In Salons**

**Bactericidal Fungicidal Tuberculocidal
Virucidal Staphylocidal Pseudomonacidal**

Active Ingredients:	
o-Benzyl-p-chlorophenol	5.25%
o-Phenylphenol	0.99%
Other Ingredients	93.76%
TOTAL	100.00%

KEEP OUT OF REACH OF CHILDREN

WARNING

First Aid

If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

If on skin or clothing: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

If swallowed: Call a poison control center immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by the poison control center or doctor. Do not give anything by mouth to an unconscious person.

HOT LINE NUMBER

Have the product container or label with you when calling a poison control center or doctor, or going for treatment. You may contact the Poison Control Center 1-800-222-1222 for emergency medical treatment information.

SEE BACK PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS

Net Contents: 32 fl. oz. (1 qt.) (946 ml.)

Manufactured by: RBR Productions Inc./Isabel Cristina
P.O. BOX 766, Hartwell, GA 30643, USA
Tel: (800) 247-4130 • 706-376-8189

EPA Reg. No. 62296-1

EPA Est. No. 3862-GA-1



LET'S TOUCH®

LET'S TOUCH® is a hospital and tuberculocidal disinfectant concentrate especially formulated for the beauty care profession. It provides broad spectrum protection: bactericidal, fungicidal, virucidal, staphylocidal, pseudomonacidal, tuberculocidal, and prevents cross contamination from all metal salon instruments when used as directed. Use it to disinfect all pre-cleaned metal salon instruments including manicurist nippers and cuticle pushers, shears and metal skin care instruments. Let's Touch® does not rust or dull quality metal instruments, is pH buffered, and is biodegradable.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Before immersion in Let's Touch®, all instruments and devices must be thoroughly cleaned to remove excess dirt, rinsed and rough dried. The removal of heavy soils and the cleaning of all surfaces prior to application of this product is required. When using Let's Touch® for the first time, thoroughly clean and rinse trays and containers in order to remove any residues. Mixing Let's Touch® in a container which had another disinfectant in it can cause a brown gel-like substance to form. If this happens, thoroughly clean and rinse all contaminated instruments, trays or containers and remix Let's Touch®.

Use one part Let's Touch® to 32 parts water for disinfection of pre-cleaned metal instruments and equipment. Let's Touch® is distributed in concentrated form in Twin-Neck containers and in 1 oz. bottles.

DIRECTIONS FOR MIXING WITH BOTTLES AND TWIN-NECK CONTAINER

Water	1 Pint (16 oz.)	1 Quart (32 oz.)	1/2 Gallon (64 oz.)	1 Gallon (128 oz.)
Let's Touch®	1/2 fl. oz.	1 fl. oz. or 1 bottle	2 fl. oz. or 2 bottles	4 fl. oz. or 4 bottles

Always add 1 bottle or 1 oz. of Let's Touch® to every 32 fl. oz. of water. When using the Twin-Neck container, squeeze the bottle until the easy-measure reservoir is filled to the desired amount.

DIRECTIONS FOR MIXING WITH THE MANICURIST AND HAIR STYLIST SYSTEMS

Manicurist System

1. Fill the tall plastic mixing jar to the 32 oz. Fill Line with water.
2. Add 1 oz. of Let's Touch® (1 Bottle). Cap & mix by turning upside down a few times.
3. Pour the mixed solution into the glass manicurist jar to the suggested use level line.
4. Store the remaining solution until needed.

Hairstylist System

1. Remove the plunger mechanism.
2. Fill the glass jar to the 32 oz. Fill Line with water.
3. Add 1 oz. of Let's Touch® (1 Bottle). Mix solution by replacing plunger and moving it up and down.

Always leave instruments in Let's Touch® solution at all times and overnight for storage. Replace solution daily, or earlier if cloudy or dirty.

FOR INSTRUMENT DISINFECTION: After cleaning instruments, place articles in Let's Touch® solution for ten minutes.

NOTE: Articles must be fully submerged in solution at all times.

FOR METAL INSTRUMENT STORAGE: A solution of one part Let's Touch® to 32 parts water (see above for mixing instructions) does not rust or dull, and does not stain or otherwise attack metal instruments such as quality salon shears, manicurist and skin care implements. The solution may be used for prolonged (overnight) storage of all but non metal instruments, plastic or rubber items.

PEEL OFF THIS PANEL FOR ADDITIONAL INFORMATION AND PRECAUTIONARY STATEMENTS

LET'S TOUCH® IS A COMPLETE PRODUCT. DO NOT MIX WITH OTHER CHEMICALS. Use only as directed. If frozen, thaw and remix before use.

Let's Touch® is recommended for use in salons, hospitals, dental offices or other facilities requiring disinfection of metal instruments.

FOR HOSPITAL AND DENTAL OFFICE USE: Let's Touch® is a concentrated disinfectant to be used at a 1:32 dilution with water for cold decontamination/disinfection of precleaned equipment. For disinfection immerse instruments for 10 minutes.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced with the bloodstream or normal sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

The 1:32 dilution kills pathogenic bacteria and fungi such as *S. Aureus*, *S. cholersuis*, *Ps. aeruginosa*, *M tuberculosis* (effective in 10 minutes at 20° C), and *C. albicans*. Let's Touch® is virucidal against Influenza A2 (Hong Kong), Herpes Simplex 1 and 2, Adenovirus Type 2 and Vaccina on hard inanimate surfaces in 10 minutes at 20° C. A 1:32 dilution is effective against *M tuberculosis* on hard, inanimate surfaces in 10 minutes at 20° C. Remove heavy soil or gross filth and thoroughly clean surfaces.

When tested by the EPA-approved Dilution Method the HIV-1 (AIDS) virus was completely inactivated by a 1:32 solution in 10 minutes at 25° C.

PRECAUTIONARY STATEMENTS

HAZARD TO HUMANS AND DOMESTIC ANIMALS

WARNING: Causes substantial but temporary eye injury. Do not get in eyes or on clothing. Wear protective eyewear (goggles, face shield or safety glasses). Harmful if swallowed. Avoid prolonged skin contact. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash clothing before reuse.

STORAGE AND DISPOSAL

Do not contaminate water food or feed by storage or disposal.

Pesticide Storage: Store only in original container away from heat or open flame. Keep container closed when not in use. Keep out of reach of children.

Pesticide Disposal: This germicide, its solutions or rinsing from empty containers should be disposed of in a toilet or service sink served by a sanitary sewer or in a landfill approved for pesticides.

Container Disposal: Triple rinse and deposit in a waste container for incineration or burial in a landfill approved for pesticide containers.

KILLS HIV (AIDS VIRUS) ON PRECLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV (AIDS VIRUS) OF SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS

Personal Protection: Disposable latex or vinyl gloves, gowns, masks and/or eye coverings, as appropriate must be worn during all cleaning and decontamination procedures of blood and other body fluids.

Cleaning Procedures: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying disinfectant.

Disinfectant Use and Contact Time: Effective against HIV-1 (AIDS virus) on hard non-porous surfaces/objects in the presence of organic soil (e.g 5% blood serum). Prepare disinfectant by mixing one part Let's Touch® to 32 parts water. Leave surfaces wet for 10 minutes.

Disposal of Infectious Material: Blood and other fluids should be autoclaved and disposed of according to Federal, State and local regulations for infectious waste disposal.