



PRECISE®

HOSPITAL FOAM CLEANER DISINFECTANT



- Cleans, Disinfects, Deodorizes
- For Hard, Non-Porous Surfaces
- Won't Scratch
- Effective against HIV-1 (AIDS Virus), Hepatitis B Virus (HBV), Hepatitis C Virus (HCV) in 1 Minute
- Effective against Multiple Microorganisms in 10 Minutes:

Mycobacterium bovis (TB); Staphylococcus aureus; Pseudomonas aeruginosa; Streptococcus pyogenes; Salmonella choleraesuis; Vaccinia Virus ATCC VR-117; Influenza A-2, strain Hong Kong (NCDC); Herpes Simplex Virus, Hominus, Type 2 (HSV-2); Poliovirus Type 1, Mahoney; Canine Parvovirus; Trichophyton mentagrophytes

Active Ingredients: o-Phenylphenol 0.37%
Other Ingredients..... 99.63%
TOTAL 100.00%



Contains no CFCs or other ozone depleting substances. Federal regulations prohibit CFC propellants in aerosols.

KEEP OUT OF REACH OF CHILDREN
WARNING

READ COMPLETE PRECAUTIONS AND DIRECTIONS ON REAR PANELS

NET WT. 6 oz.

40510-200-2

SHAKE OFTEN

FOR HOSPITAL HARD SURFACES: Quick spray and wipe for cleaning. Always properly mixed and ready to use. Put this product to work throughout the hospital. **USEFUL:** Patient room, emergency room, recovery room, operating room, laboratory, clinic, bathroom, physician office, physical therapy, nursing station, radiology, morgue, dental operator, veterinary hospital.

STORAGE AND DISPOSAL: Do not contaminate water, food, or feed by storage and disposal. **Storage:** Do not use or store near heat or open flame. Store out of reach of children. **Disposal:** This container may be recycled in the few but growing number of communities where (steel) aerosol can recycling is available. Before offering for recycling, empty the can by using the product according to the label. (DO NOT PUNCTURE!) If recycling is not available, wrap the container and discard in the trash.



PRECAUTIONARY STATEMENTS: Hazards to Humans and Domestic Animals. WARNING: Do Not Spray in Eyes. Causes mild eye irritation. Avoid contact with eyes or clothing. Wash hands before eating/drinking/using toilet. **FIRST AID: IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Have the product with you when calling a Poison Control Center or doctor or going for treatment. In case of emergency, call a doctor or National Poison Control Center at 800-222-1222. Do not puncture or incinerate container. Exposure to temperatures above 130°F may cause bursting.

©2007 Caltech Industries, Inc.
Midland, MI 48642
1-800-234-7700

®Trademark of Caltech Industries, Inc.

EPA Reg. No. 56392-1
EPA Est. Nos. 62801-MI-1
33915-MI-1

REORDER NO. 40510

DISC.

- 0 = NONE
1 = SLIGHT
2 = MODERATE
3 = HIGH
4 = EXTREME



Department of Agriculture
STATE OF HAWAII
LICENSED
2014-2016 LIC. NO. 8013.78

SHAKE OFTEN

CLEANS AND SHINES: All hospital washable surfaces - tile, tubs, bowls, basins, sinks, showers, urinals, empty bed pans, toilets, cabinets, doors and door handles, empty emesis basins, floors, walls, woodwork, workstation tables and countertops, housekeeping carts, vanity tops, vending machines, bed springs and frames, bedside furniture, dressing and surgical carts, wheelchairs, walkers, stainless steel, chrome, plastic, metal, vinyl, fiberglass, baked enamel, porcelain, examination and x-ray tables, film cassettes, chart holders, shower chairs, sitz baths, commodes and commode seats, instrument cabinets and trays, catheter trays, electrodes, bacteriology laboratory benches, centrifuges, whirlpool baths, operating tables, emergency vehicles, telephones, reception counters, dental chairs, medical equipment. **DISINFECTS:** Kills organisms commonly found on hard nonporous hospital surfaces. See directions for disinfecting. **DEODORIZES:** Kills odor causing organisms on hard surfaces. Leaves surfaces smelling clean and fresh. **LABOR SAVING:** Just spray on, wipe off. No rubbing, rinsing or scouring. **WON'T SCRATCH:** Contains no abrasive to scratch or leave gritty deposit. **REMOVES HOSPITAL SOILS:** Grime, greasy soil, soap scum, water marks, scuff marks, blood, drainage, stains, emesis, fingermarks.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. **DIRECTIONS:** Shake well. Hold can 6" from surface. **TO CLEAN:** Quick spray and wipe clean with damp sponge, cloth or paper towel. No rinsing needed. **TO CLEAN AND DISINFECT IN 1 STEP:** Spray surface thoroughly until completely covered with foam and allow to remain wet for the appropriate time, depending on the organism (See front panel). Wipe clean with damp sponge, cloth or paper towel. No rinsing needed. Kills Mycobacterium bovis (TB) in 10 minutes at room temperature (68-77°F, 20-25°C). Gross soil must be removed prior to disinfecting.

LIC. NO. 8013.78

Kills HIV-1 (AIDS Virus), HBV, and HCV on precleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV). **Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV, and HCV on Surfaces/Objects Soiled with Blood/Body Fluids: Personal Protection:** When handling items soiled with blood/body fluids, appropriate personal protective equipment such as disposable plastic gloves, gowns, masks, or eye coverings should be considered. **CLEANING PROCEDURE:** Blood/body fluids must be thoroughly cleaned from surfaces and objects prior to disinfecting. **CONTACT TIME:** This product kills HIV-1 (AIDS Virus), HBV, and HCV in 1 minute at room temperature. Contact time required to kill other organisms is 10 minutes. **DISPOSAL OF INFECTIOUS MATERIALS:** Blood/body fluids and cleaning materials should be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

THIS PRODUCT IS NOT TO BE USED AS A TERMINAL STERILANT/HIGH LEVEL DISINFECTANT ON ANY SURFACE OR INSTRUMENT THAT (1) IS INTRODUCED DIRECTLY INTO THE HUMAN BODY, EITHER INTO OR IN CONTACT WITH THE BLOODSTREAM OR NORMALLY STERILE AREAS OF THE BODY, OR (2) CONTACTS INTACT MUCOUS MEMBRANES BUT WHICH DOES NOT ORDINARILY PENETRATE THE BLOOD BARRIER OR OTHERWISE ENTER NORMALLY STERILE AREAS OF THE BODY. THIS PRODUCT MAY BE USED TO PRECLEAN OR DECONTAMINATE CRITICAL OR SEMI-CRITICAL MEDICAL DEVICES PRIOR TO STERILIZATION OR HIGH LEVEL DISINFECTION.