



Note to Customer: Weiman Products adds a revision code to all artwork items for tracking purposes. As the customer & controller of your own artwork, we ask that you either denote a specific area for this revision code and/or ask Weiman for the specific revision code that is to be added & place on artwork.



LICENSED

16.75 inches Width

PERIOD 2013-2015 LIC. NO. 8262.6

0.9593

4.25 inches Height

0.5023

When used as directed this product effectively kills the following microorganisms: *Acinetobacter baumannii*, Avian Influenza A Virus*, Hepatitis B Virus* (HBV), Hepatitis C Virus* (HCV), Herpes Simplex II Virus* (HSV 2) G Strain, Human Immunodeficiency Virus* (HIV-1) (associated with AIDS), Influenza Virus* Strain A2/Hong Kong, Methicillin-Resistant *Staphylococcus aureus* (MRSA), *Mycobacterium bovis*, BCG (tuberculosis) (TB), Pandemic 2009 H1N1 Influenza A Virus*, Polio I Virus*, *Pseudomonas aeruginosa*, Rhinovirus*, Rotavirus* Strain WA, *Salmonella enterica*, *Staphylococcus aureus*, *Trichophyton mentagrophytes* and Vancomycin-Resistant *Enterococcus faecium* (VRE).

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.
DO NOT USE THIS PRODUCT AS A BABY OR SKIN WIPE, OR ON CONTACT LENSES!
 This product is intended for use on hard, non-porous inanimate surfaces in areas where the disinfection of surfaces is required. When used as directed this product cleans and disinfects surfaces found in areas/facilities concerned with hazards of cross-contamination from infectious microorganisms: dental facilities, schools, laboratories, clinical areas, patient care areas, dental offices, hospitals, and areas/facilities concerned with the hazards of cross contamination from infectious microorganisms.
 This product may be used on dental equipment/surfaces/devices, counters, curing lights, finished surfaces, office fixtures, eye protection devices, computer keyboards, telephones, X-ray equipment/surfaces and devices/surfaces potentially contaminated with disease-causing microorganisms [and use sites made of glass, metal, Plexiglas®, plastics (i.e., polycarbonate, polypropylene, polyvinylchloride and polystyrene), baked-on painted surfaces, stainless steel and non-porous vinyl.]

DISPENSING DIRECTIONS

To start feed, remove cover and discard seal. From the center of the wipe roll, pull up a wipe corner, twist it into a point and thread it through the hole located on the container cover. Pull through about one inch. Replace cover. Pull out first wipe and tear off at an angle. Remaining wipes feed automatically, ready for the next use. When not in use, keep center cap closed to prevent evaporation.

CLEANING AND DISINFECTION INSTRUCTIONS

Use a fresh wipe to pre-clean surfaces of all gross filth and heavy soil. Repeat as necessary until surfaces are visibly clean. To effectively disinfect the pre-cleaned surfaces, use a new wipe to thoroughly wet the surfaces and allow surfaces to remain wet for 3 minutes at room temperature. Use enough wipes to ensure that the surface remains wet for the contact time. Allow surfaces to air dry before use, or if desired, rinse surface using potable water and wipe or air dry. For disinfection of Avian Influenza A, use this product only for surfaces that are not conducive to treatment by immersion or excess liquid (e.g., equipment tables, lab benches, faucet handles, door knobs and similar 'high touch' surfaces.)

Opti-Cide³

Virucidal*
 Bactericidal / Fungicidal
 Tuberculocidal
 Surface Wipes

Hospital & Dental Disinfectant, Cleaner & Decontaminant

KEEP OUT OF REACH OF CHILDREN
CAUTION

See side panels for First Aid and additional Precautionary Statements.

ACTIVE INGREDIENTS:

n-Alkyl (60% C₁₄, 30% C₁₆, 5% C₁₂, 5% C₁₈)
 dimethyl benzyl ammonium chloride 0.154%
 n-Alkyl (68% C₁₂, 32% C₁₄)
 dimethyl ethylbenzyl ammonium chloride .. 0.154%

OTHER INGREDIENTS 99.692%

TOTAL 100.000%

EPA REG. NO 70144-2-51003 EPA EST. 56952-WI-001



Net Contents: 100 Wipes
 Net Weight: 1 lb 7.8 oz. (676.6 g)
 Order No.: DOCW06-100

HBV, HCV & HIV-1 PRECAUTIONS

This product effectively inactivates HBV, HCV and HIV-1 on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids and in which the surfaces/objects can be associated with the potential for transmission of HBV, HCV and/or the human immunodeficiency virus Type 1 (HIV-1) associated with AIDS. **Special instructions for cleaning and decontamination against HBV, HCV and HIV-1 of surfaces/objects soiled with blood/body fluids: Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks, or eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces/objects before disinfection. **Contact Time:** Allow the surface to remain wet for 3 minutes prior to wiping. **Infectious Materials Disposal:** Cleaning materials used that may contain blood/body fluids must be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

This product is not to be used as a terminal sterilant or high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter sterile areas of the body. This product may be used to pre-clean or decontaminate critical and semi-critical medical devices prior to sterilization or high-level disinfection.

STORAGE AND DISPOSAL

Do not contaminate food or feed by storage or disposal.
Storage: Store in a cool, well-ventilated area. **Product and Container Disposal:** Dispose of wipe in trash after use. Do not flush. Tightly close lid between uses to retain moisture. **Nonrefillable Container:** Do not reuse or refill this empty container. Offer empty container for recycling. If recycling is not available, discard container in trash.

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals

CAUTION: May be harmful if inhaled or absorbed through skin. Causes moderate eye irritation. Avoid breathing vapors and contact with eyes, skin or clothing. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash clothing before reuse.

FIRST AID

If inhaled: Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible. **If on skin or clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. **If in eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.

Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. You may also contact 1-800-255-3924 for emergency medical treatment information.

PHYSICAL OR CHEMICAL HAZARDS

Do not use or store near heat or open flame.

Manufactured in the USA for:
 Biotrol
 13705 Shoreline Court East
 Earth City, MO 63045
 (800) 822-8550
 12002829C
 YD100-82638A-0713



No Print Area • Glue Adhesive