



STATE OF HAWAII
Department of Agriculture
ACCEPTED
LICENSE NO. **8585.4**

COEfect MinuteSpray is intended for use in:

- Hospitals and healthcare settings including: Ambulances, Anesthesia, CAT Labs, CCU, Central Supply Clinics, Dental Offices, Dialysis Clinics, Doctors' Offices, Emergency Rooms, Food Service & Food Preparation Areas, Home Health Care, Hospices, Laboratory, Laundry Rooms, Long-Term Care Facilities, Newborn Nursery, Nursing Homes, Offices, Operating Rooms, Pediatrics, Physical Therapy, Physicians' Office, Ophthalmic Offices, Radiology, Recovery Room, Rehabilitation, Respiratory Therapy, Ultrasound, Transport Vehicles, Waiting Areas, X-Ray.
- Barber Shops, Beauty Salons, Mortuaries
- Athletic Facilities, Gymsnasiums, Schools
- Transportation Terminals
- Banquet Halls, Catering Facilities
- Childcare Facilities
- Institutional Kitchens and/or Restrooms
- Food Processing & Packaging Facilities

COEfect MinuteSpray will clean and disinfect hard, non-porous surfaces of: bed railings, cabinets, carts, chairs, counters, exam tables, gurneys, isolettes, IV poles, stethoscopes, tables, telephones, toilets seats, toilets, bathtubs, bathroom fixtures, shower stalls, sinks, tiles, vanity tops, hospital beds, traction devices, MRI, CAT, scales, paddles and wheelchairs.

COEfect MinuteSpray will clean and disinfect hard, non-porous surfaces of: changing tables, cribs, diaper pails.

COEfect MinuteSpray will clean and disinfect hard, non-porous surfaces of: cabinets, chairs, computers, countertops, desks, desktops, fax machines, garbage cans, headsets, keyboards, railings, seats, sinks, tables, tabletops, appliance exteriors and faucets.

COEfect MinuteSpray will clean and disinfect hard, non-porous surfaces of: gym equipment, exercise machines, non-porous wrestling/gymnastic/exercise/yoga mats.

Use this product to clean and disinfect non-porous personal protective safety equipment, protective headgear, half mask respirators, full-face breathing apparatus, gas mask, goggles, spectacles, face shields, silicone rubber or PVC hearing protectors, and vinyl covered ear muffs. Rinse all equipment that comes in prolonged contact with skin before reuse with clean, warm water about 120F°, and allow to dry.

CAUTION: Cleaning at 120F° temperature will avoid overheating and distortion of the personal safety equipment that would necessitate replacement.

This product is formulated to disinfect hard, non-porous inanimate environmental surfaces such as floors, walls, metal surfaces, stainless steel surfaces, glazed porcelain, glazed ceramic tile, plastic surfaces, bathrooms, shower stalls, bathtubs and cabinets. This product helps to prevent cross contamination from contact with diseased surfaces.

Manufactured by
GC AMERICA INC.
3737 WEST 127TH STREET
ALSIP, IL 60803 U.S.A.
1-708-557-0900
www.gcamerica.com



COEfect™ MinuteSpray

SURFACE DISINFECTANT

One Step: Clean, Disinfect & Deodorize
One Minute: Virucidal*, Bactericidal,
Tuberculocidal, Fungicidal

ACTIVE INGREDIENTS:

Ethyl Alcohol	72.50%
Quaternary Ammonium Chloride	0.33%
OTHER INGREDIENTS	27.14%
Total	100.00%

NET CONTENTS: 1 U.S. Quart (16. 13oz)
Re-Order 557201



GC AMERICA INC.
ALSIP, IL 60803

CAUTION
KEEP OUT OF
REACH OF
CHILDREN

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. This product is not for use on critical and semi-critical medical device surfaces.

CLEANING: Wet surface sufficiently and clean with a mop, sponge, cloth, nonwoven wipe or other suitable cleaning instrument. **DISINFECTING:** To disinfect hard, non-porous surface, thoroughly wet the surface to be treated. Treated surface must remain visibly wet for one minute to achieve complete kill of all pathogens listed on this label. For visibly soiled surfaces, clean before following disinfecting instructions. Allow surfaces to air dry. ***VIRUCIDAL:** This product is an effective virucide on hard, non-porous surfaces, in the presence of an organic load, against Poliovirus serotype 1 (ATCC VR-1562), Hepatitis B Virus, Hepatitis C Virus, Rhinovirus (ATCC VR-40), Norovirus, Human Coronavirus (ATCC VR-1558), Rotavirus Strain WA, Respiratory Syncytial Virus (ATCC VR-26), Human Immunodeficiency Virus, Herpes Simplex Virus Type 1 (ATCC VR-539) and Influenza Virus type A2 (ATCC VR-544), when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

KILLS HIV-1, HBV AND HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in healthcare settings (Hospitals, Nursing Homes) or other setting in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids - and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) [associated with AIDS] Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV AND HCV FOR SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS.

PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use industry approved disposable latex gloves, gown, mask and eye protection/covering. **CLEANING PROCEDURE:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of COEfect MinuteSpray. **CONTACT TIME:** After application, allow treated surfaces to remain wet for one (1) full minute.

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

Respiratory illnesses attributed to Pandemic 2009 H1N1 are caused by Influenza A virus. This product is a broad-spectrum hard surface disinfectant that has been shown to be effective against Influenza A virus tested and listed on the label and is expected to inactivate all Influenza A viruses including Pandemic 2009 H1N1 formerly called swine flu. **BACTERICIDAL:** COEfect MinuteSpray is an effective disinfectant on hard, non-porous surfaces, in the presence of an organic load against Escherichia coli [E. coli] 0157:H7 [ATCC 35150], Pseudomonas aeruginosa [ATCC 15442], Salmonella enterica [Salmonella] [ATCC 10708], Staphylococcus aureus [ATCC 6538], Streptococcus pyogenes [Strep] [ATCC 19615], Vancomycin resistant, Enterococcus faecium [VRE] [ATCC 51559], Methicillin resistant Staphylococcus aureus [MRSA] [ATCC 33591], Klebsiella pneumoniae [ATCC 4352], Listeria monocytogenes [ATCC 9841], Acinetobacter baumannii [ACTT 19606], Campylobacter jejuni [ATCC 29428] and Burkholderia cepacia [ATCC 25416], when the treated surface is allowed to remain wet for one (1) minute. Allow surfaces to air dry. **TUBERCULOCIDAL:** COEfect MinuteSpray is an effective disinfectant on hard, non-porous surfaces, in the presence of an organic load, against Mycobacterium bovis [ATCC 35743], when the treated surface is allowed to remain wet for one (1) minute. Allow surfaces to air dry. **FUNGICIDAL:** COEfect MinuteSpray is an effective fungicide against Trichophyton mentagrophytes [ATCC 9533], causative agent of Athlete's Foot Fungus, when used as directed on hard surfaces found in bathrooms, shower stalls, locker rooms, exercise facilities or other clean, hard, non-porous, surfaces commonly contacted by bare feet. The surface must remain visibly wet for one (1) minute. Allow surfaces to air dry.

EFFICACY TESTS HAVE DEMONSTRATED THAT THIS PRODUCT IS AN EFFECTIVE VIRUCIDE, BACTERICIDE, TUBERCULOCIDE AND FUNGICIDE IN THE PRESENCE OF ORGANIC SOIL (5% BLOOD SERUM). For surfaces that may come in contact with food, a potable water rinse is required.

This product will not harm most surfaces, including stainless steel, aluminum, plastics, fiberglass, vinyl and most painted surfaces. For plastic and painted surfaces, spot test on inconspicuous area before general use. DO NOT use on unpainted wood. Not for use on people or pets.

PRECAUTIONARY STATEMENTS: Hazardous to Humans and Domestic Animals.

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

FIRST AID: IF IN EYES, hold eye(s) open and flush slowly and gently with water for 15-20 minutes. Remove contact lenses if present, after first 5 minutes, then continue rinsing. Call a Poison Control Center or doctor for treatment advice. Have the product container or label with you when calling poison control center, doctor or going for treatment. For emergency information concerning this product, call the National Pesticides Information Center (NPIC) at 1-800-858-7378 seven days a week, 6:30am to 4:30pm Pacific Time - NPIC Web site: www.npic.orst.edu

PHYSICAL OR CHEMICAL HAZARD: Flammable, keep away from heat/open flame. **STORAGE/DISPOSAL:** Do not contaminate water, food or feed by storage or disposal. **PESTICIDE STORAGE:** Always store this product and all pesticides in their original container. If a leaky container must be contained within another, mark the outer container to identify the contents. Protect from extreme heat or cold.

PESTICIDE DISPOSAL: Waste resulting from the use of this product, must be disposed of on site or at an approved waste disposal facility.

CONTAINER DISPOSAL: Nonrefillable container. DO NOT reuse or refill this container. Offer for recycling if available or discard this empty container in trash.

No warranty, either expressed or implied, is given to the usage contrary to the instructions given herein.

EPA Reg. No. 88494-1-10214
EPA Est. No. 10088-WI-2
Made in U.S.A.

