

Dieline  
will Not Print.



STATE OF HAWAII  
Department of Agriculture

**ACCEPTED**

LICENSE NO. **8754.1**

**Disc.**

30% Black

# DECON-PHASE™

## Sterile Pharmaceutical Clean Room Formula

Germicidal Detergent with Organic Soil Tolerance

**Pseudomonacidal • Tuberculocidal • Fungicidal • Virucidal \***

**ACTIVES:**

Potassium ortho-benzyl parachlorophenate . . . . .	5.03%
Potassium ortho-phenylphenate . . . . .	5.28%
Potassium para-tertiary-amylphenate . . . . .	4.30%
Inert Ingredients . . . . .	81.39%
Total . . . . .	100.00%

EPA REG. NO.211-25-68959 EPA EST. NO. 68959-PA-1

**DANGER**  
KEEP OUT OF REACH OF CHILDREN

**PRECAUTIONARY STATEMENTS**

**HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

Corrosive, Causes irreversible eye damage or skin burns. Harmful if swallowed or absorbed through skin or inhaled. Do not get in eyes, on skin, or on clothing. Avoid breathing spray mist. Wear goggles or face shield. Wear protective clothing (long-sleeve shirt and long pants, socks, plus shoes and chemical resistant gloves such as water proof gloves). Wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Remove contaminated clothing and wash before reuse.

**FIRST AID**

**IF IN EYES:** Hold eye open and rinse slowly and gently with plenty of water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

**IF INHALED:** Move person to fresh air. If person is not breathing, call 911 or an ambulance, then give artificial respiration, preferably by mouth to mouth, if possible.

**IF ON SKIN OR CLOTHING:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice.

**IF SWALLOWED:** Call a poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by a poison control center or doctor. Do not give anything by mouth to an unconscious person.

**HOTLINE NUMBER:** For emergency medical advice, call your local poison control center or a doctor. Have the product container or label with you when seeking medical advice or treatment.

**NOTE TO PHYSICIAN:** Probable mucosal damage may contraindicate the use of gastric lavage. Measure against circulatory shock, respiratory depression, and convulsion may be needed.

Product Questions? Call 1-888-478-3745 during normal business hours (8:30am - 5pm, EST)

Container and Product Distributed By:

Veltek Associates, Inc. 15 Lee Blvd. Malvern, PA 19355-1234 Tel: 610-644-8335 Fax: 610-644-8336

Made in U.S.A.



**PERSONAL PROTECTIVE EQUIPMENT (PPE)**

Handlers must wear: Coveralls over long-sleeved shirt and long pants; socks and chemical resistant shoes; goggles or face shield and chemical resistant gloves (such as barrier laminate, butyl rubber, polyvinyl chloride, or vitlon). Follow manufacturer's instructions for cleaning/maintaining Personal Protective Equipment (PPE). If no such instructions exist, use detergent and hot water. Keep and wash PPE separately from other laundry.

**USER SAFETY REQUIREMENTS**

Users should wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Users should remove PPE immediately after handling this product. Wash the outside of the gloves before removing. As soon as possible, wash thoroughly and change into clean clothing. Users should remove clothing immediately if pesticide gets inside. Then wash thoroughly and put on clean clothing.

**PHYSICAL AND CHEMICAL HAZARDS**

Combustible. Do not use or store product near heat or open flame.

**PESTICIDE STORAGE AND DISPOSAL:** Do not contaminate water, food or feed by storage or disposal.

**PESTICIDE STORAGE:** Store in securely closed original container. Avoid storage at temperature extremes or in sunlight.

**PESTICIDE DISPOSAL:** Securely wrap original container in several layers of newspaper and discard in trash.

**CONTAINER DISPOSAL:** Do not reuse container. Replace cap and discard container in trash.

**DECON-PHASE** is a germicidal detergent formulated for use where broad spectrum disinfection is required in the presence of moderate amounts of organic soil.

**DECON-PHASE** is a phosphate free, dilutable, hospital germicidal detergent effective in 10 minutes at 20° C, in the presence of 5% blood serum, used to simulate organic soil. **DECON-PHASE**'s multi-phenolic formula is designed to clean, disinfect and deodorize any washable, inanimate, non-porous surface in one easy step. When used as directed, **DECON-PHASE** is strong enough to kill a broad spectrum of pathogenic bacteria, yet is mild enough to have no harmful effects on the surface being disinfected, when used as directed.

**DECON-PHASE** is recommended for use by pharmaceutical and biotechnology industries; medical device manufacturers; hospitals and health care institutions that are dedicated to controlling the hazards of cross contamination in their establishments. For use on hard, inanimate surfaces in aseptic filling and gowning rooms, general manufacturing areas, or on machinery, tables, counters, laminar flow benches, floors, walls, stainless steel, porcelain, glass, and chrome.

**DIRECTIONS FOR USE**

It is a violation of federal law to use this product in a manner inconsistent with its labeling.

**CLEANING AND DISINFECTION**

For broad spectrum disinfection of gram negative and gram positive bacteria, including Pseudomonas aeruginosa, add 1/2 ounce **DECON-PHASE** to one gallon of water. Remove heavy soil or gross filth then apply solution with a mop, cloth or sponge so as to wet the surface thoroughly. Do not apply this product in a way that will contact workers or other persons. Only protected handlers may be in the area during application. Allow to remain wet for 10 minutes and then let air dry. If higher detergency is desired, increase dilution to 1 ounces per gallon of water. Prepare a fresh solution for each use.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

**EFFICACY PERFORMANCE**

The broad spectrum efficacy of **DECON-PHASE** has been confirmed using AOAC Protocol testing at a use-dilution of 1:256 in the presence of 5% blood serum (used to simulate organic soil load) at 20 degrees Celsius with an exposure time of 10 minutes.

**GERMICIDAL:**

**DECON-PHASE** is effective against a broad spectrum of gram positive and gram negative pathogens as confirmed against Salmonella choleraesuis, Staphylococcus aureus, and Pseudomonas aeruginosa by the AOAC Use Dilution test. In the presence of 5% organic soil load at 20 degrees Celsius with an exposure time of 10 minutes.

**FUNGICIDAL:**

At 1/2 ounce per gallon of water, **DECON-PHASE** is also effective against Trichophyton mentagrophytes in use areas as aseptic filling and gowning rooms, general manufacturing areas or on any other hard, non-porous inanimate environmental surfaces. At 1/2 ounce per gallon of water, **DECON-PHASE** is effective against Candida albicans.

**TUBERCULOCIDAL:**

**DECON-PHASE** disinfectant's effectiveness against Mycobacterium bovis (BCG) is demonstrated through testing using the AOAC Method "Confirmative In Vitro Test for Determining Tuberculocidal Activity" in 10 minutes at 20 degrees Celsius.

**\*VIRUCIDAL \*KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of the Human Immunodeficiency Virus Type-1 (HIV-1) (associated with AIDS).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 (AIDS VIRUS) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.**

**PERSONAL PROTECTION:** Clean-up should always be done wearing protective latex gloves, gowns, masks, and eye protection.

**CLEANING PROCEDURES:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

**CONTACT TIME:** Leave surface wet for 10 minutes with a, 1/2 ounce per gallon use-dilution (1:256), **DECON-PHASE** and water.

**DISPOSAL OF INFECTIOUS WASTE MATERIALS:** Blood, body fluids, cleaning materials, and clothing should be autoclaved and disposed according to local regulations for infectious waste disposal.

**PELIGRO:** Si no puede leer en ingles, pregunte a su supervisor sobre las instrucciones de uso apropiadas antes de trabajar con este producto.

MSDS #: DPH-98-01

REV10/05 DPH-02

Colors: Black

Run Direction →

### IMPORTANT

Please check ENTIRE proof carefully including ALL new, revised, repeat copy AND human readable barcode characters. Refer to Pantone Color Chips for accurate PMS colors. All colors print 100% unless noted.

Proof # 3	Date: 3/17/06	JH/ak/jc
Ticket # 46153	Product # 28999	
Customer: Veltek Associates, Inc.		
Die Size: 5" x 11"		
(.125" corner radius)		
Colors: black		

Alterations to supplied art and re-proofs will incur additional charges. This proof is provided for your inspection of copy, layout, finished size and color breakdown. Please review carefully & indicate instructions below.  
 Proof approved for printing as is. Proceed with order.  
 Proof not approved. Make corrections as indicated & reproof.

Authorized Signature

Date

RECEIVED

JAN - 6 2010

Pesticides Branch