



Department of Agriculture  
STATE OF HAWAII

LICENSED

PERIOD 2011-2013 LIC. NO.



Department of Agriculture  
STATE OF HAWAII

LICENSED

PERIOD 2012-2014 LIC. NO. **8813.1**

ARJOHUNTLEIGH  
GETINGE GROUP

Surface Disinfectant Wipes are effective against the following pathogens in 5 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed:  
**Trichophyton mentagrophytes (Athlete's foot fungus)**

**Areas of use:**

hospital	surgery	exam rooms
intensive care units	emergency medical services	X-ray
emergency rooms	respiratory therapy	Cat. Lab
patient rooms	healthcare facilities	laboratories
medical offices	physician offices	ambulances

**PREPARATION OF DISPENSER:** Remove cover from canister and discard covering seal. From center of wipes roll, pull up a wipe corner, twist into a point and thread through the tear throat in the underside of the lid cover. Pull wipe through about one (1) inch. Replace the lid cover on container, assuring a secure fit. Pull out first wipe and snap at a 90° angle. Remaining wipes feed automatically for next use. Always reseal hinged cap after each use to insure moisture retention.

**DIRECTIONS FOR USE:** It is a violation of U.S. Federal law to use this product in a manner inconsistent with its labeling.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection. This product is not for use on the inhalation therapy and/or pulmonary diagnostic equipment.

**DISINFECTION:** Thoroughly wet pre-cleaned, hard, non-porous surface with a wipe, keep wet for 2 minutes (5 minutes if fungus is suspected), and allow to air dry. Use as many wipes as needed for the treated surface to remain wet for the entire contact time. For use as a general disinfectant wipe on hard, non-porous, inanimate surfaces.

You may use a wipe to pre-clean the surface to be disinfected, but through re-wiping with a fresh wipe is required for disinfection. Use on food contact surfaces, must be followed by a potable water rinse.

**DEODORIZING:** Wipe hard, non-porous surface thoroughly to eliminate odors from bathrooms, fecal matter, urine, burns, smoke, mildew and mold, decay, stale air, tobacco, sickrooms and musty closets.

**CLEANING:** To clean surfaces, wipe thoroughly with a fresh wipe, which may be used to deodorize at the same time as cleaning. When used to clean a hard non-porous surface prior to disinfection, a second wipe must be used for disinfection.

**THE ACTIVITY OF Surface Disinfectant Wipes AGAINST THE HIV-1 (AIDS VIRUS):** Surface Disinfectant Wipes completely inactivate the HIV-1 (AIDS virus) on hard, non-porous surfaces when used at a contact time of 1 minute at room temperature (20°-25°C). Surface Disinfectant Wipes kill HIV-1 on pre-cleaned environmental surfaces/objects soiled with blood/body fluids, in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type-1 (HIV-1)(associated with AIDS).

**FOR USE ON ENVIRONMENTAL SURFACES**

Virucidal*	<b>ACTIVE INGREDIENTS:</b> N-Alkyl (60% C14, 30% C16, 5% C12, 5% C18)..... 0.12% Dimethyl Benzyl Ammonium Chloride N-Alkyl (68% C12, 32% C14)..... 0.12% Dimethyl Ethyl Benzyl Ammonium Chloride
Bactericidal	Isopropyl Alcohol ..... 41.58%
Pseudomoniacidal	
Tuberculocidal	
Fungicidal	
	<b>OTHER INGREDIENTS+:</b> ..... 58.18%
	<b>TOTAL:</b> .....100.00%

+ Does not include weight of dry wipe

PERIOD 2011-2013 LIC. NO. PERIOD 2012-2014 LIC. NO. **8813.1**

# SURFACE

## Disinfectant Wipes

**KEEP OUT OF REACH OF CHILDREN**

**WARNING:** SEE BACK PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS.

**DO NOT USE AS A BABY WIPE**

Surface Disinfectant Wipes are effective against the following pathogens in 2 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed:

<i>Staphylococcus aureus</i>	<i>Pseudomonas aeruginosa</i>
<i>Mycobacterium tuberculosis</i>	Methicillin Resistant <i>Staphylococcus aureus</i> (MRSA)
<i>Escherichia coli</i> O157:H7	Vancomycin Resistant <i>Enterococcus faecalis</i> (VRE)
* Human Hepatitis B virus	* Influenza A (A2)Japan)
* Human Hepatitis C virus	* HIV-1 (AIDS virus)**(At an exposure time of 1 min.)
<i>Salmonella enterica</i>	* Herpes simplex virus 2(genital herpes virus)

\*\* See back panel for complete HIV-1 (AIDS virus) information

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:** **Personal Protection:** Specific barrier protection items to be used when handling items soiled with blood or body fluids are disposable latex gloves, gowns, masks and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of **Surface Disinfectant Wipes.**  
**CONTACT TIME:** HIV-1 virus is inactivated in 1 minute, **however, other pathogens require 5 minutes exposure for complete kill. See disinfection directions.**  
**Disposal of Infectious Materials:** Materials that may contain blood or other body fluids should be autoclaved and disposed of according to federal, state and local regulations for infectious waste disposal.

**STORAGE AND DISPOSAL:**  
Store in cool area. Keep away from heat or open flame. Do not reuse wipes. Dispose of used wipe in trash. **Container Disposal:** Non-refillable container. Do not reuse or refill this container. Offer for recycling if available.

**PRECAUTIONARY STATEMENTS:**  
Hazards to Humans and Domestic Animals  
**WARNING:** Causes substantial, but temporary eye injury. Do not get in eyes or on clothing. Wear (specific appropriate protective eye wear such as goggles, face shield or safety glasses). Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, or using tobacco. Wear disposable latex gloves, gowns, masks and eye coverings, when handling HIV-1 infected blood or body fluids. Avoid contamination of food. Remove and wash contaminated clothing before reuse.

**FIRST AID:**  
**If in eyes:** Hold eye open and rinse slowly and gently with water for at least 15-20 minutes. Remove contact lenses, if present, after first 5 minutes, then, continue rinsing eye. Call a poison control center or doctor for treatment advice. **If on skin or clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15 to 20 minutes. Call a poison control center or doctor for treatment advice. Wash contaminated clothing before reuse.  
**WASH THOROUGHLY AFTER HANDLING WIPES.**

Distributed in the USA by ArjoHuntleigh Inc.  
2349 West Lake Street, Suite 250  
Addison, Illinois 60101  
EPA Reg. No. 59894-10 -45556  
EPA Est. No. 73027-WI-01

For product emergency,  
Call CHEMTREC: 1-800-424-9300

For an MSDS, call 800-323-1245

**NET CONTENTS: 50 Pre-saturated Wipes (5x8 in) | NET WEIGHT: .50 lbs**