

DRAFT

Note to Customer: Weiman Products adds a revision code to all artwork items for tracking purposes. As the customer & controller of your own artwork, we ask that you either denote a specific area for this revision code and/or ask Weiman for the specific revision code that is to be added & place on artwork.



LICENSED

14.5 inches Width

LICENSED

PERIOD 2014-2016 LIC. NO. **9946.7**

PERIOD 2011-2013 LIC. NO.

Sklar Disinfectant Surface Wipes effectively kills the following microorganisms at room temperature with a 3-minute contact time: *Acinetobacter baumannii*, Avian Influenza A Virus*, Hepatitis B Virus* (HBV), Hepatitis C Virus* (HCV), Herpes Simplex II Virus* (HSV 2) G Strain, Human Immunodeficiency Virus* (HIV-1) (associated with AIDS), Influenza Virus* Strain A2/Hong Kong, Methicillin-Resistant *Staphylococcus aureus* (MRSA), *Mycobacterium bovis*, BCG (tuberculosis) (TB), Polio I Virus*, *Pseudomonas aeruginosa*, Rhinovirus*, Rotavirus* Strain WA, *Salmonella enterica*, *Staphylococcus aureus* and Vancomycin-Resistant *Enterococcus faecium* (VRE).

This product has demonstrated effectiveness against Influenza A Virus* and is expected to inactivate all Influenza A Viruses including Pandemic 2009 H1N1 Influenza virus*.

Directions for Use

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

DO NOT USE THIS PRODUCT AS A BABY OR SKIN WIPE, OR ON CONTACT LENSES!

Sklar Disinfectant Surface Wipes are disposable, single-use wipes pre-saturated with Sklar disinfectant solution.

When used as directed Sklar Disinfectant Surface Wipes clean and disinfect surfaces found in blood centers, cardiac catheterization labs, dialysis units, esthetician facilities, hospitals, medical offices, surgical centers, laboratories, intensive care units, isolation rooms, neonatal units, nursing homes, ophthalmic centers, EMS units, skin care facilities, respiratory care units, and areas/facilities concerned with the hazards of cross contamination.

May be used on devices/surfaces potentially contaminated with disease-causing microorganisms and use sites made of glass, metal, Plexiglas®, plastics (i.e., polycarbonate, polypropylene, polyvinylchloride and polystyrene), baked-on painted surfaces, stainless steel and non-porous vinyl.

Dispensing Directions: To start feed, remove cover and discard seal. From the center of the wipe roll, pull up a wipe corner, twist it into a point and thread it through the hole located on the container cover. Pull through about one inch. Replace cover. Pull out first wipe and tear off at an angle. Remaining wipes feed automatically, ready for the next use. When not in use, keep center cap closed to prevent evaporation.

Cleaning and Disinfection Instructions: Use a fresh Sklar Disinfectant Surface Wipes to pre-clean surfaces of all gross filth and heavy soil. Repeat as necessary until surfaces are visibly clean. To effectively disinfect the pre-cleaned surfaces, use a new Sklar Disinfectant Surface Wipes to thoroughly wet the surfaces and allow surfaces to remain wet for 3 minutes at room temperature. Use enough wipes to ensure that the surface remains wet for the three-minute contact time. Allow surfaces to air dry before use, or if desired, rinse surface after 3 minutes using potable water and wipe or air dry. For disinfection of Avian Influenza A, use this product only for surfaces that are not conducive to treatment by immersion or excess liquid (e.g., equipment tables, lab benches and similar "high touch" surfaces). Food preparation surfaces must be rinsed with potable water prior to use.

Distributed by: Sklar Instruments
889 S. Matlack Street • West Chester, PA 19382
800.221.2166 • 610.430.3200 • fx: 610.696.9007
Made in the U.S.A.

SW06100-81513BMS-0711



Sklar Disinfectant Surface Wipes

Bactericidal • Virucidal*
Tuberculocidal Disinfectant / Cleaner

ACTIVE INGREDIENTS:	
n-Alkyl (60% C ₁₄ , 30% C ₁₆ , 5% C ₁₂ , 5% C ₁₈)	
dimethyl benzyl ammonium chloride.....	0.154%
n-Alkyl (68% C ₁₂ , 32% C ₁₄)	
dimethyl ethylbenzyl ammonium chloride.....	0.154%
OTHER INGREDIENTS	99.692%
TOTAL	100.000%

KEEP OUT OF REACH OF CHILDREN
CAUTION

See side panel for additional Precautionary Statements.

Cat. No. 10-1616

Net Contents: 100 Wipes

Net Weight: 1 lb. 8 oz. (680.4 g)

EPA REG. NO. 70144-2-31118

EPA EST. 1130-IL-1

HBV, HCV and HIV-1 Precautions: Sklar Disinfectant Surface Wipes effectively inactivates HBV, HCV and HIV-1 on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids and in which the surfaces/objects can be associated with the potential for transmission of HBV, HCV and/or the human immunodeficiency virus Type 1 (HIV-1) associated with AIDS. **Special instructions for cleaning and decontamination against HBV, HCV and HIV-1 of surfaces/objects soiled with blood/body fluids:** **Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks or eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces/objects before disinfection. **Contact Time:** Allow the surface to remain wet for 3 minutes prior to wiping. **Infectious Materials Disposal:** Cleaning materials used that may contain blood/body fluids should be autoclaved and/or disposed of in accordance with local regulations for infectious materials disposal.

This product is not to be used as a terminal sterilant or high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter sterile areas of the body. This product may be used to pre-clean or decontaminate critical and semi-critical medical devices prior to sterilization or high-level disinfection.

STORAGE AND DISPOSAL

Do not contaminate food or feed by storage or disposal.

Storage: Store in a cool, well-ventilated area. **Product and Container Disposal:** Dispose of wipe in trash after use. Do not flush. Tightly close lid between uses to retain moisture. **Nonrefillable Container:** Do not reuse or refill this empty container. Offer empty container for recycling. If recycling is not available, discard container in trash.

Precautionary Statements
Hazards to Humans and Domestic Animals

CAUTION: Harmful if inhaled or absorbed through skin. Causes moderate eye irritation. Avoid breathing vapors and contact with eyes, skin or clothing. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash clothing before reuse.

FIRST AID: If inhaled: Move person to fresh air. If person is not breathing call 911 or an ambulance, then give artificial respiration, preferably by mouth-to-mouth, if possible. **If on skin or clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. **If in eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.

Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. You may also contact 1-800-255-3924 for emergency medical treatment information.

Physical or Chemical Hazards
Do not use or store near heat or open flame.

6 inches Height

0.25

0.25