

This product is to be diluted only with the 3M™ Twist 'n Fill™ Dispenser.

DESCRIPTION

- Germicidal detergent for use in hospitals, medical & dental offices and veterinary clinics.
- Use on hard, non-porous surfaces, such as floors or walls.
- Disinfects, cleans and deodorizes in one labor-saving step.

Net Contents/Contenido Neto:

0.528 U.S. Gallon • 2 Liters
0.528 Galón de E.U.A. • 2 Litros

Directions for Use (continued)

FUNGICIDAL: This product is an effective fungicide against *Trichophyton mentagrophytes* (agent causes athlete's foot fungus) and *Candida albicans* in areas such as locker rooms, dressing rooms, shower and bath areas. Follow disinfection directions. Refer to product label hang tag for additional directions for use.

PRECAUTIONARY STATEMENTS HAZARDS TO HUMANS AND DOMESTIC ANIMALS

DANGER. Corrosive. Causes irreversible eye damage and skin burns. Harmful if swallowed or absorbed through the skin. Do not get in eyes, on skin or clothing. Avoid breathing spray mist.

PERSONAL PROTECTIVE EQUIPMENT (PPE)

Applicators and other handlers must wear: Coversalls over long-sleeved shirt and long pants, socks and chemical-resistant footwear, goggles or face shield, and chemical resistant gloves (such as barrier laminate or butyl rubber, nitrile rubber, neoprene rubber, polyvinyl chloride, or viton at ≥ 14 mils).

USER SAFETY REQUIREMENTS: Follow the manufacturer's instructions for cleaning/maintaining PPE. If no such instructions

18L



18L



18L

Phenolic Disinfecting and Cleaning Concentrate



KEEP OUT OF REACH OF CHILDREN
DANGER

ACTIVE INGREDIENTS	0.475%
o-benzyl-p-chlorophenol	4.75%
o-phenylphenol	4.75%
OTHER INGREDIENTS:	90.50%
Total	100.00%

This product contains the active ingredients 0.475% o-benzyl phenol and 4.75% o-phenyl phenol.

Use in surfaces, exterior bow windows, walls, and hard nonporous surfaces made of metal. For use in hospitals, nursing homes, schools, colleges, medical and dental offices, farms, tack shops, kennels, and animal care facilities. For use in cleaning and disinfecting and cleaning Concentrate is effective against the following Pathogenic Organisms: Pseudomonas aeruginosa, Staphylococcus aureus, Salmonella enteric, Enterococcus spp. ** Staphylococcus aureus - Methicillin resistant (MRSA), Methicillin Type V Virus/Variant, kills pandemic 2009 H1N1 Influenza A virus. Refer to product label hang tag for additional pathogens. **

DIRECTIONS FOR USE (continued)
TO DISINFECT: Use 1 oz of 3M Phenolic Disinfectant Cleaner Concentrate per gallon of water for a minimum contact time of 10 minutes in a single application. For disinfecting, remove gross filth and heavy soil deposits, then thoroughly wet surfaces. This product is extremely versatile and can be applied with a mop, sponge, or cloth as well as soaking or used as a holding solution. Thoroughly wet surfaces with the use solution for at least 10 minutes. The recommended use solution is used once and discarded on a daily basis. Rinsing is not necessary on floor surfaces unless floors are to be waxed or polished. Tuberculocidal effectiveness is in 10 minutes at 20°C.
BACTERICIDAL STABILITY OF USE-DILUTION: Tests confirm that this product, when diluted in 400 ppm hard water remains effective against *Pseudomonas aeruginosa*, *Staphylococcus aureus* and *Salmonella choleraesuis* for up to one year when stored in a sealed container such as a spray bottle. If product becomes visibly dirty or contaminated, the use-dilution must be discarded and fresh product prepared. Always use clean, properly labeled containers when diluting this product. Bactericidal stability of the use-dilution does not apply to open containers such as buckets or pails.

PPE Use (continued)

For washables, use detergent and hot water. Keep and wash PPE separately from other laundry. Users must wash hands before eating, drinking, chewing gum, using tobacco, or using the toilet. Users must remove clothing immediately if pesticide gets inside. Then wash thoroughly and put on clean clothing. Users must remove PPE immediately after handling this product. Wash the outside of gloves before removing. As soon as possible, wash thoroughly and change into clean clothing.

FIRST AID: If in Eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. **If on Skin or Clothing:** Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes. Call a poison control center or doctor for treatment advice. **If Swallowed:** Call a poison control center or doctor immediately for treatment advice. Have person sip a glass of water if able to swallow. Do not induce vomiting unless told to do so by a poison control center or doctor. Do not give anything by mouth to an unconscious person. **HOT LINE NUMBER:** In case of emergency call toll free 1-800-364-3577 (1-800-3M-HELPS). Have the product container or label with you when calling a Poison Control Center or doctor or going for treatment. **NOTE TO PHYSICIAN:** Probable mucosal damage may contraindicate the use of gastric lavage. Measures against circulatory shock, respiratory depression, and convulsion may be needed.

REFER TO MATERIAL SAFETY DATA SHEET (MSDS: 07-3145-5) FOR MORE INFORMATION. IN CASE OF EMERGENCY, CALL: (651) 737-6501 OR 1-800-364-3577 (1-800-3M-HELPS). See product label hang tag for additional product use information.

STORAGE AND DISPOSAL: Do not contaminate water, food or feed by storage or disposal. **Pesticide Storage:** Open dumping is prohibited. Store in original container in areas inaccessible to children. **Pesticide Disposal:** Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your State Pesticide or Environmental Control Agency, or the Hazardous Waste representative at the nearest EPA Regional Office for guidance. **Container Disposal:** Nonrefillable container. Do not reuse or refill this container. Clean container promptly after emptying. Triple rinse as follows: Fill container ¼ full with water and reclose the container. Agitate vigorously, and dispose of rinsate consistent with pesticide disposal instructions. Repeat two more times. Then offer for recycling if available or puncture and dispose in sanitary landfill or by other procedures approved by state and local authorities. Follow pesticide disposal instructions for rinsate. If not triple rinsed, these containers are acute hazardous wastes and must be disposed in accordance with local, state, and federal regulations. ***Typical yield

3M and exported by:
3M
3M Building and Commercial Services Division
St. Paul, MN 55144-1000
www.3M.com/building

Made in U.S.A.
© 3M 2011. All rights reserved.

3M and Twist 'n Fill are trademarks of 3M.
34-8706-8156-5-R

Phenolic Disinfecting and Cleaning Concentrate

18L

Makes 60 ready-to-use gallons.***

Net Contents/Contenido Neto: 0.528 U.S. Gallon • 2 Liters / 0.528 Galón de E.U.A. • 2 Litros

34-8706-8156-5-R